

Technical Documentation IPM

Evita 4



Warning

All servicing and/or test procedures on the device require detailed knowledge of this documentation. Use of the device requires detailed knowledge and observance of the relevant Instructions for Use.

Intensive Care Ventilator

This page has been left blank intentionally.

Contents

1	General.....	4
1.1	General notes	4
2	Abbreviations	7
2.1	List of abbreviations and glossary of terms	7
3	Function descriptions	8
3.1	General information	8
3.2	Basic principle.....	8
3.3	CPU 68332 PCB.....	10
3.4	CO ₂ Carrier PCB	13
3.5	CO ₂ measurement.....	16
3.6	Power supply unit	20
3.7	Evita 4 Graphics Controller PCB	21
3.8	Communication PCB (option)	23
3.9	Paediatric Flow PCB (optional).....	25
3.10	IFCO Carrier PCB (optional nurse call and other)	26
3.11	Pneumatic Controller PCB.....	27
3.12	HPSV Controller PCB	31
3.13	Pneumatic assembly.....	32
3.14	Gas mixture	51
3.15	Pneumatics function diagram	53
4	Parts catalog	55
5	Test instructions	163
5.1	Test Instructions / Service Card IPM	163
5.2	Result Sheet Test instructions / Service Card IPM	219

1 General

This chapter contains general notes and definitions that are important for the use of this documentation.

1.1 General notes

1.1.1 Notes on use

Read through the following notes thoroughly before applying this documentation.

The warnings set out here apply to all parts of this documentation.

Dräger reserves the right to make changes to the device and/or to this documentation without prior notice. This documentation is intended solely as an information resource for maintenance personnel or technical specialists.

1.1.2 Copyright and other protected rights

The content of this documentation, in particular its design, text, software, technical drawings, configurations, graphics, images, data and their selection and its composition and any amendments to it (content) are protected by copyright. The content must not (in whole or in part) be modified, copied, distributed, reproduced, republished, displayed, transmitted or sold without the prior written consent of Dräger.

1.1.3 Definitions

WARNING

An important advisory indicating a potentially hazardous situation which may result in death or serious injury if not prevented.

CAUTION

An important advisory indicating a potentially hazardous situation which may result in minor or moderate injury to the user or patient or in damage to the medical product or other assets if not prevented.

NOTE

A NOTE provides additional information intended to avoid inconvenience during operation and/or servicing.

Term	Definition
Maintenance	Maintaining the operative condition of a medical product by suitable means
Inspection	Assessment of the actual condition of a medical product
Servicing	Maintaining the operative condition of a medical product by recurrent, specified measures
Repair	Restoring the operative condition of a medical product after failure of a device function

1.1.4 General safety precautions

Read through each section thoroughly before beginning servicing.

CAUTION

Incorrect use of tools

The device's function may be impaired, or the device may be damaged.
Always use the correct tools and the specified test equipment.

WARNING

The device must be regularly inspected and serviced by maintenance personnel. Repairs and complex maintenance work on the medical product must be carried out by qualified specialists.

If you require a service contract, or for any necessary repair work, Dräger recommends DrägerService. Dräger recommends using original Dräger parts for servicing.

If the aforementioned instructions and recommendations are ignored, the correct functioning of the medical product may be put at risk. Pay attention to the "Servicing" section of the Instructions for Use.

WARNING

Non-conforming test values

If test values do not conform to specifications, the safety of the patient may be put at risk.

- Do not put the device into operation if test values do not conform to specifications.
- Contact your local service organization.

WARNING

Impermissible modifications to the device

If impermissible modifications are made to the device, the safety of the patient may be put at risk.

Do not modify the device without Dräger's permission.

WARNING

Risk of infection

The unit may transmit pathogens following use on the patient.

- Before carrying out any servicing, ensure that the device and its components have been handed over by the user cleaned and disinfected.
- Service only cleaned and disinfected units and unit components.

WARNING

Risk to patient.

Ensure that no patient is connected to the device before starting maintenance or repair work.

NOTE

Where reference is made to legislation, regulations and standards, in respect of devices used and serviced in Germany they are based on the laws of Germany. Users and technicians in other countries must comply with their national laws and/or international standards.

2 Abbreviations

This section contains a list of the abbreviations used in this document.

2.1 List of abbreviations and glossary of terms

Abbreviation	Term
ATC	Automatic Tube Compensation
BTPS	Measured values referring to conditions of the patient's lung, body temperature 37 °C, water vapour saturated gas, ambient pressure (Body Temperature, Pressure Saturated)
CAN	Fast serial port (Controller Area Network)
COM	Serial port
DUART	Dual Universal Asynchronous Receiver/Transmitter
EGB	Electrostatic sensitive devices
ECG-SYNC	ECG Synchronization
ESD	Electrostatic Sensitive Device
FiO ₂	Inspiratory O ₂ concentration
HPSV	High Pressure Servo Valve
ILV	Independent Lung Ventilation
NIF	Negative Inspiratory Force
NIV	Application mode for support of non-invasive ventilation therapies
NTPD	Normal Temperature, Pressure Dry
Paw	Airway pressure
PIP	Peak Inspiratory Pressure
RSB	Rapid Shallow Breathing (quotient of spontaneous breathing frequency and breathing volume)
SpO ₂	Functional oxygen saturation
VTASB	Tidal volume ASB

3 Function descriptions

This chapter contains descriptions of the device's technical functions.

3.1 General information

The Evita 4/ Evita 4 edition is a time-cycled, constant-volume long-term ventilator for adults and children. The features and ventilation modes depend on the specific device and its optional features; they are described in the Instructions for Use of the specific device.

3.1.1 Evita 4 edition

New design:

- The background foil is gray.
- The alarm silence key is yellow.
- The back of the display housing is blue.

New default option - "SW 4.0 Plus":

- Weaning parameters RSBi, NIF and VTASB.
- Measurement of an external flow source.
- Large loops.

New software version SW 4.20:

- Extended alarm management for battery, flow sensor and MediBus.
- Optimized ATC[®], NIV and NeoFlow.

By the integrated "Software 4.0. Plus" option loops can be displayed in full-screen format. A cursor can be used to determine individual values in the loop sequence. This facilitates an accurate examination of the lung mechanics.

Rapid shallow breathing index, negative inspiratory force index and VTASB contribute to an improved weaning process. Optional features such as ATC[®], PPS, NIV and NeoFlow can be added.

Evita 4 edition offers both conventional and non-invasive ventilation. Due to its new type of dynamic leakage compensation, the Evita 4 edition is ready to meet mask ventilation requirements.

3.2 Basic principle

3.2.1 Evita 4

The Evita 4 consists of three components which communicate via a CAN (fast serial interface).

- Control unit (Fig. 1/1)
- Electronics (Fig. 12)
- Pneumatics (Fig. 1/3)




Fig. 1 Basic components of the Evita 4

3.2.1.1 Control unit

The control unit is the interface between the device and the operator. The control unit serves to make adjustments, to display measured values and to generate alarms. In the control unit the display, keypad, touch-screen and Graphics Controller PCB are accommodated.

3.2.1.2 Electronic assembly

The electronics is the central control unit of the Evita 4. The electronics includes the CPU 68332 PCB, the CO₂ Carrier PCB with the Processor Board PCB and the Power Supply PCB and the power supply unit (Communication PCB, Paediatric Flow PCB, IFCO PCB and the SpO₂ PCB as option).

3.2.1.3 Pneumatic assembly

The pneumatic assembly controls the pneumatic valves following preset ventilation parameters. It includes an independent microprocessor system and the valve control. In the pneumatics the Pneumatics Controller PCB, the HPSV Controller AIR/O₂ PCB, the PEEP valve, the mixer, the pressure connection, the flow sensor and the O₂ sensor are accommodated.

3.2.2 Simplified block diagram of the Evita 4




Fig. 2 Basic principle of the Evita 4

Item	Designation	Item	Designation
A	Front panel.	12	Battery (Goldcap capacitor).
B	Electronics.	13	Supply voltages.
C	Pneumatics.	14	Power switch.
1	Keys.	15	Second inspiratory Paw sensor.
2	Control knob (press to confirm).	16	Pneumatics processor and venting reset.
3	Touchscreen.	17	Electronic processor reset and second loud-speaker alarm.
4	TFT display 640 x 480.	18	Inspiratory Paw sensor.
5	Information LEDs and Alarm LEDs.	19	O ₂ sensor.
6	CAN bus.	20	FiO ₂ (HPSV mixer).
7	Graphics processor reset.	21	AIR (HPSV mixer).
8	N/A.	22	Flow sensor.
9	Loudspeaker with sound chip.	23	Expiratory valve with PEEP.
10	Second loudspeaker (piezo).	24	Expiratory Paw sensor.
11	Voltage monitoring (activates reset of the processors and the piezo).		

3.3 CPU 68332 PCB

The CPU 68332 PCB is located in the electronic assembly of the Evita 4. The board includes an independent processor system, two external interfaces, three internal interfaces, the loudspeaker control and a serial EEPROM.




Fig. 3 Block diagram CPU 68332 PCB

3.3.1 EEPROM

The EEPROM is connected to the synchronized, serial interface of the 68332. The EEPROM characterizes the Evita 4 (enabled options, serial number, etc.).

3.3.2 Processor system

The processor system comprises a 68332 CPU, a 512 kByte RAM and a 1 MByte flash EPROM (electrically programmable and erasable read-only memory). The RAM has a battery back-up. When the battery is being replaced a Goldcap capacitor ensures voltage supply of the RAMs.

3.3.3 RS232 interface

The CPU 68332 PCB provides an RS232 interface in the device. The interface is labeled COM1. The interface is electrically isolated from the device. Electrical isolation is made by means of optocouplers.

3.3.4 ILV interface

The ILV interface is required for independent-lung ventilation with units. The ILV interface is not electrically isolated. Pin 3 of the ILV interface is provided with a filler plug. This filler plug prevents confusion with the RS232 interface.

3.3.5 Driver

The driver adjusts the access times between the 68332, the clock and the DUART (Dual Universal Asynchronous Receiver/Transmitter).

3.3.6 Clock

The clock gives the current time. When the device is off, the clock runs on battery back-up.

3.3.7 DUART

The DUART has two serial interfaces and digital inputs and outputs. The serial ports are used to connect the SpO₂ module and the CO₂ module.

3.3.8 DC/DC converter

The DC/DC converter provides the voltage supply (+5 V ISO) required for the interface. The input voltage of the DC/DC converter is +5 V.

3.3.9 CAN

The CAN interface is a fast, serial interface (Controller Area Network). The control unit, the electronics and the pneumatics communicate via a CAN interface. The transmission rate is 800 kbit/s.

3.3.10 Bus driver

The address bus, the data bus and the check-back signals are transferred by the bus driver to the motherboard. The 68332 CPU communicates with the optional PCBs on the motherboard via the bus driver. Currently, it is only the Paediatric Flow PCB (NeoFlow option).

3.3.11 Sound generator


The sound generator controls the loudspeaker in the control unit. The sound generator incorporates the volume control and sound generation for the loudspeaker. The volume is controlled by the DUART.

3.3.12 Reset logic

The CPU 68332 can reset the control unit and the pneumatics. A reset is also triggered if there is an undervoltage or overvoltage of the +5 V voltage. The pneumatics can also reset the CPU 68332 PCB. The reset logic controls and displays the resets.

3.4 CO₂ Carrier PCB

The CO₂ Carrier PCB is integrated in the electronic assembly of the device. The PCB includes the mount and the electrical isolation of the CO₂ module and the SpO₂ module, the power failure logic, the temperature measurement and the voltage monitoring.



* Optional

Fig. 4 Block diagram of the CO₂ Carrier PCB





Fig. 5 CO2 Carrier PCB block diagram, (part 1)

Fig. 6 CO₂ Carrier PCB block diagram, (part 2)

3.4.1 Electrical isolation

Electrical isolation is made by means of optocouplers. The PCB is coupled with the SpO₂ and CO₂ modules by plug-in contacts. The **X3** connector is part of the optional SpO₂ module and is not fitted.

3.4.2 Power failure logic

The power failure logic monitors the power supply. In the event that a mains voltage failure occurs while the device is operating an audible alarm will sound.

3.4.3 Voltage monitoring

The voltages ± 15 V, +10 V, +24 V, +12 V and +5 V are monitored. The relevant voltages are transmitted to a voltage divider and read out by an A/D converter. The A/D converter is read out by the CPU 68332 PCB.

3.4.4 Measurement of temperature and pressure

The temperature is measured by an NTC thermistor detector. The signals are processed in a temperature hybrid. The output signal of the temperature hybrid is transmitted to an A/D converter. The A/D converter is read out by the CPU 68332 PCB.

3.5 CO₂ measurement

The CO₂ measuring system comprises three modules:

- Sensor,
- Processor Board PCB,
- Power Supply PCB.




Fig. 7 Block diagram of the CO₂ measurement

3.5.1 CO₂ sensor

The CO₂ sensor comprises the CO₂ measuring unit and a microprocessor system. A lamp generates a light spectrum up to 4.5 μm . This light spectrum is transmitted via the cuvette and two sapphire lenses to detectors. The detectors emit electrical signals depending on the CO₂ concentration. These signals are analyzed by the microprocessor and transmitted to the Processor Board PCB via an RS232 interface. The CO₂ measuring unit is kept at a constant temperature to avoid condensation.




Fig. 8 Sectional view of the CO2 sensor




Fig. 9 Block diagram of the CO₂ sensor

3.5.2

Processor Board PCB

The Processor Board PCB controls the heating of the CO₂ measuring unit and the measured-data transfer of the CO₂ sensor.




Fig. 10 Processor Board PCB block diagram

3.5.3 Power Supply PCB

The Power Supply PCB provides the supply voltages for the Processor Board PCB and the CO₂ sensor. The supply voltages are electrically isolated from the device. The CO₂ measured values are transmitted to the CPU 68332 PCB via the Power Supply PCB.




Fig. 11 Power Board block diagram

3.6 Power supply unit

The power supply unit delivers the following output voltages:

- +24 V
- +15 V
- -15 V
- +12 V
- +5 V

The power supply unit is a switched-mode power supply unit. The output voltages are short-circuit-proof. A DC module is optionally available ensuring voltage supply in case of power failure or short duration transport of the device.




Fig. 12 Functional diagram of a switched-mode power supply unit

3.7 Evita 4 Graphics Controller PCB

The Graphics Controller PCB is located in the Evita 4 control unit. The board includes an independent processor system, the voltage supply of the control unit, the interface to the electronics and the display control.




Fig. 13 Graphics Controller PCB block diagram

3.7.1 Processor system

The processor system comprises a 68332 CPU, a 256 kByte RAM and a 2 MByte flash EPROM (electrically programmable and erasable read-only memory).

3.7.2 Uprog

Uprog generates the voltage required for programming the flash EPROMs. An enable logic prevents unintentional supply of the programming voltage to the flash EPROMs.

3.7.3 Reset

The Graphics Controller PCB operates independently and is fitted with its own reset generator. The display can be reset by the CPU 68332 PCB.

3.7.4 Loudspeaker monitoring

The loudspeaker monitoring checks the loudspeaker function. The current passing through the loudspeaker is measured by a resistor and then analyzed.

3.7.5 Keyboard port and LED driver

The keyboard and the LEDs are arranged in a **3 x 8** matrix array. The LEDs are controlled by transistors. The keyboard columns are also selected in the same way. A driver module controls the selection of lines of the LEDs in the matrix. The 68332 CPU scans the keyboard using Schmitt triggers.

3.7.6 Display control

The video processor TMS 34010 controls the image processing, generates the synchronized signals for the display and displays the image to be represented by means of the display logic. The video RAM is the memory of the video processor. The size of the video RAM is 512 kilobytes. The video RAM stores the video processor program, the video processor data and the displayed image.

3.7.7 CAN interface

The CAN interface is a fast, serial interface (Controller Area Network). The control unit, the electronics and the pneumatics communicate via a CAN interface. The transmission rate is 1 Mbit/s.

3.7.8 DC/DC converter

The DC/DC converter supplies +5 V and +12 V for the voltage supply of the control unit. The input voltage of the DC/DC converter is +24 V. Since the +5 V and +12 V voltages are generated directly on the Graphics Controller PCB there is a tolerance of $\pm 2.5\%$. The maximum power consumption of the control unit is 1.5 A for +5 V and 2.2 A for +12 V.

3.8 Communication PCB (option)

The Communication PCB is integrated in the electronic assembly of the device. The PCB includes an independent processor system, the voltage supply of the interfaces, an internal CAN interface, an external CAN interface, two RS232 interfaces and two analog outputs.




Fig. 14 Block diagram of the Communication PCB

3.8.1 Processor system

The processor system comprises a 68332 CPU, a 512 kByte RAM and a 1 MByte flash EPROM (electrically programmable and erasable read-only memory).

3.8.2 Uprog

Uprog generates the voltage required for programming the flash EPROMs. An enable logic prevents unintentional supply of the programming voltage to the flash EPROMs.

3.8.3 Reset

The reset logic generates a defined reset after power-up. The CPU 68332 can be reset by the CPU 68332 PCB.

3.8.4 CAN interfaces

The Communication PCB is not connected to the data bus of the CPU 68332 PCB. The data are transmitted via an internal CAN interface (Controller Area Network – fast, serial interface). The external CAN interface is electrically isolated from the device. Electrical isolation is made by means of optocouplers.

3.8.5 RS232 interfaces

The Communication PCB provides an RS232 interface in the device. The interfaces are labeled COM2 and COM3. The interfaces are electrically isolated from the device. Electrical isolation is made by means of optocouplers.

3.8.6 Analog outputs

The analog outputs supply voltages between 0 V and 4.095 V. The assignment of analog outputs is freely selectable. The resolution of the output voltage is 1 mV per bit.

3.8.7 DC/DC converter

The DC/DC converters supply +5 V ISO each for the voltage supply of the interfaces. The input voltages of the DC/DC converters are +5 V.

3.8.8 EEPROM

The EEPROM stores internal data of the interface. The EEPROM has a 2 kB capacity.

3.9 Paediatric Flow PCB (optional)

The Paediatric Flow PCB is integrated in the electronic unit of the device. The PCB has two flow measuring channels for connection of the Babylog flow sensor, one four-channel multiplexer, one twelve-bit A/D converter and one interface to the CPU 68332 PCB.




Fig. 15 Block diagram of the Paediatric Flow PCB

3.9.1 Flow measurement

The flow is measured with the Babylog flow sensor. The Babylog flow sensor has two measuring wires. One measuring wire is "covered" by a plastic bar; the Babylog flow sensor recognizes the direction of the flow. The flow is analyzed by a measuring bridge.

3.9.2 Multiplexer

The multiplexer consists of four analog-value selectors. The analog-value selectors are controlled by the software. The multiplexer transmits the flow sensor measurement signals, CGND and UREF LIMIT via a buffer (AMP) to an A/D converter (ADC).

3.9.3 A/D Converter

The input voltage of the A/D converter ranges from 0 V to 10 V. The A/D converter converts the flow measurement values into digital data. The CPU PCB controls the A/D converter and the multiplexer via an interface (DATA OUT; DATA IN). The voltage drop across the multiplexer, the buffer and the A/D converter is measured using the UREF LIMIT reference voltage and can be taken into account when measuring the flow.

3.10 IFCO Carrier PCB (optional nurse call and other)

3.10.1 General information

The IFCO Carrier PCB functions as carrier board for other optional features; it supports the nurse call and Remote Pad (cable remote control) functions and one additional ambient pressure sensor (optional).




Fig. 16 IFCO Carrier PCB

3.10.2 Nurse call

The nurse call transmits top-priority alarms (!!!) displayed on the screen to, e.g., a central station. An alarm is also transmitted if the internal loudspeaker for audible alarms fails. The alarm is reset automatically as soon as the cause of the alarm disappears. The alarms are suppressed during the boot phase ("start phase" of the device).

The central station alarm signaling is carried out by relay contacts. The alarm status can be determined by scanning these relay contacts. An alarm is considered to be triggered if contacts 3 and 5 of the "nurse call" connectors are closed by the relay contacts. These relay contacts are electrically isolated from the rest of the electronics. A bistable relay is used in order to keep this alarm signaling function active even when the power supply fails.

The nurse call function monitors the contacts of the power switch thus being able to detect whether the device is switched on or off. A power failure in the device can thus be detected. This alarm (power failure) can be cancelled by switching the device off. However, the device must be switched off within a defined time window. This time window is dependent on the charge status of specific capacitors. The time window is typically 2 minutes.

3.10.3 Remote Pad

The Remote Pad is a cable remote control that is connected to the IFCO Carrier PCB. The Remote Pad has 6 keys to control the ventilator. The Remote Pad is also provided with alarm LEDs. The inputs and outputs of the Remote Pads are short-circuit protected and protected against inadvertent wrong connections.

3.10.4 Ambient pressure sensor (optional)

The measured value supplied by the ambient pressure sensor on the IFCO Carrier PCB is evaluated by the safety software. This measured value is used to monitor the ambient pressure sensor on the CO₂ Carrier PCB.

3.11 Pneumatic Controller PCB

The Pneumatic Controller PCB is located in the pneumatic assembly of the device. The PCB provides the following functions:

- Supply pressure measurement.
- Inspiratory and expiratory airway pressure measurement.
- Esophagus pressure measurement.
- Flow measurement.
- Measurement of the O₂ concentration in the respiratory gas.
- Fan monitoring.
- Solenoid valve control.
- PEEP valve control.
- Interface to the HPSV Controller PCBs.
- CAN interface.
- 68332 CPU with RAM and flash memory.
- Serial EEPROM for storage of device configuration.

3.11.1 Pneumatic Controller PCB flow measurement




Fig. 17 Pneumatic Controller PCB flow measurement

3.11.2 Pneumatic Controller PCB block diagram




Fig. 18 Pneumatic Controller PCB block diagram

3.11.3 Pressure measurement

The pressure is measured with two pressure sensors. During the inspiratory phase and the expiratory phase, the pressure is measured by the respective pressure sensor. The airway pressure signals are transmitted to the multiplexer. The pressure sensors are automatically zero-calibrated every 3 minutes.

3.11.4 Flow measurement

The sensor works according to the principle of a constant-temperature hot-wire flowmeter. Respiratory gas flows along a very thin, electrically heated platinum wire in a measuring tube. The wire is heated to a temperature of 180 °C and kept at this temperature with a closed-loop control system. If gas flows past this wire, heat is dissipated. The larger the gas volume per time unit flowing past, the higher the heat dissipation. The heating power required to keep the wire at a constant temperature is a measure of the gas flow.

3.11.4.1 Cleaning the flow sensor

The actuator to clean/anneal the sensor applies a defined current to both measuring wires. During calibration the measuring wires begin to glow and burn any impurities. The sensor is cleaned when pressing the "flow calibration" key or automatically after drug nebulization. The cleaning process occurs during the next inspiratory phase or after 15 seconds at the latest.

3.11.5 O₂ measurement

The O₂ sensor works according to the galvanic cell principle. Oxygen molecules contained in the gas mixture to be measured diffuse through a plastic diaphragm into the electrochemical cell and are reduced at the noble metal electrodes. At the same time a base electrode is being oxidized. The base electrode is spent by the oxidation process and thus determines the life of a sensor. The current flowing through the cell is proportional to the oxygen partial pressure in the gas mixture to be measured.

Provided the pressure and temperature of the gas mixture to be measured are kept constant, the measured value will be directly proportional to the oxygen partial pressure. The O₂ amplifier on the O₂ Top PCB is mounted externally on the inspiratory block. The output signal is transmitted to the O₂ Contact PCB via spring contacts and transferred to the Pneumatic Controller PCB. The O₂ cell is also connected to the O₂ Top PCB via spring contacts.

3.11.6 Fan monitoring

At the front panel of the device a fan is mounted to limit the temperature and the O₂ concentration in the electronic unit of the pneumatics in case of failure; the fan is monitored by the electronics.

3.11.7 Multiplexer

The multiplexer consists of 16 analog-value selectors. The analog-value selectors are controlled by the software. The multiplexer routes the measurement signals from the pressure sensors, the O₂ amplifier, the FAN UREF and the flow sensor to an A/D converter via a buffer.

3.11.8 Solenoid valves

The solenoid valves are controlled by two power drivers. The power driver outputs for the nebulizer and the O₂/Air switchover are monitored by a comparator. The power drivers can be switched off by the electronic unit.

3.11.9 PEEP valve control

The PEEP valve is controlled by a voltage-controlled current source with power MOS FET. A quad operational amplifier serves to adapt the D/A converter output signal to the current range of the PEEP valve. The CPU controls the D/A converter. The PEEP valve control is calibrated to the PEEP valve; the calibration data are stored in the serial EEPROM.

3.11.10 HPSV interface

The status lines of the HPSV Controller PCB are led to the Pneumatic Controller PCB via the pneumatics motherboard. Two bus drivers transmit the data to the data bus of the CPU. The data are transmitted to the HPSV Controller PCB by two power drivers. The data are accepted by power swing of the respective chip-select pin.

3.11.11 CAN interface

The CAN interface comprises a CAN controller and a series-connected driver. The CAN controller is directly connected to the data bus of the CPU. The control unit, the electronics and the pneumatics communicate via a CAN interface. The transmission rate is 1 Mbit/s.

3.11.12 Processor system

The processor system on the Pneumatic Controller PCB consists of a 68332 CPU, a 256 kByte flash EPROM (electrically programmable and erasable read-only memory) and a 256 kByte RAM.

3.11.13 Serial EEPROM

The serial EEPROM stores the data of the pneumatics. The EEPROM has a 128 Byte capacity.

3.12 HPSV Controller PCB

The pneumatics of the device contain two identical HPSV Controller PCBs. The board slot determines which of the boards is assigned to O₂ and AIR. The PCB comprises the following functions:

- Microcontroller with EPROM and RAM.
- ADC for measurement of supply pressure.
- DAC for default setpoint of current.
- Control loop for current control.
- Power transistor (power source).

NOTE

The characteristic of the HPSV cartridge is stored in the cartridge and is read out by the HPSV Controller PCB.




Fig. 19 HPSV Controller PCB with HPSV

3.13 Pneumatic assembly

Air and O₂ must be available at a supply pressure of 2.7 to 6 bar to drive the device.

The pneumatic assembly consists of the following subassemblies:

- Gas connection block
- Parallel mixer or mixer block
- Pressure sensors
- PEEP/PIP valve
- Inspiratory block
- Patient system




Fig. 20 Layout of the pneumatic subassemblies

Item	Designation
1	Gas connection block
2	PEEP/PIP valve
3	Parallel mixer
4	Pressure sensors
5	Patient system (expiratory)
6	Inspiratory block

3.13.1

Gas connection block

The gas connection block comprises the **O₂** gas connection (M 12x1 female) and the **AIR** connection (M 20x1.5 male). In the future, the following gas connections will be available: NIST, DISS (USA) and DIN. The connections are fitted with filters **F1.1** and **F1.2** (metal fiber web). The non-return valves **D1.1** (AIR) and **D1.2** (O₂) prevent the gas from flowing back into the central gas supply system.

The pressure regulators **DR1.1** and **DR1.2** are set to 2 bar. The control gas flows past the **DR1.1** to the 3/2-way valve **Y1.1**, from there to the emergency valve **Y1.3**, to the PEEP/PIP valve **Y4.1** and finally to the emergency valve **Y3.1**.

The gas also flows to the expiratory pressure sensor S6.2 (purge flow) via the restrictor **R1.1** (0.08 L/min).

Gas flows to the nebulizer via the 3/2-way valve **Y1.4**, if appropriately adjusted.

In the event of AIR supply failure, the device will switch over to O₂ supply ("switchover function").




Fig. 21 Gas connections function diagram

3.13.1.1

Use of different gas connection blocks

Two different gas connection blocks are used:

- Dräger gas connection block
- FAS gas connection block




Fig. 22 Dräger gas connection block (1) and FAS gas connection block (2)

Item	Designation
1	Dräger gas connection block
2	FAS gas connection block
DR1.1	AIR pressure regulator
DR1.2	O2 pressure regulator
Y1.1	3/2-way solenoid valve, O2/AIR
Y1.2	3/2-way solenoid valve, O2 sensor calibration
Y1.3	3/2-way solenoid valve, venting
Y1.4	3/2-way solenoid valve, nebulizer

3.13.2 Characteristics of various gas connection blocks

The gas connection blocks can be differentiated by the following items:

- If the gas connection plate on the back of the device is secured with four screws, then it is a Dräger gas connection block, see also Fig. 23/1. With the FAS gas connection block, the gas connection plate is secured with two screws only - see also Fig. 23/2.
- If the pressure regulator housing on the gas connection block is made of brass, then it is a Dräger gas connection block. The pressure regulator housings of the FAS gas connection block are made of black plastic.




Fig. 23 Connection side of the gas connection blocks

Item	Designation
1	Dräger gas connection block
2	FAS gas connection block
3	Connector plate

The gas connection blocks are interchangeable since their mounts for attaching to the device are identical.

3.13.3

Parallel mixer or mixer block

The parallel mixer is a fast, electrically controllable proportional valve for gas flows between 5 and 180 L/min at supply pressures of 3 to 6 bar. Partial flows of less than 5 L/min are pulsed at a constant flow of 5 L/min. The supply gases compressed air (AIR) and oxygen (O₂) available at the parallel mixer have a supply pressure of 2.7 bar to 6 bar. In the parallel mixer the two gases are mixed in accordance with the set parameters. The parallel mixer supplies the inspiratory gas to the patient.

The parallel mixer consists of the following components:

- Mixer connection block
- Two cartridge valves with displacement sensor system for compressed air (AIR) and oxygen (O₂)
- 2 supply pressure sensors measuring the inlet pressure of the supply gases

The HPSV Controller PCB controls the parallel mixer electrically. The control signals are transmitted to the parallel mixer via the Pneumatic Motherboard PCB and the Pneumatic Controller PCB.




Fig. 24 Parallel mixer


Item	Designation
1	AIR cartridge valve
2	O ₂ cartridge valve
3	Supply pressure sensor
4	Mixer connection block

3.13.4 Mixer connection block

The two cartridge valves are mounted to the mixer connection block. The inspiratory gases in the mixer connection block are supplied to the respective cartridge valve. The respiratory gas available at the outlet of the cartridge valves is mixed in the mixer connection block and supplied to the inspiratory unit.

3.13.5 Cartridge valves (O₂ and AIR HPS valves)

The cartridge valve or HPS valve (HPSV = high-pressure servo-valve) supplies a defined amount of gas to the patient in accordance with the preset adjustment parameters for e.g. inspiration, trigger pressure, leak flow compensation.



4480

Fig. 25 Cartridge valves or O₂, AIR HPS valves

Item	Designation
1	O-rings
2	Input
3	Outlet
4	Adhesive label: O ₂ or AIR

During expiration the supply gas is available at the cartridge valve and at the supply pressure sensor. In the HPS valve ball A is pressed into the valve seat B; this action closes the valve.

In the event of a flow demand, e.g. inspiration, current is supplied to the drive system. The drive system is equipped with a coil working according to the principle of a moving coil of the type used e.g. in loudspeakers. The plunger is deflected in proportion to the supplied current. The ball A is pushed out of its valve seat B, causing an annular gap. The cartridge valve opens and gas flows to the patient.




Fig. 26 HPS valve cross-section

Item	Designation
1	Ball A
2	Displacement measuring system
3	Drive system
4	Displacement Sensor PCB
5	Plunger
6	Valve seat B

Consequently, the size of the annular gap determines the gas flow. The annular gap between the ball and the valve seat is determined by the deflection of the plunger. This deflection is measured and controlled by the displacement measuring system. At the same time a supply pressure sensor at the cartridge valve measures the supply pressure of the gas to be dosed. The supply pressure, along with the gap width, determines the gas flow. Twice the pressure (absolute) for the same gap width will produce twice the flow.

Gap width (annular gap) - or supply pressure - = gas flow -.

The gas flow in valves used to meter and mix the respiratory gas is not measured but determined indirectly by the displacement signal and the supply pressure. The valves supply a flow of 5 L/min to 180 L/min. The plunger movement (displacement) depends on the supply pressure. The supply pressure is between 3 bar and 6 bar absolute. For a flow requirement of 120 L/min and a supply pressure of 5 bar the displacement will be approx. 0.2 mm or 0.6 mm for a supply pressure of 1.5 bar.

As the mixer must permanently operate with high precision it is measured and the data (non-linearity) are stored in an EEPROM on the Displacement Sensor PCB. Therefore no calibration is required when the cartridge valves are replaced.

3.13.5.1 Displacement sensor system

The displacement sensor system measuring the position of the plunger in the valve is integrated in the cartridge valve. The displacement sensor system consists of the displacement measuring system and the Displacement Sensor PCB.

The displacement measuring system is designed as a differential transformer. The AC voltage applied to the primary winding of the transformer has a frequency of approx. 1 MHz. The two secondary windings are switched such that their output voltages balance out. If the ferrite core (plunger of the cartridge valves) moves in the differential transformer, the output voltage of the displacement sensor system will change.

As the displacement output signal is not linear to the gas flow, the characteristic of the cartridge valve is measured and stored in the EEPROM. The micro-controller on the HPSV PCB thus balances the non-linearity of the cartridge valve.

The two circuits of the cartridge valves of parallel mixers are operated asynchronously in parallel (AIR, O₂). To avoid beat interference, the frequencies of the two oscillators must differ by a minimum of 200 kHz. To do so, the two cartridge valves are measured at two different frequencies. The cartridge valves can only be measured at a special test stand.

3.13.5.2 Supply pressure sensors

The supply pressure sensors are calibrated to absolute pressure (0 bar). They measure the inlet pressure of the supplied gas. The high-pressure sensor is fitted with a P/U converter supplying a pressure-dependent output pressure.

- Measuring range: 0 bar to 7 bar.
- Sensitivity: 1.58 V/bar \pm 8 mV/bar.
- Offset voltage: 300 mV \pm 30 mV

The high-pressure sensors are linked via a flex-strip to the Displacement Sensor PCB. The Displacement Sensor PCB is installed in the HPS valves.

3.13.6 Pressure sensors

The pressure sensor mount comprises the airway pressure sensors **S6.1** for the inspiratory side and **S6.2** for the expiratory side. **S6.1** monitors the inspiratory Paw high and Paw low.

- Measuring range: 140 mbar.
- Sensitivity: 36.5 mV/mbar \pm 0.3 mV/mbar.
- Offset voltage: 1.74 V \pm 0.04 V.

3.13.6.1 Pressure sensor calibration

The solenoid valves "Y6.1" and "Y6.2" expose the relevant pressure sensors to atmospheric pressure at specific time intervals, consequently, the sensors are automatically calibrated. The pressure sensors "S6.1" and "S6.2" are zero-calibrated every 3 minutes. (Calibrate the electronic zero drift as necessary). To do so, the solenoid valves "Y6.1" and "Y6.2" are subsequently exposed to atmospheric pressure and the pressure sensors automatically calibrated.




Fig. 27 Pressure sensors function diagram

3.13.7

Inspiratory block

The safety valve "D3.3" limits the pressure in the inspiratory line to 100 mbar max.

In the event of compressed air failure or power failure the pneumatically controlled emergency air valve "Y3.1" will open so that the patient can breathe ambient air passing the filter "F3.1". The non-return valve "D3.1" prevents rebreathing of the air through the inspiratory line. When the valve "Y3.1" opens, the spring-loaded non-return valve "D3.2" allows the pressure to be relieved.

In the case of emergency air spontaneous breathing the patient can exhale through the expiratory valve "Y5.1" on account of the spring loading (5 mbar), thus preventing rebreathing.

The inspiration block is provided with the plug-in connection for the oxygen sensor.

The restrictor "R1.2" limits the drug nebulizer flow to 9 L/min.




Fig. 28 Emergency air function diagram

Item	Designation
1	Nebulizer
2	Inspiration




Fig. 29 Inspiratory block

Item	Designation	Item	Designation
1-1a	Patient connection (inspiration)	9	O2 calibration control
2-2a	Nebulizer port	10	O2 sensor chamber (behind the O2 amplifier)
3-3a	Esophagus pressure P1	11	O2 calibration purge flow outlet
4-4a	Esophagus pressure P2	12	Emergency air non-return valve D3.1 in O2 amplifier
5	10 mbar pressure relief valve D3.2	13	O2 calibration diaphragm lattice Y3.3 with R3.1
6	100 mbar pressure relief valve D3.3 (adjustable; Note: until middle of '96 valve slot milled in the block, after '96 separate piece screwed onto the block)	14	Valve Y3.1
7	Emergency venting control	15	O2 amplifier
8	Pressure measurement (inspiration)	-	

3.13.8 Patient system

The expiratory gas flows from the patient directly to the expiratory valve **Y5.1**. The copper measuring line at the **8a** connection has a germicidal effect and connects the expiratory side to the pressure sensor **S6.2**.

The expiratory valve has a transmission ratio of approx. 1:1. The non-return valve **D5.1** allows flow in one direction only and makes sure that gases do not travel backwards. The expiratory flow is measured with flow sensor **S5.1**.




Fig. 30 Patient system function diagram

Item	Designation
1	Expiration

The ratio between the control pressure at the **7a** connection of the PEEP/PIP valve and the resulting pressure at the expiratory port is linear to the following values:

- Control pressure of 3 mbar => expiratory pressure of 0 mbar.
- Control pressure of 33 mbar => expiratory pressure of 33 mbar.




Fig. 31 Sectional view of the patient system

3.13.8.1

Patient system with RS water trap 84 13 125

The water trap avoids flow measurement faults caused by water droplets. Such faults may occur if the water traps on the ventilation tubes are not positioned at the lowest possible point. In this case the condensation water is collected in the water trap of the patient system.

The collector jar of the water trap can be removed during operation. The opening to the patient system is sealed automatically.




Fig. 32 Function diagram of the patient system with water trap

Item	Designation
1	Expiration
2	Water trap




Fig. 33 Sectional view of the patient system with water trap

Item	Designation
1	Water trap conversion kit
2	Angle socket
3	Locking mechanism
4	Connector

3.13.9

PEEP/PIP valve

The PEEP valve **Y4.1** consists of a diaphragm valve acting as a flow-control device and the linear drive whose plunger closes the diaphragm valve. A coil drives the PEEP valve **Y4.1**. The values are set via the ventilation settings. These settings are processed by a computer program and the coil is driven by an appropriate current. The PEEP valve opens and adjusts a pressure proportional to the adjusted electric current. 0 mA will correspond to -1 mbar, 500 mA to 120 mbar.

The valve **Y4.1** controls the expiratory valve **Y5.1** in the patient system via a servo-line. The solenoid valve **Y1.3** and the restrictor **R4.1** supplies the patient system with control gas. The non-adjustable restrictor **R4.1** is set to a flow of 3.5 L/min.

Depending on the setting the plunger coil is activated causing an appropriate servo-pressure to be applied to the expiratory valve.

The software compares the preset and measured airway pressures. This comparison is a measure of the Pneumatic Controller PCB's control action on the PEEP/PIP valve. The PEEP/PIP valve is calibrated to the electronics. The calibration data are stored on the Pneumatic Controller PCB.




Fig. 34 PEEP valve function diagram




Fig. 35 Section view of the PEEP valve

3.13.10 AIR supply

AIR flows through the filter **F1.1** via the non-return valve **D1.1** to the mixing and dosage unit (pressure sensor **S2.1** and HPSV **Y2.1**), and at the same time via the pressure regulator **DR1.1**, which is set to 2 bar, to the 3/2-way solenoid valve **Y1.1**. From here the gas flows through the 3/2-way solenoid valve **Y1.3** to the emergency air valve **Y3.1** which closes. Furthermore, AIR passes the restrictor **R4.1** to reach the PEEP/PIP valve **Y4.1** and from there – depending on the setting – to the expiratory valve **Y5.1**. Finally, AIR passes the restrictor **R1.1** to flow to the expiratory pressure sensor **S6.2** connecting line on the patient side. At this point, expiratory humidity is prevented from reaching the pressure sensor **S6.2**.




Fig. 36 AIR supply function diagram

Item	Designation
1	Nebulizer
2	Inspiration
3	Expiration

3.13.11 O₂ supply

Compressed oxygen flows through the filter **F1.2** via the non-return valve **D1.2** to the mixing and dosage unit (pressure sensor **S2.2** and HPSV **Y2.2**), at the same time via the **DR1.2**, which is set to a downstream pressure of 2 bar, to the 3/2-way solenoid valve **Y1.1**.




Fig. 37 O₂ supply function diagram

3.13.11.1 O2/AIR changeover valve (switchover)

In the event that compressed air fails, in the course of O₂ sensor calibration and when controlling the nebulizer function (depending on the FiO₂ setting) the 3/2-way solenoid valve **Y1.1** switches over. The servo-system will then be supplied with O₂.

3.13.12 Inspiration

Depending on the settings (O₂ concentration, inspiratory volume, frequency, T_i, inspiratory flow, inspiratory pressure) the HPSVs **Y2.1** and **Y2.2** open. The gas flows via the inspiratory connector to the patient. At the same time, gas flows to the O₂ sensor **S3.1** and to the safety valve **D3.3**; from there, it flows through the 3/2-way solenoid valve **Y6.1** to the inspiratory pressure sensor **S6.1**.

The safety valve **D3.3** is fixed to 100 mbar and serves as an additional safety device in the event of a complete failure of the electronic control.

When calibrating the O₂ sensor **S3.1** the O₂ sensor **S3.1** will be disconnected with valve **Y3.3** from the inspiratory gas. The O₂ sensor **S3.1** is flushed with calibration gas via valve **Y1.2**, restrictor **R1.3**, restrictor **R3.1** and valve **Y3.2**. The O₂ concentration and the inspiratory gas flow are not affected.

The pressure sensors **S6.1** and **S6.2** monitor the inspiratory pressure. During the entire inspiratory time a pressure is applied to the expiratory valve **Y5.1** via the PEEP/PIP valve **Y4.1**.




Fig. 38 Inspiration function diagram

Item	Designation
1	Nebulizer
2	Inspiration
3	Expiration

3.13.13 Paw high

If the set Paw high alarm limit is exceeded during inspiration, the HPSVs **Y2.1** and **Y2.2** interrupt the gas flow. The PEEP/PIP valve **Y4.1** is switched to expiration and the patient can exhale.

3.13.14 Emergency venting

If the Paw high alarm limit is exceeded by 5 mbar, the 3/2-way valve **Y1.3** additionally opens. This will cause the emergency air valve **Y3.1** to open and the pressure can be reduced via the non-return valve **D3.2**.

3.13.15 Expiration

Upon start of expiration (Fig. 39/1) the HPSVs **Y2.1** and **Y2.2** are closed. No gas flows to the patient. The PEEP/PIP valve **Y4.1** is switched to the set PEEP value. The expiratory valve **Y5.1** will also be relieved and the patient can exhale via non-return valve **D5.1** and the flow sensor **S5.1**. The flow sensor **S5.1** measures the expiratory volume.




Fig. 39 Expiration function diagram

3.13.16 Drug nebulizer

After pressing the "Drug nebulizer" key, the drug nebulizer is switched on for 30 minutes. At the same time the solenoid valve **Y1.4** is switched through in the flow-active inspiratory phase. The drug nebulizer is supplied with drive gas through the restrictor **R1.2**. After completion of the inspiratory gas supply phase, the solenoid valve **Y1.4** is also switched back. The minute volume remains constant while the flow setting is being corrected. When drug nebulizing is complete, the flow sensor **S5.1** is automatically annealed clean.

The drug nebulizer needs a minimum inspiratory flow of 16 L/min.




Fig. 40 Nebulizer function diagram

Item	Designation
1	Nebulizer
2	Inspiration

3.14 Gas mixture

3.14.1 Functional principle

The AIR flow and the O₂ flow are metered directly out of the supply lines with one HPSV each; both sub-flows flow together to the inspiration block in the so-called parallel mixer. Depending on the set FiO₂ (0.21 to 1.00), the total flow to be metered is split into an air sub-flow and an oxygen sub-flow. Sub-flows of less than 5 L/min are no longer metered continuously, but in pulses lasting at least 8 ms and with a constant flow value of 5 L/min. This results in a pulse/pause ratio that corresponds to the sub-flow value.

3.14.2 Correcting the oxygen HPSV

Corresponding to its operating principle, the HPSV meters a mass flow. Owing to the differing gas density values of the compressed air and the oxygen, different volume flows would therefore be applied if these values were not corrected.

Owing to the different gas density values, the portion of the inspiration flow demanded of the oxygen HPSV is therefore boosted by 5%.

3.14.3 Dependence on the supply pressure

With regard to the gas supply, the operating range of the device is specified from 2.7 bar to 6 bar overpressure. The device monitors this operating range with the aid of the absolute pressure sensors on the HPSVs by the following criteria:

- Supply pressure always higher than 1.2 bar absolute
- At no flow, supply pressure greater than 2.5 bar relative.

At low supply pressures (below 3 bar overpressure), the HPSV can no longer apply high flow demands without errors, i.e. the actually supplied flow is less than the demanded flow. The inspiratory flow must be limited for reasons of flow accuracy and thus also for oxygen concentration accuracy reasons.

Conversion inside device. Prior to the start of an inspiration, the maximum inspiration flow for spontaneous inspiration of the HPSV is defined on the basis of the measured supplied overpressure:

- P_{supply} > 4 bar absolute > max. inspiratory flow = 180 L/min.
- P_{supply} < 4 bar absolute > max. inspiratory flow = 150 L/min.

For mandatory inspirations, the inspiration flow is limited by the setting to 120 L/min.

3.14.4

BTPS/NTPD conversion

BTPS = Body Temperature, Pressure, Saturated = body temperature (37 °C), lung pressure, 100% humidity

NTPD = Normal Temperature, Pressure, Dry = normal temperature (20 °C), ambient pressure, 0% humidity

The HPSV is a mass source and therefore supplies the flow under NTPD conditions. Nevertheless, the user setting of the inspiration flow or the demand by the demand flow generator is based on BTPS conditions.

Prior to control of the HPSVs, the BTPS flow demands must be converted to an NTPD flow supply.

$$\text{Flow_NTPD} = \text{Flow_BTPS} * \frac{P_{\text{umgebung}} + P_{\text{Insp_End}} - 63 \text{ mbar}}{1013 \text{ mbar}} * \frac{293\text{K}}{310\text{K}}$$

Fig. 41 Conversion

Item	Designation
P_umgebung	Value measured by the ambient pressure sensor.
P_Insp_End	Measured airway pressure at the end of the last inspiration.

3.15 Pneumatics function diagram




Fig. 42 Pneumatics diagram

Item	Designation
1	Gas connection.
2	Parallel mixer.
3	Inspiration block.
4	Nebulizer.
5	Inspiration.
6	Pressure measurement
7	Expiration.
8	Patient system.
9	PEEP/PIP valve.
AIR	Gas connection, compressed air.
D1.1	Non-return valve.
D1.2	Non-return valve.
D3.1	Non-return valve.
D3.2	Non-return valve 10 mbar.
D3.3	Non-return valve 100 mbar.
D5.1	Non-return valve.
DR1.1	AIR pressure regulator.
DR1.2	O ₂ pressure regulator.
F1.1	Filter.

Item	Designation
F1.2	Filter.
F3.2	Filter.
O2	Gas connection, compressed oxygen.
R1.1	Restrictor 0.08 L/min/2 bar.
R1.2	Restrictor 9 L/min/2 bar.
R1.3	Restrictor 0.4 L/min/2 bar.
R3.1	Restrictor (hole in the diaphragm in Y3.3) 0.25 L/min/1.4 bar.
R4.1	Restrictor 3.5 L/min/2 bar.
S2.1	AIR pressure sensor (HPSV).
S2.2	O2 pressure sensor (HPSV).
S3.1	O2 sensor.
S5.1	Flow sensor.
S6.1	Inspiratory pressure sensor.
S6.2	Expiratory pressure sensor.
Y1.1	3/2-way solenoid valve, O2/AIR.
Y1.2	3/2-way solenoid valve calibration O2 sensor.
Y1.3	3/2-way solenoid valve, venting.
Y1.4	3/2-way solenoid valve, nebulizer.
Y2.1	HPSV AIR (high-pressure servo-valve) parallel mixer.
Y2.2	HPSV O2 (high-pressure servo-valve) parallel mixer.
Y3.1	Emergency air valve.
Y3.3	Inspiratory valve.
Y4.1	PEEP/PIP valve.
Y5.1	Expiratory valve.
Y6.1	3/2-way solenoid valve, inspiration.
Y6.2	3/2-way solenoid valve, expiration.

4 Parts catalog

Parts catalog

This chapter contains a list of the device's orderable parts.



Parts catalog

Evita 4/ Evita 4 edition

Revision: 22
5664.500

Parts catalog
Evita 4/ Evita 4 edition


Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
	MX44200	<input type="checkbox"/>	Tools		St	
	MX44201	<input type="checkbox"/>	Products concerned		St	
	MX44068	<input type="checkbox"/>	Instructions for use		St	
	MX44205	<input type="checkbox"/>	Basic Unit		St	
	MX44203	<input type="checkbox"/>	Modification Kits/Options		St	
	MX44212	<input type="checkbox"/>	Maintenance Parts/Service Kits		St	
	MX44204	<input type="checkbox"/>	Accessories/Consumables		St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

2/107

Parts catalog
Tools




Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
	MX44355	<input type="checkbox"/>	Test equipment,non-calibrated	1.000	St	
	MX44356	<input type="checkbox"/>	Test equipment, calibrated	1.000	St	
	MX44357	<input type="checkbox"/>	Service tools	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

3/107

Parts catalog
Test equipment, non-calibrated




Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
	8290271	<input checked="" type="checkbox"/>	Testing gas for CO2 Monitor	1.000	St	
	7910364	<input checked="" type="checkbox"/>	Leakage curr. test lead, 2-pin, temperat	1.000	St	
	7901204	<input checked="" type="checkbox"/>	Wrench, special open-ended, WAF 17	1.000	St	
	7901607	<input checked="" type="checkbox"/>	27MM Spanner	1.000	St	
	7901541	<input checked="" type="checkbox"/>	Syring 50ml, empty	1.000	St	
	8400904	<input checked="" type="checkbox"/>	Taper adapter	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

4/107

Parts catalog
Test equipment, calibrated




Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
	7910724	<input checked="" type="checkbox"/>	Pres. meas. device, dig. 20bar	1.000	St	
	7910722	<input checked="" type="checkbox"/>	Pres. meas. device, dig. 1bar	1.000	St	
	7901161	<input checked="" type="checkbox"/>	Flowm., bl. , 0.02 - 14 L/min.	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

5/107

Parts catalog
Service tools

Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
	8412710	<input checked="" type="checkbox"/>	Calibration set CO2	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

6/107

Parts catalog
Products concerned




Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
	8411740	<input type="checkbox"/>	Evita 4, CE	1.000	St	
	8411900	<input type="checkbox"/>	Evita 4, Construction Kit	1.000	St	
	8412800	<input type="checkbox"/>	Evita 4 international	1.000	St	
	8412970	<input type="checkbox"/>	Evita 4, F	1.000	St	
	8412980	<input type="checkbox"/>	Evita 4, USA	1.000	St	
	8412990	<input type="checkbox"/>	Evita 4, J	1.000	St	
	8413454	<input type="checkbox"/>	ICU respirator Evita 4	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

7/107

Parts catalog
Instructions for use


Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
----------	-----------	----------	-------------	------	----------	--------

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

8/107

Parts catalog
Basic Unit




Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
	MX44207	<input type="checkbox"/>	Electronic	1.000	St	Servicedrawing
	MX44206	<input type="checkbox"/>	Pneumatic	1.000	St	Servicedrawing
	MX44208	<input type="checkbox"/>	Operating Device	1.000	St	Servicedrawing
1	8412380	<input checked="" type="checkbox"/>	Stripe	1.000	St	
2	8412858	<input checked="" type="checkbox"/>	Vent plug	2.000	St	
3	2M17288	<input checked="" type="checkbox"/>	Screw	8.000	St	
4	8412857	<input checked="" type="checkbox"/>	Vent plug	2.000	St	
5	8412384	<input checked="" type="checkbox"/>	Filter, single	2.000	St	
6	8406763	<input checked="" type="checkbox"/>	Bar	2.000	St	
7	G60455	<input checked="" type="checkbox"/>	Vent plug	4.000	St	
8	8406504	<input checked="" type="checkbox"/>	Socket	3.000	St	
9	1824481	<input checked="" type="checkbox"/>	Power cable CE,3m,10A,C13L,bk	1.000	St	
10	8411784	<input checked="" type="checkbox"/>	Parking support	1.000	St	
11	8412375	<input checked="" type="checkbox"/>	Hood	1.000	St	
12	8419618	<input checked="" type="checkbox"/>	Plastic screw	1.000	St	
13	1839780	<input checked="" type="checkbox"/>	Connecting cable	1.000	St	
14	8410575	<input checked="" type="checkbox"/>	Angular porcelain Bush	1.000	St	
15	8410580	<input checked="" type="checkbox"/>	Expiration valve, Patientsys.	1.000	St	
16	8416980	<input checked="" type="checkbox"/>	Op.Device Evita 4 edition	1.000	St	
17	8414828	<input checked="" type="checkbox"/>	Tray	1.000	St	
18	8412379	<input checked="" type="checkbox"/>	Hood 4H,1B	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

9/107

Parts catalog Pneumatic



Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
1	8412361	<input checked="" type="checkbox"/>	Socket	1.000	St	
2	8712065	<input checked="" type="checkbox"/>	Cable clip	7.000	St	
3	8412339	<input checked="" type="checkbox"/>	Chassis pneumaticPNEUMATIC	1.000	St	
4	8306545	<input checked="" type="checkbox"/>	Cable harness	2.000	St	
5	8306225	<input checked="" type="checkbox"/>	Screw	2.000	St	
6	8306560	<input checked="" type="checkbox"/>	Fan pneumatic, compl.	1.000	St	
7	MX44218	<input type="checkbox"/>	Inspiration Block	1.000	St	
8	MX44216	<input type="checkbox"/>	Housing	1.000	St	
9	MX44215	<input type="checkbox"/>	Parallel Mixer 4	1.000	St	
10	MX44225	<input type="checkbox"/>	Patient System	1.000	St	
11	8306601	<input type="checkbox"/>	PCB pneumatic controller	1.000	St	sale
11	8420855	<input checked="" type="checkbox"/>	Rep-Kit LP Pneumatikcont 42dXL	1.000	St	
12	8306561	<input checked="" type="checkbox"/>	PCB HPSV-controller	2.000	St	for repair only one piece necessary
13	MX44214	<input type="checkbox"/>	Pressure Metering Block	1.000	St	
14	8410717	<input checked="" type="checkbox"/>	Peep Valve	1.000	St	
15	2M04819	<input checked="" type="checkbox"/>	Rubber buffer	2.000	St	
35	8412718	<input checked="" type="checkbox"/>	Cover plate, backpanel	1.000	St	
36	8306571	<input checked="" type="checkbox"/>	PCB pneumatic motherboard	1.000	St	


Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

10/107

Parts catalog

Pneumatic



Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
16	MX44217	<input type="checkbox"/>	Patientsystem heater	1.000	St	
17	8412605	<input checked="" type="checkbox"/>	Plate	1.000	St	
18	8412363	<input checked="" type="checkbox"/>	Latch	2.000	St	
19	8408902	<input checked="" type="checkbox"/>	Pin	1.000	St	
20	2M17303	<input checked="" type="checkbox"/>	Catch	2.000	St	
21	8412354	<input checked="" type="checkbox"/>	Side part 2H	1.000	St	
22	8415709	<input checked="" type="checkbox"/>	Cable harness spiolog sensor	1.000	St	
23	8410663	<input checked="" type="checkbox"/>	Plug accommodation	1.000	St	
24	8410570	<input type="checkbox"/>	Connecting housing	1.000	St	
25	1190520	<input checked="" type="checkbox"/>	Hose 4X1,5-SI 50 SH A NF	1.000	m	
26	8410631	<input checked="" type="checkbox"/>	Spring	1.000	St	
27	8306556	<input checked="" type="checkbox"/>	Cable harness flowswitcher	1.000	St	
28	8410917	<input checked="" type="checkbox"/>	Rubber element	1.000	St	
29	8407689	<input checked="" type="checkbox"/>	Lip seal	2.000	St	
30	8412707	<input checked="" type="checkbox"/>	Hinge,cpl.	2.000	St	
37	8412382	<input checked="" type="checkbox"/>	Sheet metal,cpl	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

11/107



Parts catalog
Pneumatic

Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
30	8412707	<input checked="" type="checkbox"/>	Hinge,cpl.	2.000	St	
31	2M17302	<input checked="" type="checkbox"/>	Base	2.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

12/107

Parts catalog
Pneumatic

Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
1	8412361	<input checked="" type="checkbox"/>	Socket	1.000	St	
33	8414466	<input checked="" type="checkbox"/>	90° connector Air	1.000	St	
34	8413641	<input checked="" type="checkbox"/>	90 degrees connector O2	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

13/107

Parts catalog
Pneumatic


Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
	MX44221	<input type="checkbox"/>	Tubing Plan GS FAS	1.000	St	
	MX44223	<input type="checkbox"/>	Gas Supply FAS	1.000	St	
	1301489	<input checked="" type="checkbox"/>	Lock washer 6 DIN 6799	1.000	St	
	MX44222	<input type="checkbox"/>	Gas Supply Draeger	1.000	St	
	MX44220	<input type="checkbox"/>	Tubing Plan GS Draeger	1.000	St	
32	6800478	<input checked="" type="checkbox"/>	Buffer	4.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

14/107

Parts catalog
Tubing Plan GS FAS




Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
1	8401083	<input checked="" type="checkbox"/>	T-piece	4.000	St	
2	8412998	<input checked="" type="checkbox"/>	Dosage	1.000	St	
3	8412201	<input checked="" type="checkbox"/>	Dosage	1.000	St	
4	2M17645	<input checked="" type="checkbox"/>	Tube	1.000	St	
5	1190520	<input checked="" type="checkbox"/>	Hose 4X1,5-SI 50 SH A NF	1.000	m	
6	6800187	<input checked="" type="checkbox"/>	T-Socket	1.000	St	
7	1198343	<input checked="" type="checkbox"/>	Hose 7X2,5 SI NF M29908	1.000	m	
8	8412892	<input checked="" type="checkbox"/>	Dosage	1.000	St	
9	8414789	<input checked="" type="checkbox"/>	Dosage	1.000	St	
10	1203622	<input checked="" type="checkbox"/>	Hose 2X1,5 SI NF 8403323	1.000	m	
11	1204807	<input checked="" type="checkbox"/>	Hose 2X1-SIGN green	1.000	m	
12	1204823	<input checked="" type="checkbox"/>	Hose 2X1-sign orange	1.000	m	
13	1204793	<input checked="" type="checkbox"/>	Hose 2X1-SIGN red	1.000	m	
14	1204815	<input checked="" type="checkbox"/>	Hose 2X1-sign brown	1.000	m	
15	1204785	<input checked="" type="checkbox"/>	Hose 2X1-SIGN BL	1.000	m	
16	1210203	<input checked="" type="checkbox"/>	Hose 4X1PAE SW/WS M31989	1.000	m	
17	1210165	<input checked="" type="checkbox"/>	Hose 4X1 PAE WS	1.000	m	
18	8419520	<input checked="" type="checkbox"/>	Kit Volume for PEEP Valve	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

Parts catalog
Gas Supply Draeger




Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
1	8412209	<input checked="" type="checkbox"/>	Micro-electrovalve	4.000	St	
2	8408197	<input checked="" type="checkbox"/>	Socket	3.000	St	
3	M31601	<input checked="" type="checkbox"/>	Thrust collar 6 black	1.000	St	
4	M31603	<input checked="" type="checkbox"/>	Thrust collar 6, white	1.000	St	
5	M30960	<input checked="" type="checkbox"/>	Plug-type connection	2.000	St	
6	8303296	<input checked="" type="checkbox"/>	Adhesive label	1.000	St	
7	D21138	<input checked="" type="checkbox"/>	Check valve	2.000	St	
8	8408204	<input checked="" type="checkbox"/>	Gasket	1.000	St	
9	8408208	<input checked="" type="checkbox"/>	Filter (for 8411848)	2.000	St	
10	8412871	<input checked="" type="checkbox"/>	Receiver	1.000	St	
11	M19311	<input checked="" type="checkbox"/>	Washer	2.000	St	
12	M32489	<input checked="" type="checkbox"/>	Adapter O2	1.000	St	
13	8412872	<input checked="" type="checkbox"/>	Anschluss Air (Nist)	1.000	St	
14	8412873	<input checked="" type="checkbox"/>	Connection Air (DIN)	1.000	St	
15	8408205	<input checked="" type="checkbox"/>	Pressure regulator	2.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

Parts catalog
Tubing Plan GS Draeger




Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
1	8401083	<input checked="" type="checkbox"/>	T-piece	4.000	St	
2	8412998	<input checked="" type="checkbox"/>	Dosage	1.000	St	
3	8412201	<input checked="" type="checkbox"/>	Dosage	1.000	St	
4	2M17645	<input checked="" type="checkbox"/>	Tube	1.000	St	
5	1190520	<input checked="" type="checkbox"/>	Hose 4X1,5-SI 50 SH A NF	1.000	m	
6	6800187	<input checked="" type="checkbox"/>	T-Socket	1.000	St	
7	1198343	<input checked="" type="checkbox"/>	Hose 7X2,5 SI NF M29908	1.000	m	
8	8412892	<input checked="" type="checkbox"/>	Dosage	1.000	St	
9	8414789	<input checked="" type="checkbox"/>	Dosage	1.000	St	
10	1203622	<input checked="" type="checkbox"/>	Hose 2X1,5 SI NF 8403323	1.000	m	
11	1204807	<input checked="" type="checkbox"/>	Hose 2X1-SIGN green	1.000	m	
12	1204823	<input checked="" type="checkbox"/>	Hose 2X1-sign orange	1.000	m	
13	1204793	<input checked="" type="checkbox"/>	Hose 2X1-SIGN red	1.000	m	
14	1204815	<input checked="" type="checkbox"/>	Hose 2X1-sign brown	1.000	m	
15	1204785	<input checked="" type="checkbox"/>	Hose 2X1-SIGN BL	1.000	m	
16	1210203	<input checked="" type="checkbox"/>	Hose 4X1PAE SW/WS M31989	1.000	m	
17	1210165	<input checked="" type="checkbox"/>	Hose 4X1 PAE WS	1.000	m	
18	8419520	<input checked="" type="checkbox"/>	Kit Volume for PEEP Valve	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

17/107

Parts catalog
Gas Supply FAS




Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
1	8412209	<input checked="" type="checkbox"/>	Micro-electrovalve	1.000	St	Y1.1;Y1.2;Y1.3;Y1.4
2	8411515	<input checked="" type="checkbox"/>	Nonreturn valve (for 8411848)	2.000	St	
3	8411516	<input checked="" type="checkbox"/>	O-ring for 8411848	2.000	St	
4	8419505	<input checked="" type="checkbox"/>	FAS connection board, complete	1.000	St	
5	8408208	<input checked="" type="checkbox"/>	Filter (for 8411848)	2.000	St	
6	8411848	<input checked="" type="checkbox"/>	Gas Input Unit	1.000	St	
7	8411514	<input checked="" type="checkbox"/>	Valve tappet (for 8411848)	2.000	St	
8	8411513	<input checked="" type="checkbox"/>	Diaphragm FPM (for 8411848)	2.000	St	
9	8419255	<input checked="" type="checkbox"/>	FAS pressure regulator kit	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

18/107

Parts catalog
Inspiration Block




Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
1	8412955	<input type="checkbox"/>	Cover 1	1.000	St	
2	8414179	<input checked="" type="checkbox"/>	Diaphragm connection SV	1.000	St	
3	D13509	<input checked="" type="checkbox"/>	Closing spring	1.000	St	
4	R27696	<input checked="" type="checkbox"/>	O-ring 19x1	1.000	St	
5	8412951	<input checked="" type="checkbox"/>	Receiver	1.000	St	
6	8414597	<input checked="" type="checkbox"/>	O-ring seal 2,8x1,6	8.000	St	
7	8412981	<input checked="" type="checkbox"/>	Inspiration Unit,complete	1.000	St	
8	8410682	<input checked="" type="checkbox"/>	Diaphragm	1.000	St	
9	8411749	<input checked="" type="checkbox"/>	Cover, cpl. (Evita 4)	1.000	St	
10	8413031	<input checked="" type="checkbox"/>	Gasket	1.000	St	
11	8306671	<input checked="" type="checkbox"/>	Pba O2-contact	1.000	St	
12	8412938	<input checked="" type="checkbox"/>	Lip seal	1.000	St	
13	8600213	<input checked="" type="checkbox"/>	LUER-LOCK,female	2.000	St	
14	M20622	<input checked="" type="checkbox"/>	O-ring	1.000	St	
15	8410676	<input checked="" type="checkbox"/>	Angular porcelain Bush	1.000	St	
16	8414178	<input checked="" type="checkbox"/>	Diaphragm connection O2	1.000	St	
17	8412952	<input checked="" type="checkbox"/>	Screw	1.000	St	
18	2M12034	<input checked="" type="checkbox"/>	Spring	1.000	St	
19	8410307	<input checked="" type="checkbox"/>	Sealing washer	1.000	St	
20	2M08777	<input checked="" type="checkbox"/>	O-ring seal	1.000	St	
21	8411745	<input checked="" type="checkbox"/>	Pressure control valve	1.000	St	
22	R26807	<input checked="" type="checkbox"/>	O-ring, 18.77x1.78mm	1.000	St	
23	8412939	<input checked="" type="checkbox"/>	Socket	1.000	St	glue in with 7911820

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

19/107

Parts catalog
Inspiration Block




Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
24	7911820	<input checked="" type="checkbox"/>	UHU plus self-setting 300	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

Parts catalog
Housing



Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
1	8409177	✓	Sheet metal	1.000	St	
2	8407265	✓	Microswitch	1.000	St	Dreloba
3	8412356	✓	Housing, cpl	1.000	St	
4	8408197	✓	Socket	8.000	St	
5	E20567	✓	O-ring	1.000	St	
6	8407827	✓	Sintered filter	1.000	St	
7	2M08777	✓	O-ring seal	1.000	St	


Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

21/107

Parts catalog

Parallel Mixer 4




Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
1	8306225	<input checked="" type="checkbox"/>	Screw	4.000	St	
2	8412128	<input checked="" type="checkbox"/>	Valve Air	1.000	St	
3	R16442	<input checked="" type="checkbox"/>	O-ring	2.000	St	
4	8414375	<input checked="" type="checkbox"/>	O-ring seal	1.000	St	
5	E09548	<input checked="" type="checkbox"/>	O-ring 13.5x1.5	2.000	St	
6	8306248	<input checked="" type="checkbox"/>	Rep.set high pressure sensor	1.000	St	
7	8410713	<input checked="" type="checkbox"/>	O-ring seal	2.000	St	
8	M30961	<input checked="" type="checkbox"/>	Angle connection	2.000	St	
9	M31601	<input checked="" type="checkbox"/>	Thrust collar 6 black	1.000	St	
10	M31603	<input checked="" type="checkbox"/>	Thrust collar 6, white	1.000	St	
11	8412126	<input checked="" type="checkbox"/>	Valve O2	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

22/107

Parts catalog
Patient System



Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
1	8410576	<input checked="" type="checkbox"/>	Gasket	1.000	St	
2	8412202	<input checked="" type="checkbox"/>	Sieve	1.000	St	
3	8410682	<input checked="" type="checkbox"/>	Diaphragm	1.000	St	
4	8410580	<input checked="" type="checkbox"/>	Expiration valve, Patientsys.	1.000	St	Complete Unit
5	8410571	<input checked="" type="checkbox"/>	Valve housing	1.000	St	
6	M29902	<input checked="" type="checkbox"/>	O-Ring	2.000	St	without Watertrap
7	8410575	<input checked="" type="checkbox"/>	Angular porcelain Bush	1.000	St	
8	8413125	<input checked="" type="checkbox"/>	KIT Water trap, for exp. valve	1.000	St	
9	8403976	<input checked="" type="checkbox"/>	Pot	1.000	St	
10	8410572	<input checked="" type="checkbox"/>	Valve cover	1.000	St	
11	E22651	<input checked="" type="checkbox"/>	O-ring	1.000	St	
12	8410181	<input checked="" type="checkbox"/>	Diaphragm	1.000	St	
13	8407979	<input checked="" type="checkbox"/>	Sealing washer	1.000	St	
14	8412015	<input checked="" type="checkbox"/>	Diaphragm	1.000	St	
15	8410574	<input checked="" type="checkbox"/>	Tube	1.000	St	
16	8410672	<input checked="" type="checkbox"/>	Connecting cap	1.000	St	
17	M20622	<input checked="" type="checkbox"/>	O-ring	1.000	St	only with "Kit Watertrap"

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

23/107

Parts catalog
Pressure Metering Block




Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
1	8408197	<input checked="" type="checkbox"/>	Socket	2.000	St	
2	8412210	<input checked="" type="checkbox"/>	Pressure Measurement Unit	1.000	St	Complete Unit
3	8410713	<input checked="" type="checkbox"/>	O-ring seal	2.000	St	
4	6804141	<input checked="" type="checkbox"/>	Isolate socket	4.000	St	
5	1865870	<input checked="" type="checkbox"/>	Pres.sensor 120mbar diff. RoHS	1.000	St	sale, replaced through 8419627 + 8419628
5	8419628	<input checked="" type="checkbox"/>	Kit pressure sensor EX	1.000	St	replacement for 1865870
5	8419627	<input checked="" type="checkbox"/>	Kit pressure sensor IN	1.000	St	replacement for 1865870
6	D04766	<input checked="" type="checkbox"/>	Packing ring	2.000	St	
7	8306538	<input checked="" type="checkbox"/>	Cable harness,press.absorption	2.000	St	
8	8412993	<input checked="" type="checkbox"/>	Electrovalve	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

24/107

Parts catalog
Patientsystem heater




Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
1	8415975	<input checked="" type="checkbox"/>	Heater	1.000	St	(mit Lüfter u. Ansteuerung/ with fan and controller)
2	8411739	<input checked="" type="checkbox"/>	Insulation foil	1.000	St	
3	8410631	<input checked="" type="checkbox"/>	Spring	1.000	St	
4	8414198	<input checked="" type="checkbox"/>	Exhalation outlet protector	1.000	St	
5	8415973	<input checked="" type="checkbox"/>	Cap	1.000	St	
6	8415971	<input checked="" type="checkbox"/>	Connecting housing	1.000	St	
7	8407689	<input checked="" type="checkbox"/>	Lip seal	1.000	St	
8	8415978	<input checked="" type="checkbox"/>	Dust filter	1.000	St	
9	8415979	<input type="checkbox"/>	Protective guard	1.000	St	
10	8416804	<input checked="" type="checkbox"/>	Kit heating patient part	1.000	St	For patientsystem heater 2005 only

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

25/107

Parts catalog
Electronic



Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
1	8306546	<input checked="" type="checkbox"/>	Cable harness,pneumatics drive	1.000	St	
2	8306225	<input checked="" type="checkbox"/>	Screw	6.000	St	
3	2M17304	<input checked="" type="checkbox"/>	Support	2.000	St	
4	8304117	<input checked="" type="checkbox"/>	Sideparts 4H	2.000	St	
5	1342142	<input checked="" type="checkbox"/>	Screw M3X12 DIN7985	4.000	St	
6	MX44231	<input type="checkbox"/>	PCB Communication	1.000	St	Service drawing
7	8306581	<input checked="" type="checkbox"/>	PCB motherboard	1.000	St	
8	1340727	<input checked="" type="checkbox"/>	pan head screw ISO7045 M3x6-A2-H	21.000	St	


Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

26/107

Parts catalog

Electronic



Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
1	8306546	<input checked="" type="checkbox"/>	Cable harness,pneumatics drive	1.000	St	
2	8306225	<input checked="" type="checkbox"/>	Screw	6.000	St	
8	1340727	<input checked="" type="checkbox"/>	pan head screw ISO7045 M3x6-A2-H	21.000	St	
9	8412384	<input checked="" type="checkbox"/>	Filter, single	1.000	St	
10	MX44232	<input type="checkbox"/>	PCB IFCO-Carrier	1.000	St	Servicedrawing
11	MX44230	<input type="checkbox"/>	PCB Paed.Flow	1.000	St	Servicedrawing
12	MX44226	<input type="checkbox"/>	PCB-CPU	1.000	St	Servicedrawing
13	MX44234	<input type="checkbox"/>	Kit SpO2	1.000	St	
14	MX44228	<input type="checkbox"/>	PCB CO2-Carrier (discontinued)	1.000	St	
15	MX44229	<input type="checkbox"/>	Modification Kit-CO2 Carrier	1.000	St	
16	MX44233	<input type="checkbox"/>	CO2-Accessories	1.000	St	
17	8412374	<input checked="" type="checkbox"/>	Chassis elektronik	1.000	St	
18	MX44662	<input type="checkbox"/>	Power supply Delta AC/DC 2013	1.000	St	
18	MX44001	<input type="checkbox"/>	Power Supply Delta (Frako)	1.000	St	
19	MX44002	<input type="checkbox"/>	Power Supply MB	1.000	St	
20	8421988	<input checked="" type="checkbox"/>	kit 2 batteries Evita internal	1.000	St	
21	1853945	<input checked="" type="checkbox"/>	Cap	1.000	St	
22	1822357	<input checked="" type="checkbox"/>	DSUB-Threaded boltM3 RoHS	2.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

27/107

Parts catalog
Electronic


Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
	8306547	<input checked="" type="checkbox"/>	Plate	4.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

28/107

Parts catalog
PCB Communication




Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
1	8306621	<input type="checkbox"/>	PCB COMMUNICATION	1.000	St	sale
1	8420856	<input checked="" type="checkbox"/>	Rep-Kit PCB Communication42dXL	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

29/107

Parts catalog
PCB IFCO-Carrier


Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
1	8415611	<input checked="" type="checkbox"/>	P.C.B. nursecall	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

30/107

Parts catalog
PCB Paed.Flow




Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
1	8411130	<input checked="" type="checkbox"/>	Neonatal Flow Sensor ISO15	1.000	St	
2	8409626	<input checked="" type="checkbox"/>	Connector cable flow sensor	1.000	St	
3	8306631	<input checked="" type="checkbox"/>	Pba paediatric-flow	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

31/107

Parts catalog
PCB-CPU




Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
1	1835343	<input checked="" type="checkbox"/>	Lithium stor.battery 3V/1400	1.000	St	
2	1852345	<input checked="" type="checkbox"/>	Cap	1.000	St	
3	1837087	<input checked="" type="checkbox"/>	Realtimelock DIL24	1.000	St	
4	8420854	<input checked="" type="checkbox"/>	Rep-Kit LP CPU 42dXL	1.000	St	please D22 with take over
4	8306591	<input type="checkbox"/>	PCB CPU 68332	1.000	St	please D22 with take over/sale

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

32/107

Parts catalog
Kit SpO2




Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
1	8201001	<input type="checkbox"/>	Durasensor DS-100A	1.000	St	
2	8600859	<input checked="" type="checkbox"/>	Extension cable sensor 2,6M	1.000	St	
3	8712065	<input checked="" type="checkbox"/>	Cable clip	2.000	St	
4	8600677	<input type="checkbox"/>	Flat/mounting SPO2-cable	1.000	St	
5	8600787	<input checked="" type="checkbox"/>	Cable harness SPO2	1.000	St	
6	D04766	<input checked="" type="checkbox"/>	Packing ring	3.000	St	
7	1837869	<input checked="" type="checkbox"/>	DC/DC-converter +-15V	1.000	St	
8	8411728	<input type="checkbox"/>	Cable harness SPO2-power	1.000	St	
9	8600481	<input checked="" type="checkbox"/>	SAO2 module	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

33/107

Parts catalog
PCB CO2-Carrier (discontinued)




Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
1	8306611	<input type="checkbox"/>	PCB CO2-carrier	1.000	St	
2	8306567	<input checked="" type="checkbox"/>	Wiring harness CO2-power RS232	1.000	St	
3	8350411	<input checked="" type="checkbox"/>	PCB power	1.000	St	
4	D04766	<input checked="" type="checkbox"/>	Packing ring	1.000	St	
5	8305182	<input checked="" type="checkbox"/>	Spacer bolt	1.000	St	
6	6870296	<input checked="" type="checkbox"/>	PCB signal processor (CO2-main)	1.000	St	
7	6871930	<input checked="" type="checkbox"/>	Software 1.14 (CO2-Main)	1.000	St	Substitute for(6870508 discontinued)
8	8306766	<input checked="" type="checkbox"/>	Cable for CO2 sensor	1.000	St	
9	1837125	<input checked="" type="checkbox"/>	Al. DL Capacitor 0.22F 5V RoHS	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

34/107

Parts catalog
Modification Kit-CO2 Carrier




Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
1	8306567	<input checked="" type="checkbox"/>	Wiring harness CO2-power RS232	1.000	St	
2	8350411	<input checked="" type="checkbox"/>	PCB power	1.000	St	
3	D04766	<input checked="" type="checkbox"/>	Packing ring	1.000	St	
4	8305182	<input checked="" type="checkbox"/>	Spacer bolt	1.000	St	
5	6870296	<input checked="" type="checkbox"/>	PCB signal processor (CO2-main)	1.000	St	
6	6871930	<input checked="" type="checkbox"/>	Software 1.14 (CO2-Main)	1.000	St	
7	8306766	<input checked="" type="checkbox"/>	Cable for CO2 sensor	1.000	St	
8	8414361	<input checked="" type="checkbox"/>	Horn, complete	1.000	St	Mounted after the filter on the back wall
9	8414362	<input checked="" type="checkbox"/>	KIT CO2-Carrier	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

35/107

Parts catalog
CO2-Accessories




Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
1	6870280	<input checked="" type="checkbox"/>	CO2-Cuvette Paed.	1.000	St	
2	6870279	<input checked="" type="checkbox"/>	CO2-Cuvette ERW	1.000	St	
3	MP01062	<input checked="" type="checkbox"/>	Disposable CO2 Cuvette Adult	1.000	St	
4	8412840	<input checked="" type="checkbox"/>	Parking support, cpl.	1.000	St	
5	6871500	<input checked="" type="checkbox"/>	Mainstream sensor	1.000	St	New CO2 Sensor
6	6870343	<input checked="" type="checkbox"/>	Set of inscripted clips	1.000	St	All possible clips with all values
7	6870281	<input checked="" type="checkbox"/>	Test filter	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

36/107

Parts catalog
Power supply Delta AC/DC 2013




Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
1	8421079	<input checked="" type="checkbox"/>	Repair kit AC/DC power supply unit	1.000	St	
2	8415955	<input checked="" type="checkbox"/>	Kit flap for power switch	1.000	St	
3	8421988	<input checked="" type="checkbox"/>	kit 2 batteries Evita internal	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

37/107

Parts catalog
Power supply Delta AC/DC 2013




Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
4	1851810	<input checked="" type="checkbox"/>	Cabel Australia,3m,10A,C13W	1.000	St	
5	1856553	<input checked="" type="checkbox"/>	Power cord CE,3m,10A,C13W	1.000	St	
6	1866788	<input checked="" type="checkbox"/>	Cable Switzerland,3m,C13 RoHS	1.000	St	
7	1856626	<input checked="" type="checkbox"/>	Power cord USA,3m,10A,C13W	1.000	St	
8	1856561	<input checked="" type="checkbox"/>	Power cord CH,3m,10A,C13W	1.000	St	
9	1868950	<input checked="" type="checkbox"/>	PWR Cord DK,3m,10A,C13,HG RoHS	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

38/107

Parts catalog
Power supply Delta AC/DC 2013



Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
10	1856596	<input checked="" type="checkbox"/>	Cable Great Britain 3m,C13W	1.000	St	
11	1878840	<input checked="" type="checkbox"/>	Pwr. cord IL,3m,10A,C13W RoHS	1.000	St	
12	1875973	<input checked="" type="checkbox"/>	Mains cable SA,3m,C13W,GY RoHS	1.000	St	
13	1859706	<input checked="" type="checkbox"/>	Power cable 10A,3m,black,China	1.000	St	
14	1866915	<input checked="" type="checkbox"/>	Power cable 10A 3m bk RoHS	1.000	St	
15	1875531	<input checked="" type="checkbox"/>	Mains cable BR,3m,C13W,BK RoHS	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

39/107

Parts catalog
Power supply Delta AC/DC 2013


Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
0	8421078	<input checked="" type="checkbox"/>	Kit connecting housing f. ext.	1.000	St	only for option external battery
0	8421072	<input checked="" type="checkbox"/>	Stain relief, 30 degrees	1.000	St	only for CH and DK
0	8422125	<input checked="" type="checkbox"/>	kit 2 batteries Evita external	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

40/107

Parts catalog
Power Supply Delta (Frako)




Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
1	8415955	<input checked="" type="checkbox"/>	Kit flap for power switch	1.000	St	
2	1843680	<input checked="" type="checkbox"/>	Pull relief	1.000	St	
3	1834967	<input checked="" type="checkbox"/>	Ferrite core DM6-7 LV01583	1.000	St	
4	1851683	<input checked="" type="checkbox"/>	Supply main, 3m, 10A	1.000	St	
4	1841793	<input checked="" type="checkbox"/>	PWR Cord 10A,3m,gr,USA/J RoHS	1.000	St	
4	1851713	<input checked="" type="checkbox"/>	Cable Great Britian,3m,10A	1.000	St	
4	1851691	<input checked="" type="checkbox"/>	Mains cable Swiss ,3m,10A	1.000	St	
4	1851705	<input checked="" type="checkbox"/>	Cable Australia,3m,10A,C13	1.000	St	
4	1851721	<input checked="" type="checkbox"/>	Power cable DK, 3 m, 10 A	1.000	St	
5	8306520	<input type="checkbox"/>	AC-modul power supply Evita 4	1.000	St	without DC-module option/ sale
6	8419059	<input type="checkbox"/>	Opt. AC/DC pwrpack w/o batt.	1.000	St	with DC-module option

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

41/107

Parts catalog
Power Supply Delta (Frako)




Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
8	8421988	<input checked="" type="checkbox"/>	kit 2 batteries Evita internal	1.000	St	
9	8300358	<input checked="" type="checkbox"/>	Cable clip (3,6X140)	1.000	St	
10	8419275	<input checked="" type="checkbox"/>	Kit PCB Battery Type B (Evita)	1.000	St	
11	8417298	<input checked="" type="checkbox"/>	Kit PCB Battery Type B cpl.	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

42/107

Parts catalog
Power Supply MB




Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
1	8415638	<input checked="" type="checkbox"/>	Drop Flap	1.000	St	
2	1855549	<input type="checkbox"/>	Fuse 6,3A IEC60127-2/I	2.000	St	
3	1843680	<input checked="" type="checkbox"/>	Pull relief	1.000	St	
4	8415639	<input checked="" type="checkbox"/>	Battery Drawer	1.000	St	
5	8421988	<input checked="" type="checkbox"/>	kit 2 batteries Evita internal	1.000	St	
6	8415636	<input type="checkbox"/>	Set of cable	1.000	St	
7	1855530	<input checked="" type="checkbox"/>	Fuse 20A/80V	1.000	St	
8	8415572	<input checked="" type="checkbox"/>	Airfilter	1.000	St	
9	8415637	<input checked="" type="checkbox"/>	Filter covering	2.000	St	
10	8306225	<input checked="" type="checkbox"/>	Screw	2.000	St	
11	8414577	<input type="checkbox"/>	Power pack USA	1.000	St	nur für Bundeswehr/ only for German army

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

43/107

Parts catalog
Operating Device




Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
1	8416980	<input checked="" type="checkbox"/>	Op.Device Evita 4 edition	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts
Evita 4/ Evita 4 edition
Revision: 22

Parts catalog

Operating Device




Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
2	8412380	<input checked="" type="checkbox"/>	Stripe	1.000	St	
3	M29655	<input checked="" type="checkbox"/>	Control knop	1.000	St	
4	M29952	<input checked="" type="checkbox"/>	Plate spring	1.000	St	
5	8416960	<input checked="" type="checkbox"/>	Housing, blue	1.000	St	
6	8412602	<input checked="" type="checkbox"/>	Outer frame	1.000	St	
7	8412626	<input checked="" type="checkbox"/>	Insert nut	4.000	St	
8	M09360	<input checked="" type="checkbox"/>	Spring	4.000	St	
9	8412625	<input type="checkbox"/>	Guide tappet	4.000	St	
10	8412632	<input checked="" type="checkbox"/>	Release taste	2.000	St	
11	M14362	<input checked="" type="checkbox"/>	Spring	2.000	St	
12	8412884	<input checked="" type="checkbox"/>	Clamp	1.000	St	
13	1340727	<input checked="" type="checkbox"/>	pan head screw ISO7045 M3x6-A2-H	6.000	St	
14	D04766	<input checked="" type="checkbox"/>	Packing ring	4.000	St	
15	1341790	<input checked="" type="checkbox"/>	Oval head screw DIN7985-M4X12-A2	4.000	St	
16	8416965	<input checked="" type="checkbox"/>	Keyboard Evita 4 edition	1.000	St	
16A	8412137	<input type="checkbox"/>	Display, cpl.	1.000	St	keyboard Evita 4
17	8306638	<input checked="" type="checkbox"/>	Touchscreen	1.000	St	
18	8411746	<input checked="" type="checkbox"/>	Touchrahmen, kompl.	1.000	St	8412941-seal does not belong to the scope of delivery
20	8412941	<input checked="" type="checkbox"/>	O-ring seal	1.000	St	
21	1339974	<input checked="" type="checkbox"/>	Screw M4X16 DIN965	6.000	St	
22	1822357	<input checked="" type="checkbox"/>	DSUB-Threaded boltM3 RoHS	1.000	St	
23	1839780	<input checked="" type="checkbox"/>	Connecting cable	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

Parts catalog
Operating Device




Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
24	8420857	<input checked="" type="checkbox"/>	Rep-Kit Grafikcontroller 4	1.000	St	
24	8306641	<input type="checkbox"/>	PCB graphic controller	1.000	St	
43	8412631	<input checked="" type="checkbox"/>	Roller	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

46/107

Parts catalog
Operating Device




Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
4	M29952	<input checked="" type="checkbox"/>	Plate spring	1.000	St	
13	1340727	<input checked="" type="checkbox"/>	pan head screw ISO7045 M3x6-A2-H	6.000	St	
15	1341790	<input checked="" type="checkbox"/>	Oval head screw DIN7985-M4X12-A2	4.000	St	
21	1339974	<input checked="" type="checkbox"/>	Screw M4X16 DIN965	6.000	St	
25	8412885	<input checked="" type="checkbox"/>	Plate	4.000	St	
26	8412844	<input checked="" type="checkbox"/>	socket	1.000	St	
27	8412637	<input checked="" type="checkbox"/>	Toggle	1.000	St	
28	E24916	<input checked="" type="checkbox"/>	Spring	1.000	St	
29	8412843	<input checked="" type="checkbox"/>	Sheet metal	1.000	St	
30	6800478	<input checked="" type="checkbox"/>	Buffer	1.000	St	
31	8412703	<input checked="" type="checkbox"/>	Springplate	1.000	St	
32	1327542	<input checked="" type="checkbox"/>	Washer ISO 7089-4-200HV-A4	4.000	St	
33	8412627	<input checked="" type="checkbox"/>	Bar	1.000	St	
34	8412622	<input checked="" type="checkbox"/>	Lock	2.000	St	
35	M14164	<input checked="" type="checkbox"/>	Closing spring	4.000	St	
36	1263056	<input checked="" type="checkbox"/>	Hexagon socket head cap screw	1.000	St	
37	8412611	<input checked="" type="checkbox"/>	Housing (hinge)	2.000	St	
38	8412898	<input checked="" type="checkbox"/>	Cable harness touch-SC(Evita 4	1.000	St	
39	8306565	<input checked="" type="checkbox"/>	Shaft encoder,cpl.	1.000	St	
40	1334913	<input checked="" type="checkbox"/>	Hexagon nut ISO 4032-M3-A4	3.000	St	
41	8306558	<input type="checkbox"/>	Loudspeaker, cpl.	1.000	St	
42	8412808	<input checked="" type="checkbox"/>	Spacer	4.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

47/107

Parts catalog
Operating Device

Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
17	8306638	<input checked="" type="checkbox"/>	Touchscreen	1.000	St	
18	8411746	<input checked="" type="checkbox"/>	Touchrahmen, kompl.	1.000	St	8412941-seal does not belong to the scope of delivery

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

48/107

Parts catalog
Operating Device


Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
	MX44236	<input type="checkbox"/>	Display Toshiba Typ1	1.000	St	sale
	MX44238	<input type="checkbox"/>	Display NEC	1.000	St	sale
	MX44237	<input type="checkbox"/>	Display Toshiba Typ0	1.000	St	
	MX44240	<input type="checkbox"/>	Plate Front Evita 4	1.000	St	Servicedrawing
	MX44239	<input type="checkbox"/>	Display Sharp	1.000	St	sale

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

49/107

Parts catalog
Plate Front Evita 4



Pos. 0 ohne Abbildung/ Item 0 without illustration/ Pos. 0 sin ilustración/ Pos. 0 sans illustration


Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
1	1311522	<input checked="" type="checkbox"/>	Screw B2,9X13 DIN7973	1.000	St	
2	8412884	<input checked="" type="checkbox"/>	Clamp	1.000	St	
3	1822357	<input checked="" type="checkbox"/>	DSUB-Threaded boltM3 RoHS	2.000	St	
4	8414815	<input checked="" type="checkbox"/>	Front plate	1.000	St	
5	8306775	<input checked="" type="checkbox"/>	Cable harness,front (Evita 4)	1.000	St	
6	8412839	<input checked="" type="checkbox"/>	Front plate, cpl	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

50/107

Parts catalog
Display Toshiba Typ1




Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
0	8420794	<input checked="" type="checkbox"/>	Rep-Kit Display AUO (Evita 4)	1.000	St	consider please conversion instructions
1	D04766	<input checked="" type="checkbox"/>	Packing ring	6.000	St	
2	8306640	<input checked="" type="checkbox"/>	Multiplexer color Toshiba	1.000	St	
3	8306795	<input checked="" type="checkbox"/>	Lookup TFT-Toshiba	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

51/107

Parts catalog
Display Toshiba Typ0




Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
0	8416801	<input type="checkbox"/>	Kit Display Toshiba E4	1.000	St	TSB/IDM Nr.57
0	8309126	<input checked="" type="checkbox"/>	Cable harness Toshiba Evita 4	1.000	St	
0	8420794	<input checked="" type="checkbox"/>	Rep-Kit Display AUO (Evita 4)	1.000	St	consider please conversion instructions
1	8712065	<input checked="" type="checkbox"/>	Cable clip	1.000	St	
2	D04766	<input checked="" type="checkbox"/>	Packing ring	4.000	St	
3	8411859	<input checked="" type="checkbox"/>	Multiplexer color Toshiba	1.000	St	
4	8306795	<input checked="" type="checkbox"/>	Lookup TFT-Toshiba	1.000	St	
5	1843559	<input checked="" type="checkbox"/>	DC/AC-inverter 12V/550V RoHS	1.000	St	
6	1868594	<input checked="" type="checkbox"/>	TFT-Display 10,4" 640X480 RoHS	1.000	St	
7	8411765	<input checked="" type="checkbox"/>	Distance bolt	4.000	St	
8	D20682	<input checked="" type="checkbox"/>	Gasket	4.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

52/107

Parts catalog
Display NEC




Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
0	8420794	<input checked="" type="checkbox"/>	Rep-Kit Display AUO (Evita 4)	1.000	St	consider please conversion instructions
1	8300428	<input checked="" type="checkbox"/>	Spacer	4.000	St	
2	1841610	<input type="checkbox"/>	LCD back lighting (Nec 2.vers.	1.000	St	
3	8411715	<input type="checkbox"/>	Holder, left	1.000	St	
4	8411686	<input checked="" type="checkbox"/>	Multiplexer NEC	1.000	St	
5	8306795	<input checked="" type="checkbox"/>	Lookup TFT-Toshiba	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

53/107

Parts catalog
Display Sharp




Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
0	8420794	<input checked="" type="checkbox"/>	Rep-Kit Display AUO (Evita 4)	1.000	St	consider please conversion instructions
1	8411732	<input checked="" type="checkbox"/>	Inverter Sharp	1.000	St	
2	1843095	<input type="checkbox"/>	Back lighting (Sharp)	1.000	St	
3	8712065	<input checked="" type="checkbox"/>	Cable clip	3.000	St	
4	8411765	<input checked="" type="checkbox"/>	Distance bolt	2.000	St	
5	8411753	<input type="checkbox"/>	Screw	1.000	St	
6	8300428	<input checked="" type="checkbox"/>	Spacer	4.000	St	
7	8411776	<input checked="" type="checkbox"/>	Multiplexer Sharp	1.000	St	
8	8306795	<input checked="" type="checkbox"/>	Lookup TFT-Toshiba	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

54/107

Parts catalog
Modification Kits/Options




Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
8413780		<input type="checkbox"/>	Kit Capno Plus	1.000	St	Bestellung über Konfi.-Nr: 8414240
8411735		<input type="checkbox"/>	Kit communication	1.000	St	Bestellung über Konfi.-Nr: 8414240
8413562		<input type="checkbox"/>	Kit Breathing Support Package	1.000	St	Bestellung über Konfi.-Nr: 8414240
8413563		<input type="checkbox"/>	Kit Neoflow	1.000	St	Bestellung über Konfi.-Nr: 8414240
8414474		<input type="checkbox"/>	Kit mask ventilation	1.000	St	Bestellung über Konfi.-Nr: 8414240
8414472		<input type="checkbox"/>	Kit Evita Remote	1.000	St	Bestellung über Konfi.-Nr: 8414240
8414476		<input type="checkbox"/>	Kit nurse call	1.000	St	Bestellung über Konfi.-Nr: 8414240
8414879		<input type="checkbox"/>	Kit Evita 4 zu EvitaXL	1.000	St	
8414556		<input type="checkbox"/>	Kit ATC	1.000	St	Bestellung über Konfi.-Nr: 8414240
8413034		<input type="checkbox"/>	Kit DC-modul	1.000	St	Bestellung über Konfi.-Nr: 8414240
8414359		<input checked="" type="checkbox"/>	Kit muted loudspeaker	1.000	St	
MX44336		<input type="checkbox"/>	Evita 4 edition	1.000	St	
MX44209		<input type="checkbox"/>	Evita mobil with Transport	1.000	St	
MX44210		<input type="checkbox"/>	Mobile Trolley MSU	1.000	St	
MX44252		<input type="checkbox"/>	Hinged Bracket 240-Degree	1.000	St	
MX44255		<input type="checkbox"/>	Water Seperator, cpl.	1.000	St	
MX44254		<input type="checkbox"/>	High Pressure Water Trap AIR	1.000	St	
MX44253		<input type="checkbox"/>	Pneum. Med. Nebulizer Evita 4	1.000	St	
8414875		<input type="checkbox"/>	Kit PPS	1.000	St	order via Confi.-No: 8414240

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

55/107

Parts catalog
Modification Kits/Options




Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
8415570		<input type="checkbox"/>	Kit nurse call/pressure sensor	1.000	St	order via Confi.-No: 8414240
8419417		<input type="checkbox"/>	Kit SW 4.22 Evita 4	1.000	St	order via Confi.-No: 8414240
MX44241		<input type="checkbox"/>	Two-Column Trolley to 1997	1.000	St	discontinued
8419419		<input checked="" type="checkbox"/>	Kit SW4.22 Evita 4 w/o UM	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

56/107

Parts catalog
 Evita 4 edition



Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
1	8414828	<input checked="" type="checkbox"/>	Tray	1.000	St	
2	8416960	<input checked="" type="checkbox"/>	Housing, blue	1.000	St	
3	8416965	<input checked="" type="checkbox"/>	Keyboard Evita 4 edition	1.000	St	
4	1868594	<input checked="" type="checkbox"/>	TFT-Display 10,4" 640X480 RoHS	1.000	St	
5	8419064	<input checked="" type="checkbox"/>	Evita 4 label	1.000	St	nur für Mittel- und Südamerika/ only for Middle- and Southamerika
6	8416969	<input checked="" type="checkbox"/>	Keypad strip kit, EL	1.000	St	
6	8416974	<input checked="" type="checkbox"/>	Keypad strip kit, nl	1.000	St	
6	8416968	<input checked="" type="checkbox"/>	set insertion strip ZH	1.000	St	
6	8416973	<input checked="" type="checkbox"/>	Keypad strip kit, es	1.000	St	
6	8416972	<input checked="" type="checkbox"/>	Keypad strip kit, it	1.000	St	
6	8416970	<input checked="" type="checkbox"/>	Keypad strip kit, de	1.000	St	
6	8416978	<input checked="" type="checkbox"/>	Keypad strip kit, ja	1.000	St	
6	8416975	<input checked="" type="checkbox"/>	Keypad strip kit, sv	1.000	St	
6	8416976	<input checked="" type="checkbox"/>	Keypad strip kit, fr	1.000	St	
6	8416977	<input checked="" type="checkbox"/>	Keypad strip kit, enUS	1.000	St	
6	8416967	<input checked="" type="checkbox"/>	Keypad strip kit, PT	1.000	St	
6	8416966	<input checked="" type="checkbox"/>	Keypad strip kit, RU	1.000	St	
6	8416971	<input checked="" type="checkbox"/>	Keypad strip kit, enGB	1.000	St	
7	8416980	<input checked="" type="checkbox"/>	Op.Device Evita 4 edition	1.000	St	
8	MX44257	<input type="checkbox"/>	Frontplate	1.000	St	auch für Evita 4 edition/also for Evita 4 edition

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

 Evita 4/ Evita 4 edition
 Revision: 22

57/107

Parts catalog
Frontplate




Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
1	1340727	<input checked="" type="checkbox"/>	pan head screw ISO7045 M3x6-A2-H	1.000	St	
2	1309749	<input checked="" type="checkbox"/>	Split washer B3 DIN127-X12CRNI	1.000	St	
3	8306775	<input checked="" type="checkbox"/>	Cable harness,front (Evita 4)	1.000	St	
4	1822357	<input checked="" type="checkbox"/>	DSUB-Threaded boltM3 RoHS	2.000	St	
5	8414815	<input checked="" type="checkbox"/>	Front plate	1.000	St	
6	8412884	<input checked="" type="checkbox"/>	Clamp	1.000	St	
7	1311522	<input checked="" type="checkbox"/>	Screw B2,9X13 DIN7973	1.000	St	
8	8412839	<input checked="" type="checkbox"/>	Front plate, cpl	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

58/107

Parts catalog
Evita mobil with Transport




Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
1	8411954	<input checked="" type="checkbox"/>	Cover plate for top	1.000	St	
2	8411957	<input checked="" type="checkbox"/>	Handle	1.000	St	
3	8414229	<input checked="" type="checkbox"/>	Cylinder rack, compl.	1.000	St	
4	8411972	<input checked="" type="checkbox"/>	Hook and loop fastener	1.000	St	
5	8411973	<input checked="" type="checkbox"/>	Inlet	1.000	St	
6	8411955	<input checked="" type="checkbox"/>	cover for battery	1.000	St	
7	8414226	<input checked="" type="checkbox"/>	Castor with brake	1.000	St	
8	8414228	<input checked="" type="checkbox"/>	Roller, unlocked	1.000	St	
9	8414227	<input checked="" type="checkbox"/>	Castor	1.000	St	
10	8411518	<input checked="" type="checkbox"/>	Slide, cpl.	1.000	St	
11	M16452	<input checked="" type="checkbox"/>	Tension spring	1.000	St	
12	8411822	<input checked="" type="checkbox"/>	KIT battery cable Evita Mobil	1.000	St	
13	8414424	<input checked="" type="checkbox"/>	Handle	1.000	St	
14	1314440	<input checked="" type="checkbox"/>	Split washer 8 DIN7980-X12CRNI	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

59/107

Parts catalog
Evita mobil with Transport




Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
15	8414484	<input type="checkbox"/>	Evita TransportMobil	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

60/107

Parts catalog
Evita mobil with Transport

Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
16	1330705	<input checked="" type="checkbox"/>	Washer ISO 7089-6-200HV-A4	1.000	St	
17	8419252	<input checked="" type="checkbox"/>	Screw	1.000	St	
18	8414495	<input checked="" type="checkbox"/>	Base (Transportmobile)	1.000	St	
19	8414496	<input checked="" type="checkbox"/>	Top section (Transportmobile)	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

61/107

Parts catalog
Evita mobil with Transport


Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
1326279		<input checked="" type="checkbox"/>	Washer ISO 7089-5-200HV-A4	1.000	St	
1283081		<input checked="" type="checkbox"/>	Hexagon nut ISO 4032-M5-A4	1.000	St	
8414437		<input checked="" type="checkbox"/>	Kit Evita TransportMobil	1.000	St	see modification instructions
1335820		<input checked="" type="checkbox"/>	Ch.head screw AM5X16 DIN84-A2	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

62/107

Parts catalog
Mobile Trolley MSU





Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
1	8411965	<input checked="" type="checkbox"/>	Evita mobil, short	1.000	St	
2	M16452	<input checked="" type="checkbox"/>	Tension spring	2.000	St	
3	8411954	<input checked="" type="checkbox"/>	Cover plate for top	1.000	St	
4	8411518	<input checked="" type="checkbox"/>	Slide, cpl.	1.000	St	
5	1282123	<input checked="" type="checkbox"/>	Hexagon socket head cap screw	1.000	St	
6	8411957	<input checked="" type="checkbox"/>	Handle	1.000	St	
7	8411973	<input checked="" type="checkbox"/>	Inlet	1.000	St	
8	8411955	<input checked="" type="checkbox"/>	cover for battery	1.000	St	
9	8411945	<input checked="" type="checkbox"/>	Castor with fixing	2.000	St	
10	8411946	<input checked="" type="checkbox"/>	Castor	2.000	St	
11	8419252	<input checked="" type="checkbox"/>	Screw	1.000	St	
12	1330705	<input checked="" type="checkbox"/>	Washer ISO 7089-6-200HV-A4	1.000	St	
13	8411950	<input checked="" type="checkbox"/>	Evita mobil, tall	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

63/107

Parts catalog
Mobile Trolley MSU

Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
14	8414495	<input checked="" type="checkbox"/>	Base (Transportmobile)	1.000	St	
15	8414496	<input checked="" type="checkbox"/>	Top section (Transportmobile)	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

64/107

Parts catalog
Mobile Trolley MSU


Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
	MX44242	<input type="checkbox"/>	Cylinder Support	1.000	St	
	MX44243	<input type="checkbox"/>	Humidifier Rack	1.000	St	
	MX44244	<input type="checkbox"/>	Kit Socket Strip Evita Mobil	1.000	St	Servicedrawing
	MX44245	<input type="checkbox"/>	Kit Battery Cable Evita Mobil	1.000	St	Servicedrawing

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

65/107

Parts catalog
Kit Socket Strip Evita Mobil




Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
1	8712065	<input checked="" type="checkbox"/>	Cable clip	1.000	St	
2	8409310	<input checked="" type="checkbox"/>	Socket	3.000	St	
3	1815148	<input checked="" type="checkbox"/>	Fuse link DIN 41662 T3,15A	6.000	St	
4	8411814	<input checked="" type="checkbox"/>	Kit Power socket strip Evita Mobil	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

66/107

Parts catalog
Kit Battery Cable Evita Mobil


Pos. 0 (wenn vorhanden) ohne Abbildung/ Item 0 (if available) without illustration/
Pos. 0 (si existe) sin ilustración/ Pos. 0 (si existent) sans illustration

Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
1	8422125	<input checked="" type="checkbox"/>	kit 2 batteries Evita external	1.000	St	
2	8411822	<input checked="" type="checkbox"/>	KIT battery cable Evita Mobil	1.000	St	
3	1843133	<input checked="" type="checkbox"/>	Blade fuse 25A 32V RoHS	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

67/107

Parts catalog
Cylinder Support


Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
1	8411970	<input checked="" type="checkbox"/>	Set cylinder bracket ev.mobil	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

68/107

Parts catalog
Humidifier Rack




Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
1	8411956	<input checked="" type="checkbox"/>	Humidifier bracket Evita mobil	1.000	St	
2	MP01002	<input checked="" type="checkbox"/>	Mounting bracket F&P	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

69/107

Parts catalog
Hinged Bracket 240-Degree




Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
1	8410383	<input checked="" type="checkbox"/>	Set clamping lever	1.000	St	
2	8411424	<input checked="" type="checkbox"/>	Rail catch	1.000	St	
3	8410382	<input checked="" type="checkbox"/>	Star grip with disc	2.000	St	
4	8410075	<input checked="" type="checkbox"/>	Binder, cpl.	1.000	St	
5	8409746	<input checked="" type="checkbox"/>	Support for hinged arm A	1.000	St	
6	8409841	<input checked="" type="checkbox"/>	Hose clamp for hinged arm	1.000	St	
7	8409609	<input checked="" type="checkbox"/>	Hinged arm	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

70/107

Parts catalog
Water Seperator, cpl.


Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
1	M13142	<input checked="" type="checkbox"/>	o-ring	1.000	St	
2	M32027	<input checked="" type="checkbox"/>	Packing ring	1.000	St	
3	8414271	<input checked="" type="checkbox"/>	Air-Watertrap big NIST/NIST	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

71/107

Parts catalog
High Pressure Water Trap AIR




Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
1	8414274	<input checked="" type="checkbox"/>	Air water trap medium	1.000	St	
2	M13142	<input checked="" type="checkbox"/>	o-ring	1.000	St	
3	M32027	<input checked="" type="checkbox"/>	Packing ring	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

72/107

Parts catalog
Pneum. Med. Nebulizer Evita 4




Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
1	8405046	<input checked="" type="checkbox"/>	Cap	1.000	St	
2	2M10633	<input checked="" type="checkbox"/>	O-ring	1.000	St	
3	8412934	<input checked="" type="checkbox"/>	Atomizer housing, cpl.	1.000	St	
4	8407291	<input checked="" type="checkbox"/>	Gasket	1.000	St	
5	R27810	<input checked="" type="checkbox"/>	Sealing ring	1.000	St	
6	8412932	<input checked="" type="checkbox"/>	Nozzle	1.000	St	
7	E20566	<input checked="" type="checkbox"/>	O-ring 7x1,5	1.000	St	
8	8404979	<input checked="" type="checkbox"/>	Atomizer	1.000	St	
9	8406584	<input checked="" type="checkbox"/>	Receptacle	1.000	St	
10	8412985	<input checked="" type="checkbox"/>	Nebulizer hose	1.000	St	
11	8412935	<input checked="" type="checkbox"/>	Pneum. Medication Nebulizer	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

73/107

Parts catalog
Two-Column Trolley to 1997




Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
	MX44246	<input type="checkbox"/>	Support	1.000	St	Servicedrawing
	MX44248	<input type="checkbox"/>	O2-Distributor Mod. Kit Evita	1.000	St	Servicedrawing
	MX44247	<input type="checkbox"/>	Modification Kit Socket	1.000	St	Servicedrawing
	MX44249	<input type="checkbox"/>	Fixing Set Cupboard Evita	1.000	St	Servicedrawing
	MX44250	<input type="checkbox"/>	Receiver	1.000	St	Servicedrawing
1	8408904	<input type="checkbox"/>	Sheet metal	1.000	St	
2	M29200	<input checked="" type="checkbox"/>	Tray	2.000	St	
3	2M19130	<input type="checkbox"/>	Guide roll	2.000	St	
4	2M19140	<input checked="" type="checkbox"/>	Fixing guide roll	2.000	St	
5	8408902	<input checked="" type="checkbox"/>	Pin	1.000	St	
6	M29336	<input checked="" type="checkbox"/>	Profile	2.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

74/107

Parts catalog
Support


Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
1	M29110	<input checked="" type="checkbox"/>	Plate	1.000	St	
2	8412219	<input checked="" type="checkbox"/>	Support	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

75/107

Parts catalog
O2-Distributor Mod. Kit Evita




Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
1	G12065	<input checked="" type="checkbox"/>	L-Swivel connection	4.000	St	
2	8410859	<input checked="" type="checkbox"/>	Modification kit drop signal	1.000	St	
3	8403058	<input checked="" type="checkbox"/>	Bow	2.000	St	
4	8405629	<input type="checkbox"/>	Connecting line	1.000	St	
5	1190520	<input checked="" type="checkbox"/>	Hose 4X1,5-SI 50 SH A NF	1.000	m	
6	M26349	<input checked="" type="checkbox"/>	Hook	1.000	St	
7	RM07823	<input checked="" type="checkbox"/>	O-ring	2.000	St	
8	8403322	<input type="checkbox"/>	Screw	1.000	St	
9	M25210	<input checked="" type="checkbox"/>	Screw-in plug-type connection	2.000	St	
10	8402441	<input checked="" type="checkbox"/>	3/2-Port distributing valve	1.000	St	
11	D11521	<input checked="" type="checkbox"/>	Mesh bottom	1.000	St	
12	8409010	<input type="checkbox"/>	O2 distributor inst.set Evita	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

77/107

Parts catalog
Fixing Set Cupboard Evita




Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
1	8409018	<input type="checkbox"/>	Cabinet attachment set Evita	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

78/107

Parts catalog
Receiver


Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
1	MX08108	<input type="checkbox"/>	Hinge	2.000	St	
2	8408959	<input checked="" type="checkbox"/>	Fixing device	1.000	St	
3	8409015	<input type="checkbox"/>	Drawer, cpl.	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

79/107

Parts catalog
Maintenance Parts/Service Kits




Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
0	MP03904	<input type="checkbox"/>	Filter Bowl	1.000	St	
0	MP03901	<input type="checkbox"/>	AIR-High Pressure Line Filter	1.000	St	only for USA
0	MP03902	<input type="checkbox"/>	O2-High Pressure Line Filter	1.000	St	only for USA
0	MP03903	<input type="checkbox"/>	Filter Element	1.000	St	
1	MX40018	<input type="checkbox"/>	Evita XL/4/2d/2/S2, yearly	1.000	St	
2	MX44260	<input type="checkbox"/>	Evita 2 dura/4, 2-yearly	1.000	St	
3	MX40035	<input type="checkbox"/>	Evita XL/4/2d Draeger,2-yearly	1.000	St	
4	MX40034	<input type="checkbox"/>	Evita XL/4/2d FAS, 2-yearly	1.000	St	
5	MX41894	<input type="checkbox"/>	Evita XL/4/2d Dräger, 6-yearly	1.000	St	
6	MX40061	<input type="checkbox"/>	Evita XL/4/2d FAS, 6-yearly	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

80/107

Parts catalog
Evita XL/4/2d/2/S2, yearly




Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
1	MX08220	<input checked="" type="checkbox"/>	Evita 2/S2/2Dura/4/XL Set 1y	1.000	St	Set kompl.
2	8410181	<input checked="" type="checkbox"/>	Diaphragm	1.000	St	
3	8407979	<input checked="" type="checkbox"/>	Sealing washer	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

81/107

Parts catalog
Evita 2 dura/4, 2-yearly




Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
2	8407689	<input checked="" type="checkbox"/>	Lip seal	2.000	St	
3	M20622	<input checked="" type="checkbox"/>	O-ring	3.000	St	
4	8408208	<input checked="" type="checkbox"/>	Filter (for 8411848)	2.000	St	
5	1835343	<input checked="" type="checkbox"/>	Lithium stor.battery 3V/1400	1.000	St	
6	8408204	<input checked="" type="checkbox"/>	Gasket	1.000	St	
7	8411516	<input checked="" type="checkbox"/>	O-ring for 8411848	2.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

82/107

Parts catalog
Evita XL/4/2d Draeger, 2-yearly




Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
1	MX08226	<input checked="" type="checkbox"/>	EvitaXL/4/2d Gasinlet Draeg.2y	1.000	St	Set kompl.
2	8410181	<input checked="" type="checkbox"/>	Diaphragm	1.000	St	
3	8407979	<input checked="" type="checkbox"/>	Sealing washer	1.000	St	
4	1835343	<input checked="" type="checkbox"/>	Lithium stor.battery 3V/1400	1.000	St	
5	8408204	<input checked="" type="checkbox"/>	Gasket	1.000	St	
6	8408208	<input checked="" type="checkbox"/>	Filter (for 8411848)	2.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

83/107

Parts catalog
Evita XL/4/2d FAS, 2-yearly




Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
1	MX08225	<input checked="" type="checkbox"/>	EvitaXL/4/2d Gasinlet FAS 2Y	1.000	St	Set kompl.
2	8410181	<input checked="" type="checkbox"/>	Diaphragm	1.000	St	
3	8407979	<input checked="" type="checkbox"/>	Sealing washer	1.000	St	
4	1835343	<input checked="" type="checkbox"/>	Lithium stor.battery 3V/1400	1.000	St	
5	8411516	<input checked="" type="checkbox"/>	O-ring for 8411848	2.000	St	
6	8408208	<input checked="" type="checkbox"/>	Filter (for 8411848)	2.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

84/107

Parts catalog
Evita XL/4/2d Dräger, 6-yearly




Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
1	MX08228	<input checked="" type="checkbox"/>	EvitaXL/4/2d Gasinl.Draeg. 6y	1.000	St	Set kompl.
2	8410181	<input checked="" type="checkbox"/>	Diaphragm	1.000	St	
3	8407979	<input checked="" type="checkbox"/>	Sealing washer	1.000	St	
4	1835343	<input checked="" type="checkbox"/>	Lithium stor.battery 3V/1400	1.000	St	
5	8414179	<input checked="" type="checkbox"/>	Diaphragm connection SV	1.000	St	
6	8306560	<input checked="" type="checkbox"/>	Fan pneumatic, compl.	1.000	St	
7	8408205	<input checked="" type="checkbox"/>	Pressure regulator	2.000	St	
8	8415709	<input checked="" type="checkbox"/>	Cable harness spirolog sensor	1.000	St	
9	8408204	<input checked="" type="checkbox"/>	Gasket	1.000	St	
10	1837087	<input checked="" type="checkbox"/>	Realtimeclock DIL24	1.000	St	
11	8408208	<input checked="" type="checkbox"/>	Filter (for 8411848)	2.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

85/107

Parts catalog
Evita XL/4/2d FAS, 6-yearly




Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
1	MX08227	<input checked="" type="checkbox"/>	EvitaXL/4/2D Gasinlet FAS 6y	1.000	St	Set kompl.
2	8410181	<input checked="" type="checkbox"/>	Diaphragm	1.000	St	
3	8407979	<input checked="" type="checkbox"/>	Sealing washer	1.000	St	
4	1835343	<input checked="" type="checkbox"/>	Lithium stor.battery 3V/1400	1.000	St	
5	8306560	<input checked="" type="checkbox"/>	Fan pneumatic, compl.	1.000	St	
6	8411513	<input checked="" type="checkbox"/>	Diaphragm FPM (for 8411848)	2.000	St	
7	8414179	<input checked="" type="checkbox"/>	Diaphragm connection SV	1.000	St	
8	8411514	<input checked="" type="checkbox"/>	Valve tappet (for 8411848)	2.000	St	
9	1837087	<input checked="" type="checkbox"/>	Realtimeclock DIL24	1.000	St	
10	8415709	<input checked="" type="checkbox"/>	Cable harness spiolog sensor	1.000	St	
11	8411516	<input checked="" type="checkbox"/>	O-ring for 8411848	2.000	St	
12	8408208	<input checked="" type="checkbox"/>	Filter (for 8411848)	2.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

86/107

Parts catalog
Accessories/Consumables




Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
	MX44202	<input type="checkbox"/>	Order List	1.000	St	Servicedrawing
	MX44263	<input type="checkbox"/>	Accessories Adult Ventilation	1.000	St	Servicedrawing
	MX44264	<input type="checkbox"/>	Accessories Paed. Ventilation	1.000	St	Servicedrawing
	MX44089	<input type="checkbox"/>	Calibration Set CO2	1.000	St	Servicedrawing
	MX44358	<input type="checkbox"/>	Aquapor EL	1.000	St	
	MX44649	<input type="checkbox"/>	CP parts	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

87/107

Parts catalog
Order List



Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
	8411900	<input type="checkbox"/>	Evita 4, Construction Kit	1.000	St	
	M29231	<input type="checkbox"/>	O2 connecting hose 3M	1.000	St	
	M29251	<input type="checkbox"/>	O2 connecting hose 5M	1.000	St	
	M34402	<input type="checkbox"/>	O2-hose NIST 3M DIN probe	1.000	St	
	M34404	<input type="checkbox"/>	ZV-Schlauch N20 1,5M,DIN-ST.	1.000	St	
	M29239	<input type="checkbox"/>	Air-connecting hose 3M	1.000	St	
	M29259	<input type="checkbox"/>	Air-connecting hose 5M	1.000	St	
	M34409	<input type="checkbox"/>	Air-hose NIST 5M DIN probe	1.000	St	
	8414700	<input type="checkbox"/>	Aquapor EL, Humidifier	1.000	St	
	8409716	<input type="checkbox"/>	Bacteria filter EV800/801	1.000	St	
	6870281	<input checked="" type="checkbox"/>	Test filter	1.000	St	
	8412710	<input checked="" type="checkbox"/>	Calibration set CO2	1.000	St	
	6850435	<input checked="" type="checkbox"/>	CO2-cylinder (LINDE) 5%	1.000	St	
	8412840	<input checked="" type="checkbox"/>	Parking support, cpl.	1.000	St	
	8412723	<input checked="" type="checkbox"/>	Tray Evita 4	1.000	St	
	2120046	<input type="checkbox"/>	Resutator 2000 adults	1.000	St	
	5702321	<input checked="" type="checkbox"/>	Resuscitator Child MR_100	1.000	St	
	2120941	<input type="checkbox"/>	Baby-Resutator 2000	1.000	St	
	M26349	<input checked="" type="checkbox"/>	Hook	1.000	St	
	8411950	<input checked="" type="checkbox"/>	Evita mobil, tall	1.000	St	
	8411965	<input checked="" type="checkbox"/>	Evita mobil, short	1.000	St	
	M31796	<input type="checkbox"/>	Cupboard 8H (4 drawers)	1.000	St	
	M31795	<input checked="" type="checkbox"/>	Cupboard 4H (2 drawers)	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

88/107

Parts catalog
Order List




Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
	8409018	<input type="checkbox"/>	Cabinet attachment set Evita	1.000	St	
	8411970	<input checked="" type="checkbox"/>	Set cylinder bracket ev.mobil	1.000	St	
	8411814	<input checked="" type="checkbox"/>	Kit Power socket strip Evita Mobil	1.000	St	
	8413125	<input checked="" type="checkbox"/>	KIT Water trap, for exp. valve	1.000	St	
	8306489	<input checked="" type="checkbox"/>	Printer cable bl.8000	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

89/107

Parts catalog
Accessories Adult Ventilation



Pos. 0 (wenn vorhanden) ohne Abbildung/ Item 0 (if available) without illustration/
Pos. 0 (si existe) sin ilustración/ Pos. 0 (si existent) sans illustration


Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
1	8412723	<input checked="" type="checkbox"/>	Tray Evita 4	1.000	St	
2	8412840	<input checked="" type="checkbox"/>	Parking support, cpl.	1.000	St	
3	8411900	<input type="checkbox"/>	Evita 4, Construction Kit	1.000	St	
4	6850645	<input checked="" type="checkbox"/>	O2-Sensor (capsule)	1.000	St	
5	8412384	<input checked="" type="checkbox"/>	Filter, single	1.000	St	
6	8411950	<input checked="" type="checkbox"/>	Evita mobil, tall	1.000	St	
7	8411965	<input checked="" type="checkbox"/>	Evita mobil, short	1.000	St	
8	M31796	<input type="checkbox"/>	Cupboard 8H (4 drawers)	1.000	St	
9	M31795	<input checked="" type="checkbox"/>	Cupboard 4H (2 drawers)	1.000	St	
10	8411970	<input checked="" type="checkbox"/>	Set cylinder bracket ev.mobil	1.000	St	
12	8411956	<input checked="" type="checkbox"/>	Humidifier bracket Evita mobil	1.000	St	
13	8410580	<input checked="" type="checkbox"/>	Expiration valve, Patientsys.	1.000	St	
14	8403735	<input checked="" type="checkbox"/>	Spirolong Flow Sensor (5x)	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

90/107

Parts catalog
Accessories Adult Ventilation




Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
15	8405371	<input checked="" type="checkbox"/>	Temperature sensor	1.000	St	
16	8405435	<input checked="" type="checkbox"/>	Y-piece evita	1.000	St	
17	8403201	<input checked="" type="checkbox"/>	Test lung	1.000	St	
18	6870279	<input checked="" type="checkbox"/>	CO2-Cuvette ERW	1.000	St	
19	M23840	<input checked="" type="checkbox"/>	Standard conn.12,5 plastic	1.000	St	
20	8402041	<input checked="" type="checkbox"/>	Corrugated hose	1.000	St	
21	8401644	<input checked="" type="checkbox"/>	Cap (11 mm)	1.000	St	
22	8403076	<input checked="" type="checkbox"/>	Adapter E	1.000	St	
23	8403685	<input checked="" type="checkbox"/>	Set 12 stand. cone a plastic	1.000	St	
24	6871500	<input checked="" type="checkbox"/>	Mainstream sensor	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

91/107

Parts catalog
Accessories Adult Ventilation




Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
25	2165619	<input checked="" type="checkbox"/>	Breathing hose 22/22 (A) 35	2.000	St	
26	M25647	<input checked="" type="checkbox"/>	Double adapter 22/22	1.000	St	
27	M25649	<input checked="" type="checkbox"/>	Maskenkrummer	1.000	St	
28	8414700	<input type="checkbox"/>	Aquapor EL, Humidifier	1.000	St	
29	8403345	<input checked="" type="checkbox"/>	Set of spare brackets	1.000	St	
30	8405029	<input checked="" type="checkbox"/>	Patient part (Aquapor)	1.000	St	
31	8403566	<input checked="" type="checkbox"/>	Tube clamps	1.000	St	
32	8403976	<input checked="" type="checkbox"/>	Pot	1.000	St	
33	8404985	<input checked="" type="checkbox"/>	Water trap	1.000	St	
34	2165627	<input checked="" type="checkbox"/>	Breathing hose 22/22 (A) 60	3.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

92/107

Parts catalog
Accessories Adult Ventilation




Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
35	8409746	<input checked="" type="checkbox"/>	Support for hinged arm A	1.000	St	
36	8409841	<input checked="" type="checkbox"/>	Hose clamp for hinged arm	1.000	St	
37	8409609	<input checked="" type="checkbox"/>	Hinged arm	1.000	St	
38	2M85706	<input checked="" type="checkbox"/>	Quickstop hinged arm 2	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

93/107

Parts catalog
Accessories Adult Ventilation




Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
16	8405435	<input checked="" type="checkbox"/>	Y-piece evita	1.000	St	
19	M23840	<input checked="" type="checkbox"/>	Standard conn.12,5 plastic	1.000	St	
20	8402041	<input checked="" type="checkbox"/>	Corrugated hose	1.000	St	
21	8401644	<input checked="" type="checkbox"/>	Cap (11 mm)	1.000	St	
22	8403076	<input checked="" type="checkbox"/>	Adapter E	1.000	St	
23	8403685	<input checked="" type="checkbox"/>	Set 12 stand. cone a plastic	1.000	St	
25	2165619	<input checked="" type="checkbox"/>	Breathing hose 22/22 (A) 35	2.000	St	
26	M25647	<input checked="" type="checkbox"/>	Double adapter 22/22	1.000	St	
27	M25649	<input checked="" type="checkbox"/>	Maskenkrummer	1.000	St	
31	8403566	<input checked="" type="checkbox"/>	Tube clamps	1.000	St	
33	8404985	<input checked="" type="checkbox"/>	Water trap	1.000	St	
34	2165627	<input checked="" type="checkbox"/>	Breathing hose 22/22 (A) 60	3.000	St	
39	8412092	<input checked="" type="checkbox"/>	Hose kit adult, Aquapor	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

94/107

Parts catalog
Accessories Paed. Ventilation



Pos. 0 (wenn vorhanden) ohne Abbildung/ Item 0 (if available) without illustration/
Pos. 0 (si existe) sin ilustración/ Pos. 0 (si existent) sans illustration


Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
1	8412723	<input checked="" type="checkbox"/>	Tray Evita 4	1.000	St	
2	8412840	<input checked="" type="checkbox"/>	Parking support, cpl.	1.000	St	
3	8411900	<input type="checkbox"/>	Evita 4, Construction Kit	1.000	St	
4	6850645	<input checked="" type="checkbox"/>	O2-Sensor (capsule)	1.000	St	
5	8412384	<input checked="" type="checkbox"/>	Filter, single	1.000	St	
6	8411950	<input checked="" type="checkbox"/>	Evita mobil, tall	1.000	St	
7	8411965	<input checked="" type="checkbox"/>	Evita mobil, short	1.000	St	
8	M31796	<input type="checkbox"/>	Cupboard 8H (4 drawers)	1.000	St	
9	M31795	<input checked="" type="checkbox"/>	Cupboard 4H (2 drawers)	1.000	St	
10	8411970	<input checked="" type="checkbox"/>	Set cylinder bracket ev.mobil	1.000	St	
12	8411956	<input checked="" type="checkbox"/>	Humidifier bracket Evita mobil	1.000	St	
13	8410580	<input checked="" type="checkbox"/>	Expiration valve, Patientsys.	1.000	St	
14	8403735	<input checked="" type="checkbox"/>	Spirolog Flow Sensor (5x)	1.000	St	
15	6871500	<input checked="" type="checkbox"/>	Mainstream sensor	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

95/107

Parts catalog
Accessories Paed. Ventilation




Pos. 0 (wenn vorhanden) ohne Abbildung/ Item 0 (if available) without illustration/
Pos. 0 (si existe) sin ilustración/ Pos. 0 (si existant) sans illustration

Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
16	8409716	<input type="checkbox"/>	Bacteria filter EV800/801	1.000	St	
17	2165856	<input checked="" type="checkbox"/>	Breathing hose 22/10 (N/P)40	1.000	St	
18	M19351	<input checked="" type="checkbox"/>	Catheter 11	1.000	St	
19	8411048	<input checked="" type="checkbox"/>	Dual airway temp.sensor 1,45M	1.000	St	
20	8411074	<input checked="" type="checkbox"/>	Mounting kit	1.000	St	
23	8411047	<input checked="" type="checkbox"/>	Humidifier chamber MR 340 (E)	1.000	St	
24	8411050	<input checked="" type="checkbox"/>	Draw wire	1.000	St	
25	8411045	<input checked="" type="checkbox"/>	Heating 1,10m	1.000	St	
26	2165651	<input checked="" type="checkbox"/>	Breathing hose 22/10 (N/P) 110	1.000	St	
27	8411044	<input checked="" type="checkbox"/>	Probe-thermometer housing drae	1.000	St	
28	8409634	<input checked="" type="checkbox"/>	accordion hose flex.	1.000	St	
29	8401645	<input checked="" type="checkbox"/>	Cap (7 mm)	1.000	St	
30	8403075	<input checked="" type="checkbox"/>	Y-piece 90 degree	1.000	St	
31	6870280	<input checked="" type="checkbox"/>	CO2-Cuvette Paed.	1.000	St	
32	8409742	<input checked="" type="checkbox"/>	Test lung, cpl.	1.000	St	
33	8409897	<input checked="" type="checkbox"/>	Double cone 11A	1.000	St	
34	2165848	<input checked="" type="checkbox"/>	Breathing hose 10/10 (N/P) 60	1.000	St	
35	8403976	<input checked="" type="checkbox"/>	Pot	1.000	St	
36	8409627	<input checked="" type="checkbox"/>	Condensate trap ex.	1.000	St	
37	2165821	<input checked="" type="checkbox"/>	Breathing hose 22/10 (N/P) 60	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

96/107

Parts catalog
Accessories Paed. Ventilation


Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
38	8409746	<input checked="" type="checkbox"/>	Support for hinged arm A	1.000	St	
39	8409841	<input checked="" type="checkbox"/>	Hose clamp for hinged arm	1.000	St	
40	8409609	<input checked="" type="checkbox"/>	Hinged arm	1.000	St	
41	2M85706	<input checked="" type="checkbox"/>	Quickstop hinged arm 2	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

97/107

Parts catalog
Accessories Paed. Ventilation




Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
17	2165856	<input checked="" type="checkbox"/>	Breathing hose 22/10 (N/P) 40	1.000	St	
18	M19351	<input checked="" type="checkbox"/>	Catheter 11	1.000	St	
26	2165651	<input checked="" type="checkbox"/>	Breathing hose 22/10 (N/P) 110	1.000	St	
27	8411044	<input checked="" type="checkbox"/>	Probe-thermometer housing drae	1.000	St	
28	8409634	<input checked="" type="checkbox"/>	accordion hose flex.	1.000	St	
29	8401645	<input checked="" type="checkbox"/>	Cap (7 mm)	1.000	St	
30	8403075	<input checked="" type="checkbox"/>	Y-piece 90 degree	1.000	St	
32	8409742	<input checked="" type="checkbox"/>	Test lung, cpl.	1.000	St	
33	8409897	<input checked="" type="checkbox"/>	Double cone 11A	1.000	St	
34	2165848	<input checked="" type="checkbox"/>	Breathing hose 10/10 (N/P) 60	1.000	St	
36	8409627	<input checked="" type="checkbox"/>	Condensate trap ex.	1.000	St	
37	2165821	<input checked="" type="checkbox"/>	Breathing hose 22/10 (N/P) 60	1.000	St	
42	8412081	<input checked="" type="checkbox"/>	Hose kit ped., F&P	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

98/107

Parts catalog
Calibration Set CO2



Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
1	8290272	<input checked="" type="checkbox"/>	Valve for CO2-monitor	1.000	St	
2	1180614	<input checked="" type="checkbox"/>	Hose 2X1-SI NF M17749	1.000	m	
3	1198343	<input checked="" type="checkbox"/>	Hose 7X2,5 SI NF M29908	1.000	m	
4	M07582	<input checked="" type="checkbox"/>	Secretion control glass	1.000	St	
5	6850435	<input checked="" type="checkbox"/>	CO2-cylinder (LINDE) 5%	1.000	St	
6	8412710	<input checked="" type="checkbox"/>	Calibration set CO2	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

99/107

Parts catalog
Aquapor EL

Zum Öffnen des Teilekatalogs hier doppelklicken.
Double-click here to open the parts catalog.
Hacer dobleclick aquí para abrir el catálogo de piezas.
Double-cliquer ici pour ouvrir le catalogue des pièces.


Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
----------	-----------	----------	-------------	------	----------	--------

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

100/107

Parts catalog
CP parts



Pos. 0 ohne Abbildung/ Item 0 without illustration/ Pos. 0 sin ilustración/ Pos. 0 sans illustration

Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
1	8414986	<input checked="" type="checkbox"/>	BlueSet Heated (P/A)	1.000	St	
2	8414987	<input checked="" type="checkbox"/>	BlueSet Heated (N)	1.000	St	
3	6733895	<input checked="" type="checkbox"/>	Set Mic.filter 654ST-ISOCCLICK	1.000	St	
4	6850645	<input checked="" type="checkbox"/>	O2-Sensor (capsule)	1.000	St	
5	8403735	<input checked="" type="checkbox"/>	Spirolog Flow Sensor (5x)	1.000	St	
6	8414776	<input checked="" type="checkbox"/>	Set single use valve, 10 pcs.	1.000	St	
7	MP01579	<input checked="" type="checkbox"/>	NovaStar® TS NIV-Mask,SE,S	1.000	St	
8	MP01580	<input checked="" type="checkbox"/>	NovaStar® TS NIV-Mask,SE,M	1.000	St	
9	MP01581	<input checked="" type="checkbox"/>	NovaStar® TS NIV-Mask,SE,L	1.000	St	
10	MP01573	<input checked="" type="checkbox"/>	Mask ClassicStar,NIV,SE,S	1.000	St	
11	MP01574	<input checked="" type="checkbox"/>	Mask ClassicStar,NIV,SE,M	1.000	St	
12	MP01575	<input checked="" type="checkbox"/>	Mask ClassicStar,NIV,SE,L	1.000	St	
13	MP01623	<input checked="" type="checkbox"/>	Mask ClassicStar Nasal NV, S	1.000	St	
14	MP01624	<input checked="" type="checkbox"/>	Mask ClassicStar Nasal NV, M	1.000	St	
15	MP01625	<input checked="" type="checkbox"/>	Mask ClassicStar Nasal NV, L	1.000	St	
16	MP01063	<input checked="" type="checkbox"/>	Disposable CO2 Cuvette paed.	1.000	St	
17	MP01062	<input checked="" type="checkbox"/>	Disposable CO2 Cuvette Adult	1.000	St	
18	8414991	<input checked="" type="checkbox"/>	Hose kit for MR850, RT 212	1.000	St	
19	MP01650	<input type="checkbox"/>	Evita VPack w ClassicS Fullf S	1.000	St	
20	MP01651	<input type="checkbox"/>	Evita VPack w ClassicS Fullf M	1.000	St	
21	MP01652	<input type="checkbox"/>	Evita VPack w ClassicS Fullf L	1.000	St	
22	MP01653	<input type="checkbox"/>	Evita VPack w NovaS Fullf S	1.000	St	
23	MP01654	<input type="checkbox"/>	Evita VPack w NovaS Fullf M	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

101/107

Parts catalog
CP parts



Pos. 0 ohne Abbildung/ Item 0 without illustration/ Pos. 0 sin ilustración/ Pos. 0 sans illustration


Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
24	MP01655	<input type="checkbox"/>	Evita VPack w NovaS Fullf L	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

102/107

Parts catalog
CP parts



Pos. 0 ohne Abbildung/ Item 0 without illustration/ Pos. 0 sin ilustración/ Pos. 0 sans illustration


Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
25	MP00306	<input checked="" type="checkbox"/>	VentStar dual heated	1.000	St	
26	MP00307	<input checked="" type="checkbox"/>	VentStar heated	1.000	St	
27	MP00308	<input checked="" type="checkbox"/>	VentStar heated (N)	1.000	St	
28	MP00309	<input checked="" type="checkbox"/>	VentStar heated (N) basic	1.000	St	
29	MP00305	<input checked="" type="checkbox"/>	Vent Set Flex	1.000	St	
30	MP00315	<input checked="" type="checkbox"/>	Vent Set COAX	1.000	St	
31	MP00316	<input checked="" type="checkbox"/>	Vent Set COAX with Holder	1.000	St	
32	MP00338	<input checked="" type="checkbox"/>	VentStar Basic P 150 w/o LuerL	1.000	St	
33	MP00349	<input checked="" type="checkbox"/>	VentStar Basic 180	1.000	St	
34	MP00350	<input checked="" type="checkbox"/>	VentStar Basic 250	1.000	St	
35	MP00351	<input checked="" type="checkbox"/>	VentStar Basic (P) 180	1.000	St	
36	MP00355	<input checked="" type="checkbox"/>	VentStar Flex 220	1.000	St	
37	MP00356	<input checked="" type="checkbox"/>	VentStar Coax 180	1.000	St	
38	MP00361	<input checked="" type="checkbox"/>	VentStar Watertrap (P) 180	1.000	St	
39	MP00362	<input checked="" type="checkbox"/>	VentStar Watertrap (P) 180 w/oLL	1.000	St	
40	MP00363	<input checked="" type="checkbox"/>	VentStar Watertrap (N) 180	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

103/107

Parts catalog
CP parts



Pos. 0 ohne Abbildung/ Item 0 without illustration/ Pos. 0 sin ilustración/ Pos. 0 sans illustration


Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
41	MP00234	<input checked="" type="checkbox"/>	Disp. humi.chamber f.F&P MR850	1.000	St	
42	MP01590	<input checked="" type="checkbox"/>	2side pumpball f.airmanagement	1.000	St	
43	MP00334	<input checked="" type="checkbox"/>	VentStar Watert. 180 w/o LuerL	1.000	St	
44	MP00337	<input checked="" type="checkbox"/>	VentStar Watertrap 180	1.000	St	
45	MP01855	<input checked="" type="checkbox"/>	ErgoStar CM 55	1.000	St	
46	MP01840	<input checked="" type="checkbox"/>	ErgoStar CM 40	1.000	St	
47	MP01845	<input checked="" type="checkbox"/>	ErgoStar CM 45	1.000	St	
48	MP01850	<input checked="" type="checkbox"/>	ErgoStar CM 50	1.000	St	
49	MP01890	<input checked="" type="checkbox"/>	ErgoStar AC 90	1.000	St	
50	MP01895	<input checked="" type="checkbox"/>	ErgoStar AC 95	1.000	St	
51	MP01860	<input checked="" type="checkbox"/>	ErgoStar CM 60	1.000	St	
52	8418435	<input type="checkbox"/>	BabyFlow Starterp. reusable	1.000	St	
53	8418410	<input type="checkbox"/>	BabyFlow reusable	1.000	St	
54	8418414	<input type="checkbox"/>	Silicon hoses with connectors	1.000	St	
55	8418411	<input type="checkbox"/>	Prongadapter	1.000	St	
56	8418540	<input type="checkbox"/>	Headband XXL plus	1.000	St	
57	8418423	<input type="checkbox"/>	Headband XXL (extra ex. large)	1.000	St	
58	8418422	<input type="checkbox"/>	Headband XL (extra large)	1.000	St	
59	8418418	<input type="checkbox"/>	Stirnband XS (extra klein)	1.000	St	
60	8418419	<input type="checkbox"/>	Headband S (small)	1.000	St	
61	8418420	<input type="checkbox"/>	Headband M (medium)	1.000	St	
62	8418421	<input type="checkbox"/>	Headband L (large)	1.000	St	
63	MP01755	<input checked="" type="checkbox"/>	Filter CareStar 45	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

104/107

Parts catalog
CP parts



Pos. 0 ohne Abbildung/ Item 0 without illustration/ Pos. 0 sin ilustración/ Pos. 0 sans illustration


Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
64	MP01770	<input checked="" type="checkbox"/>	Filter CareStar 30	1.000	St	
65	MP01765	<input checked="" type="checkbox"/>	Filter CareStar 40A	1.000	St	
66	MP01785	<input checked="" type="checkbox"/>	Filter SafeStar 80	1.000	St	
67	MP01790	<input checked="" type="checkbox"/>	Filter SafeStar 55	1.000	St	
68	MP01795	<input checked="" type="checkbox"/>	Filter SafeStar 60A	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

105/107

Parts catalog
CP parts



Pos. 0 ohne Abbildung/ Item 0 without illustration/ Pos. 0 sin ilustración/ Pos. 0 sans illustration


Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
69	MP01810	<input checked="" type="checkbox"/>	Filter/HME TwinStar 65A	1.000	St	
70	MP01800	<input checked="" type="checkbox"/>	Filter/HME TwinStar 90	1.000	St	
71	MP01825	<input checked="" type="checkbox"/>	Filter/HME TwinStar 10A	1.000	St	
72	MP01820	<input checked="" type="checkbox"/>	Filter/HME TwinStar 8	1.000	St	
73	MP01815	<input checked="" type="checkbox"/>	Filter/HME TwinStar 25	1.000	St	
74	MP01801	<input checked="" type="checkbox"/>	Filter/HME TwinStar HEPA	1.000	St	
75	MP01805	<input checked="" type="checkbox"/>	Filter/HME TwinStar 55	1.000	St	
76	8418550	<input checked="" type="checkbox"/>	BabyFlow Demop. disposable	1.000	St	
77	8418583	<input checked="" type="checkbox"/>	BabyFlow disposable (20x)	1.000	St	
78	8418415	<input checked="" type="checkbox"/>	Prong XS (10pcs)	1.000	St	
79	8418605	<input checked="" type="checkbox"/>	Prong S (10pcs)	1.000	St	
80	8418416	<input checked="" type="checkbox"/>	Prong M (10pcs)	1.000	St	
81	8418531	<input checked="" type="checkbox"/>	Prong L (10pcs)	1.000	St	
82	8418417	<input checked="" type="checkbox"/>	Prong XL (10pcs)	1.000	St	
83	8418491	<input checked="" type="checkbox"/>	NeoMask S (10pcs)	1.000	St	
84	8418490	<input checked="" type="checkbox"/>	NeoMask M (10pcs)	1.000	St	
85	8418619	<input checked="" type="checkbox"/>	NeoMask L (10pcs)	1.000	St	
86	8418534	<input checked="" type="checkbox"/>	Cap S (smal 5pcs)	1.000	St	
87	8418535	<input checked="" type="checkbox"/>	Cap M (medium 5pcs)	1.000	St	
88	8418536	<input checked="" type="checkbox"/>	Cap L (large 5pcs)	1.000	St	
89	8418537	<input checked="" type="checkbox"/>	Cap XL (extra large 5 pcs)	1.000	St	
90	8418538	<input checked="" type="checkbox"/>	cap XXL(ex. ex. large 5pcs)	1.000	St	
91	8418539	<input checked="" type="checkbox"/>	Cap XXL plus (5pcs)	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

106/107

Parts catalog
CP parts



Pos. 0 ohne Abbildung/ Item 0 without illustration/ Pos. 0 sin ilustración/ Pos. 0 sans illustration

Item No.	Order No.	Oderable	Description	Qty.	Qty.unit	Remark
92	MP02400	<input checked="" type="checkbox"/>	SelfTestLung	1.000	St	
93	8412935	<input checked="" type="checkbox"/>	Pneum. Medication Nebulizer	1.000	St	
94	8411030	<input checked="" type="checkbox"/>	Pneum. Medication Nebulizer N	1.000	St	
95	8415978	<input checked="" type="checkbox"/>	Dust filter	1.000	St	

Items that are shown in the illustration but are not listed below the illustration are not available as spare parts

Evita 4/ Evita 4 edition
Revision: 22

107/107

5 Test instructions

Test Instructions

This chapter contains the measures required to determine the actual condition of the device.

5.1 Test Instructions / Service Card IPM

Dräger

Test Instructions / Service Card IPM

Evita 4 / Evita 4 edition



Warning

All servicing and/or test procedures on the device require detailed knowledge of this documentation. Use of the device requires detailed knowledge and observance of the relevant Instructions for Use.

Revision 13.0

Table of contents

Important notes	5
1 Device configuration	5
1.1 Device configuration	6
1.1.1 Serial number (if not otherwise recorded)	6
1.1.2 Software version	6
1.1.3 Operating hours	7
1.1.4 Device ID	7
2 Maintenance parts	7
2.1 Maintenance intervals, overview	7
2.1.1 Maintenance intervals and required sets and parts	7
2.2 Maintenance parts by specified interval	8
2.2.1 Service set, once a year, MX08220	8
2.2.2 Option for Canada and USA, filter element MP03903	8
2.2.3 Gas inlet block	8
2.2.3.1 Service set, gas inlet block FAS, MX08225	9
2.2.3.2 Service set, gas inlet block Dräger, MX08226	10
2.2.3.3 Service set, gas inlet block FAS, MX08227	11
2.2.3.4 Service set, gas inlet block Dräger, MX08228	12
2.2.4 Optional rechargeable battery (DC module), 8421988	13
2.2.5 Optional rechargeable battery (trolley), 8422125	13
2.2.6 Replacing the LCD backlight after 25000 operating hours	14
2.2.6.1 Replacement of no longer available displays and backlights	14
2.2.6.2 Maintenance-free display	15
2.2.7 Maintenance parts as required	16
2.2.7.1 Maintenance parts as per instructions for use	16
3 Electrical safety	17
3.1 Electrical safety to DIN EN 62353 (IEC 62353)	17
3.1.1 Visual check	17
3.1.2 Protective earth resistance	18
3.1.3 Protective earth resistance measuring points	19
3.1.4 Test adapters for equipment leakage current	19
3.1.5 Equipment leakage current	19
3.1.6 Test adapters for patient leakage current	21
3.1.7 Leakage current, mains on applied part with airway temperature sensor AWT 01, if installed	21
3.1.8 Leakage current on applied part with test adapter (normal condition)	22
3.1.9 Optional SpO2, leakage current, mains on applied part	24
3.2 Electrical safety according to IEC 60601-1	25
3.2.1 Visual check	26
3.2.2 Protective earth resistance	27
3.2.3 Protective earth resistance measuring points	27
3.2.4 Earth leakage current	28
3.2.5 Patient leakage current	28
4 Function and condition test	30

4.1	Condition of basic device.....	30
4.1.1	Labels and instructions for use.....	30
4.1.2	Base unit	31
4.1.3	Trolley.....	31
4.1.4	Compressed gas connecting hoses	31
4.2	Function check	31
4.2.1	PEEP/PIP valve calibration	31
4.2.2	Power-on test	33
4.2.3	Test according to internal device checklist as per instructions for use.....	33
4.2.4	Non-return valve (D5.1) in expiratory valve	34
4.2.5	Non-return valve D3.2 (10 mbar).....	35
4.2.6	Emergency air valve Y3.1/non-return valve D3.1	35
4.2.7	Safety valve D3.3 (120 mbar).....	37
4.2.8	DC power supply unit with LEDs	38
4.2.9	DC module (if present)	39
4.3	Testing ventilation in adult mode	39
4.3.1	Calibrating the sensors.....	41
4.3.2	O ₂ comparative measurement	41
4.3.3	Volume metering/flow measurement.....	42
4.3.4	BIPAP/PCV+ ventilation mode	43
4.4	Additional function tests.....	44
4.4.1	Power failure and data backup (ventilators without DC module).....	44
4.5	DC module or Benning power supply unit (if fitted), additional tests	44
4.5.1	Internal rechargeable battery	44
4.5.2	Mains power failure and data backup.....	45
4.5.3	"External rechargeable battery" option	45
4.6	CO ₂ sensor option.....	45
4.6.1	Condition	45
4.6.2	Measurement accuracy with test filter	46
4.7	Options	46
4.7.1	Evita SAT (SpO ₂ , hardware).....	46
4.7.2	NeoFlow (hardware).....	47
4.7.3	Breathing Support Package (software)	47
4.7.4	ATC (Automatic Tubus Compensation, software)	47
4.7.5	CapnoPlus (hardware)	48
4.7.6	SW 4.00 Plus (software).....	48
4.7.7	NIV (Non-Invasive Ventilation, software).....	48
4.7.8	IFCO Carrier/NurseCall (hardware).....	48
4.7.9	Remote Pad (hardware)	49
4.7.10	Communication between ventilator and monitor (hardware).....	49
4.8	Final action	50
4.8.1	Device handover	50
5	Test equipment	50
5.1	Test equipment.....	50
5.1.1	Test equipment subject to mandatory calibration.....	50
5.1.2	Test equipment not subject to mandatory calibration.....	51

Evita 4 / Evita 4 edition
Table of contents

5.1.3	Additionally required items	54
6	Annex	55
6.1	Evita 4 service mode start screen.....	55
6.1.1	Access	55

No.0306_0000007424

Important notes

NOTE

Do not use these test instructions for testing after a repair procedure.

NOTE

Prior to using these test instructions, check that they are the latest revision (compare revision with latest service documentation).

All results and inputs must be documented in the "Test Report" and "Result Sheet".

NOTE

Clean and disinfect the device and device components prior to each maintenance procedure.

NOTE

For the function test the rechargeable batteries must be fully charged.

These test instructions apply from software version 4.n.

For the devices listed below, use the test instructions specific to the relevant device:

- Respiratory air compressor
- Humidifier



No.0306_0000007424

1 Device configuration

This section records the device configuration.

1.1 Device configuration

The basic device configuration comprises the following assemblies:

- Base unit
- Expiratory valve
- Ventilation hose system
- Trolley
- Compressed gas connecting hoses

1.1.1 Serial number (if not otherwise recorded)

Action • Read the serial number on the basic device.

Result **Serial number of the basic device**

[_____]txt]

Action • Document the expiratory valve used.

Result **Document the use of a disposable expiratory valve with "Yes" or "No".**

[_____]txt]

NOTE

Disposable expiratory valves have no serial number.

Action • Read the serial number of the expiratory valve(s).

Result **Serial number of expiratory valve 1**

[_____]txt]

Result **Serial number of expiratory valve 2**

[_____]txt]

Action • Read the serial number on the control unit of the ventilator.

Result **Serial number of the control unit**

[_____]txt]

Action • Read the serial number on the CO₂ sensor (option).

Result **Serial number of the CO₂ sensor**

[_____]txt]

1.1.2 Software version

- Action • Connect unit to power supply.
- Switch on the device.
 - Press the "Standby" key.
 - Press "Configuration" key.
 - Read off the software version shown on the display.

Result **Software version**

[_____]txt]

No.0306_000007424

1.1.3 Operating hours

- Action
- as described under [1.1.2 Software version](#).
 - Read off the operating hours shown on the display.

Result **Operating hours**

[_____ h]

1.1.4 Device ID

- Action
- as described under [1.1.2 Software version](#).
 - Read off the device ID shown on the display.

Result **Device ID**

[_____ txt]

2 Maintenance parts

This chapter contains interval-related maintenance parts, measures, and tests that can only be performed on an open device.

2.1 Maintenance intervals, overview**2.1.1 Maintenance intervals and required sets and parts**

The following table presents an overview of the maintenance procedures and the required sets and maintenance parts over 78 months. After 78 months the table is worked through again starting from the first column.

Designation	Part number	Service interval in months										
		12	18	24	30	36	42	48	54	60	66	72
Service set, every 1 years	MX08220	x		x		x		x		x		x
Option: Filters for Canada and USA	MP03903	x		x		x		x		x		x
Devices with gas supply block FAS or Devices with gas supply block Dräger	MX08225 or MX08226			x				x				
Devices with gas supply block FAS or Devices with gas supply block Dräger	MX08227 or MX08228											x
Set 2 batteries Evita internal	8421988			x				x				x
Set 2 batteries Evita external	8422125			x				x				x

2.2 Maintenance parts by specified interval

2.2.1 Service set, once a year, MX08220

NOTE

When using disposable expiratory valves, the diaphragm 8410181 and the sealing washer 8407979 are not required.

The service set (quantity 1) with the number MX08220 includes the following items:

Quantity	Designation	Number	Location/Remark
1	Diaphragm	8410181	Expiratory valve / Replacement by specialist
1	Sealing washer	8407979	Expiratory valve / Replacement by specialist

Result **Service set, once a year, MX08220**

Next replacement: [_____] dat

2.2.2 Option for Canada and USA, filter element MP03903

or if the above filters are no longer available:

Quantity	Designation	Number	Location/Remark
1	AIR filter	MP03903	Gas inlet block/Replacement by specialist

Result **AIR filter**

Next replacement: [_____] dat

Quantity	Designation	Number	Location/Remark
1	O ₂ filter	MP03903	Gas inlet block/Replacement by specialist

Result **O₂ filter**

Next replacement: [_____] dat

2.2.3 Gas inlet block

NOTE

Different maintenance parts are required depending on the gas inlet block used (FAS or Dräger), as explained below.

2.2.3.1 Service set, gas inlet block FAS, MX08225**Fig. 1** Rear of device; FAS gas inlet block

The service set (quantity 1) with the number MX08225 for Evita with FAS gas inlet block includes the following items:

Quantity	Designation	Number	Location/Remark
1	Diaphragm	8410181	Expiratory valve / Replacement by specialist
1	Sealing washer	8407979	Expiratory valve / Replacement by specialist
1	Battery, lithium, 3 V / 1400 mAh	1835343	CPU / Replacement by specialist
2	Filter	8408208	Gas inlet block FAS / Replacement by specialist
2	O-ring	8411516	Gas inlet block FAS / Replacement by specialist

NOTE

When using disposable expiratory valves, the diaphragm 8410181 and the sealing washer 8407979 are not required.

Result **Devices with FAS gas inlet block**

Next replacement: [_____]dat]

Evita 4 / Evita 4 edition
Maintenance parts

2.2.3.2 Service set, gas inlet block Dräger, MX08226




Fig. 2 Rear of device; Dräger gas inlet block

The service set (quantity 1) with the number MX08226 for Evita with Dräger gas inlet block includes the following items:

Quantity	Designation	Number	Location/Remark
1	Diaphragm	8410181	Expiratory valve / Replacement by specialist
1	Sealing washer	8407979	Expiratory valve / Replacement by specialist
1	Battery, lithium, 3 V / 1400 mAh	1835343	CPU / Replacement by specialist
2	Filter	8408208	Gas inlet block Dräger / Replacement by specialist
1	Flat seal	8408204	Gas inlet block Dräger / Replacement by specialist


NOTE

When using disposable expiratory valves, the diaphragm 8410181 and the sealing washer 8407979 are not required.

Result **Devices with Dräger gas inlet block**

Next replacement: [_____] dat]

No.0306_000007424

2.2.3.3 Service set, gas inlet block FAS, MX08227**Fig. 3** Rear of device; FAS gas inlet block

The 6-yearly service set (quantity 1) with the number MX08227 for Evita with FAS gas inlet block includes the following items:

Quantity	Designation	Number	Location/Remark
1	Diaphragm	8410181	Expiratory valve / Replacement by specialist
1	Sealing washer	8407979	Expiratory valve / Replacement by specialist
1	Battery, lithium, 3 V / 1400 mAh	1835343	CPU / Replacement by specialist
2	Filter	8408208	Gas inlet block FAS / Replacement by specialist
2	O-ring	8411516	Gas inlet block FAS / Replacement by specialist
2	Diaphragm	8411513	Gas inlet block FAS / Replacement by specialist
2	Valve tappet	8411514	Gas inlet block FAS / Replacement by specialist
1	Real-time clock DIL24	1837087	CPU / Replacement by specialist
1	Fan	8306560	Pneumatic assembly/Replacement by specialist
1	Spirolog sensor cable harness	8415709	Pneumatic assembly/Replacement by specialist
1	SV diaphragm assembly	8414179	Safety valve / Replacement by specialist

NOTE

When using disposable expiratory valves, the diaphragm 8410181 and the sealing washer 8407979 are not required.

Result **Devices with FAS gas inlet block**

Next replacement: [_____] dat]

No.0306_0000007424

Evita 4 / Evita 4 edition
Maintenance parts

2.2.3.4 Service set, gas inlet block Dräger, MX08228




Fig. 4 Rear of device; Dräger gas inlet block

The 6-yearly service set (quantity 1) with the number MX08228 for Evita with Dräger gas inlet block includes the following items:

Quantity	Designation	Number	Location/Remark
1	Diaphragm	8410181	Expiratory valve / Replacement by specialist
1	Sealing washer	8407979	Expiratory valve / Replacement by specialist
1	Battery, lithium, 3 V / 1400 mAh	1835343	CPU / Replacement by specialist
2	Filter	8408208	Gas inlet block Dräger / Replacement by specialist
1	Flat seal	8408204	Gas inlet block Dräger / Replacement by specialist
2	Pressure regulator	8408205	Gas inlet block Dräger / Replacement by specialist
1	Real-time clock DIL24	1837087	CPU / Replacement by specialist
1	Fan	8306560	Pneumatic assembly/Replacement by specialist
1	Spirolog sensor cable harness	8415709	Pneumatic assembly/Replacement by specialist
1	SV diaphragm assembly	8414179	Safety valve / Replacement by specialist

NOTE

When using disposable expiratory valves, the diaphragm 8410181 and the sealing washer 8407979 are not required.

Result **Devices with Dräger gas inlet block**

Next replacement: [_____] dat]

No.0306_000007424

2.2.4 Optional rechargeable battery (DC module), 8421988**CAUTION****Risk of overcharging of internal rechargeable batteries!**

- ▶ When using power supply unit 8419059 (8306520+8306530), the "Battery PCB upgrade kit type B" must always be installed (8419275). When the PCB is installed, the label +/- REV B should be fixed from the outside on the top of the screen frame.
- ▶ When using power supply unit 8419616 (identifiable by the LEDs next to the power switch), the "Battery PCB upgrade kit type B" (8419275) and the label are not needed.

**Fig. 5** Power Supply 8419275**"Panasonic" rechargeable battery (power supply unit)**

Quantity	Designation	Number	Location/Remark
1	Set 2 batteries Evita internal	8421988	DC module (power supply unit) / Replacement by specialist

Result **Optional rechargeable battery (DC module)**

Next replacement: [_____] dat]

2.2.5 Optional rechargeable battery (trolley), 8422125

Quantity	Designation	Number	Location/Remark
1	Set 2 batteries Evita external	8422125	Trolley

Result **Optional rechargeable battery (trolley)**

Next replacement: [_____] dat]

Evita 4 / Evita 4 edition
Maintenance parts

2.2.6 Replacing the LCD backlight after 25000 operating hours

NOTE

Display 8411841 is always indicated in Service mode regardless of the display installed.

Different displays have been used over the years. There is often a label on the rear of the control panel identifying the display. If there is no label on the rear of the display, the control panel must be opened. The type designation is given on the rear of the display.

2.2.6.1 Replacement of no longer available displays and backlights

Display Toshiba type 1, if fitted

- The display was in production in 1995.
- The display type designation is LTM10C015K.
- The LCD backlight is no longer available.
- The **replacement** is the repair kit Display AUO (Evita 4).

Display NEC, if fitted

- The display was in production in 1996.
- The LCD backlight is no longer available.
- The **replacement** is the repair kit Display AUO (Evita 4).

Display Sharp, if fitted

- The display was in production in 1996.
- The LCD backlight is no longer available.
- The **replacement** is the repair kit Display AUO (Evita 4).

Display Toshiba 8411841, if fitted

The display was in production from 1997 to 2007, and is identifiable by the label under the control unit rating plate. The label indicates the number 8411841. If the label is not present, open the control unit. On the rear the Toshiba display type LTM 10C209A is specified.

The LCD backlight is no longer available. The **replacement** is the repair kit Display AUO (Evita 4).

Repair kit Display AUO (Evita 4)

The repair kit Display AUO (Evita 4) replaces all the displays and backlights listed above. The repair kit Display AUO (Evita 4) includes the maintenance-free Toshiba 8418948 display cited in the following test item.

Quantity	Designation	Number	Location/Remark
1	Repair kit Display AUO (Evita 4)	8418948	Front panel / Replacement by specialist

Result **Date of installing repair kit 8420794**

[txt]

No.0306_000007424

2.2.6.2 Maintenance-free display

Display Toshiba 8418948, if fitted

The display entered production as from mid 2011, or has been retrofitted with repair kit Display AUO (Evita 4). The LCD backlight is maintenance-free.

The display is identifiable by the label beneath the type plate on the control unit (Fig. 7). The label indicates the number "8418948".



Fig. 6 Labels

If the label is not present, open the control unit. On the rear of the fixing plate the display type "FPDS-104KID" is specified.

Display Toshiba 8417028, if fitted

The display has been in production since 2008. The LCD backlight is maintenance-free.

The display is identifiable by the label beneath the type plate on the control unit (Fig. 7). The label indicates the number "8417028".

Evita 4 / Evita 4 edition
Maintenance parts



Fig. 7 Labels

If the label is not present, open the control unit. On the rear the Toshiba display type "LTA 104A261F" is specified.

Result **No maintenance required**

[☐ OK]

2.2.7 Maintenance parts as required

2.2.7.1 Maintenance parts as per instructions for use

NOTE

Replacement of maintenance parts is detailed in the instructions for use, and is carried out by the users.

Quantity	Designation	Part number	Location/Remark
1	O ₂ sensor	6850645	Replace O ₂ sensor when the display "O ₂ measurement fault !!!" appears and when calibration is no longer possible.
1	Room air filter	8412384	Replace room air filter after no more than one year.
1	Cooling air filter	8412384	Replace cooling air filter after no more than one year.
1	Dust protection filter	8415978	Replace dust protection filter after no more than one year. The dust protection filter is fitted in devices from date of manufacture 2005.

Action • Check whether the above-mentioned maintenance parts have been replaced.

Result **The above-mentioned maintenance parts have been replaced.**

[☐ OK]

No.0306_000007424

3 Electrical safety

This section contains tests which have to be performed in order to determine the operational readiness of the medical-electrical system.

3.1 Electrical safety to DIN EN 62353 (IEC 62353)

NOTE

The device conforms to the conditions of protection class I, type BF.

NOTE

The tester, e.g. SECUTEST, must be correctly configured for all measurements. If implausible measurement results are obtained, such as a leakage current of 0.0 μ A, check the tester configuration in addition to the test setup!

NOTE

In testing to IEC 62353, the medical electrical device (ME device) or the medical electrical system (ME system) must be tested.

ME systems must be treated like ME devices.

An ME system is a combination of several devices, as specified by the manufacturer, of which at least one must be an ME device, which are interconnected by a functional connection or by means of a multiple socket outlet.

NOTE

With devices that are connected to other devices by means of a data cable, this connection must be disconnected prior to performing the electrical safety check, in order to avoid incorrect measurements.

3.1.1 Visual check

Prerequisites The tester and the device under test are switched off.

WARNING

Hazardous voltage.

Touching live components can lead to serious injury or death.

- Disconnect the power cord from the AC outlet before checking the power fuse-links.

- | | |
|--------|---|
| Action | • Disconnect the power plug from the mains socket. |
| Test | – The power fuse-links of the device under test match the specifications on the rating plate. |
| | – The power cable and plug are not dirty or damaged. |

Result **Condition checked.**

[_____ OK]

No.0306_0000007424

Evita 4 / Evita 4 edition
Electrical safety

3.1.2 Protective earth resistance

Test set-up




Fig. 8 Protective earth conductor resistance

Item	Designation
1	Tester
2	Device under test
3	Tester probe cable
L	Conductor
N	Neutral conductor
PE	Protective earth

- Action
- Prepare the test setup.
 - Switch the tester on.
 - Configure the tester appropriately and follow the instructions on the tester.
 - Using the tip of the probe cable, scan each of the measuring points on the device under test listed under "Protective earth resistance measuring points" one after the other, moving the mains power supply cord along the entire length during the measurement. The resistance must not change when you do so.

Test The protective earth resistance of devices with detachable but connected mains power supply cords must not exceed **0.3 Ohm** in each case.

Result **Maximum measured value of device with power supply cord.**

[_____] Ω

Test If other optional power supply cords are fitted, the respective protective earth resistance must not exceed **0.1 Ohm**. Move the power supply cord along the entire length during the measurement. The resistance must not change when you do so.

Result **Maximum measured value of optional power supply cord.**

[_____] Ω

No.0306_000007424

3.1.3 Protective earth resistance measuring points

- Action
- Scan the following measurement points for protective conductor resistance measurement one after the other using the tip of the probe cable:
 - Equipotential equilization pin on power pack
 - Gas connection and metal connectors on the pneumatic assembly
 - Inspiratory block
 - Side rails on the device under test
- Result Measurement points scanned.

[_____ OK]

3.1.4 Test adapters for equipment leakage current

- Prerequisites The following device-specific test adapters are required for the device under test:
- 7910195, Patient adapter
 - 7910364, Measuring line, leakage current, 2-pin, VDE test connector, temperature
 - 7901772, VDE test connector Babylog (optional)
 - 7901068, Shorting cable Nellcor, type L (optional)
- Action
- Connect all adapter cables to the device under test.
 - Connect all adapter cables to the sockets for applied parts on the tester.
 - Perform the device leakage current test.
- Result **Test adapters for device leakage current**

[_____ OK]

3.1.5 Equipment leakage current**NOTE**

The equipment leakage current can be tested by the differential measurement method or the direct measurement method.

In direct measurement, set up the device under test with insulation and scan all accessible conductive components using the probe (the protective earth is internally interrupted in the tester).

- Prerequisites The tester is switched on.

Evita 4 / Evita 4 edition
Electrical safety

Test set-up




Fig. 9 Equipment leakage current

Item	Designation
1	Tester
2	Device under test
3	(Applied part) Device-specific test adapter for tester
4	(Applied part) Configurable ports for applied part
L	Conductor
N	Neutral conductor
PE	Protective earth

- Action
- Prepare the test setup.
 - (Applied part) Connect the device-specific test adapter on one end to the device under test and on the other end to the tester's configurable port "A" for applied parts (paying attention to the configuration!).
 - Follow the instructions on the tester.

NOTE

For symmetrical mains plugs that have no preferential position in the socket-outlet, perform the test twice! Perform the second test with the plug rotated 180° in the socket. In many test devices the mains plug rotation is simulated by means of a built-in selector switch. Document the higher measured value.

NOTE

Always enter the reference value (initial value measured) in the "Test Report" or "Result Sheet" document.

NOTE

If the measured values are between 90% and 100% of the permissible limit value, apply the reference value and the previously measured values of the recurrent test to assess electrical safety!

Test The reference value must not exceed **500 μA**.

Result **Reference value**

[_____ μA]

Test The recurrent test value must not exceed **500 μA**.

Result **Recurrent test**

[_____ μA]

No.0306_000007424

3.1.6 Test adapters for patient leakage current

- Prerequisites** The following device-specific test adapters are required for the device under test:
- 7910195, Patient adapter
 - 7910364, Measuring line, leakage current, 2-pin, VDE test connector, temperature
 - 7901772, VDE test connector Babylog (optional)
 - 7901068, Shorting cable Nellcor, type L (optional)
- Action**
- Connect all adapter cables to the device under test.
 - Connect all adapter cables to the sockets for applied parts on the tester.

NOTE

Measure applied parts individually!

- Perform patient leakage current test.

Result **Test adapters for patient leakage current**

[_____ OK]

3.1.7 Leakage current, mains on applied part with airway temperature sensor AWT 01, if installed**NOTE**

In the following test the leakage current is measured at the airway temperature sensor AWT01. The expected value is very low (the typical measured value is 1.5 μ A to 2 μ A).

- Prerequisites** The tester is switched on.
- Action**
- Prepare the following test setup.
- Test set-up**




Fig. 10 Leakage current of airway temperature sensor

Item	Designation
1	Tester (test device)
2	Device under test
3	Connection port of device under test
4	Airway temperature sensor AWT 01
5	Test clip with measuring lead
6	Configurable ports for applied part

Evita 4 / Evita 4 edition
Electrical safety

Item	Designation
L	Conductor
N	Neutral conductor
PE	Protective earth

- Action
- (Applied part) Connect the airway temperature sensor on one end to the device under test and on the other end by a test clip with measuring lead and the tester, configurable port "A" for applied parts (paying attention to the configuration!).
 - Follow the instructions on the tester.

NOTE

For symmetrical mains plugs that have no preferential position in the socket-outlet, perform the test twice! Perform the second test with the plug rotated 180° in the socket. In many test devices the mains plug rotation is simulated by means of a built-in selector switch. Document the higher measured value.

NOTE

Always enter the reference value (initial value measured) in the "Test Report" or "Result Sheet" document.

NOTE

If the measured values are between 90% and 100% of the permissible limit value, apply the reference value and the previously measured values of the recurrent test to assess electrical safety!

Test The reference value must not exceed **5000** μA .

Result **Reference value**

[_____ μA]

Test The recurrent test value must not exceed **5000** μA .

Result **Recurrent test**

[_____ μA]

3.1.8 Leakage current on applied part with test adapter (normal condition)

NOTE

The following measurement is performed in the "normal condition" setting on the tester.

Prerequisites The tester is switched on.

No.0306_000007424

Test set-up




Fig. 11 Leakage current on applied part

Item	Designation
1	Tester (test device)
2	Device under test
3	Applied part of device under test
4	Test adapter
5	Configurable ports for application components
L	Conductor
N	Neutral conductor
PE	Protective earth

- Action
- Prepare the test setup.
 - (Applied part) Connect the device-specific test adapter on one end to the device under test and on the other end to the tester's configurable port "A" for applied parts (paying attention to the configuration!).
 - Follow the instructions on the tester.

NOTE

For symmetrical mains plugs that have no preferential position in the socket-outlet, perform the test twice! Perform the second test with the plug rotated 180° in the socket. In many test devices the mains plug rotation is simulated by means of a built-in selector switch. Document the higher measured value.

NOTE

Always enter the reference value (initial value measured) in the "Test Report" or "Result Sheet" document.

NOTE

If the measured values are between 90% and 100% of the permissible limit value, apply the reference value and the previously measured values of the recurrent test to assess electrical safety!

Airway temperature sensor

- Test The initial value must not exceed 100 µA "IAC".
- Result **Initial measured value "IAC", airway temperature sensor**

[_____ µA]

Evita 4 / Evita 4 edition
Electrical safety

Test	The initial value must not exceed 10 μA "IDC".	
Result	Initial measured value "IDC", airway temperature sensor	[] μA
Test	The recurrent test value must not exceed 100 μA "IAC".	
Result	Recurrent test "IAC", airway temperature sensor	[] μA
Test	The recurrent test value must not exceed 10 μA "IDC".	
Result	Recurrent test "IDC", airway temperature sensor	[] μA
NeoFlow option		
Test	The initial value must not exceed 100 μA "IAC".	
Result	Initial measured value "IAC", NeoFlow sensor	[] μA
Test	The initial value must not exceed 10 μA "IDC".	
Result	Initial measured value "IDC", NeoFlow sensor	[] μA
Test	The recurrent test value must not exceed 100 μA "IAC".	
Result	Recurrent test "IAC", NeoFlow sensor	[] μA
Test	The recurrent test value must not exceed 10 μA "IDC".	
Result	Recurrent test "IDC", NeoFlow sensor	[] μA

3.1.9 Optional SpO₂, leakage current, mains on applied part

Prerequisites The tester is switched on.

Test set-up




Fig. 12 Leakage current on applied part

Item	Designation
1	Tester (test device)
2	Device under test
3	Application component, device-specific test adapter for tester
4	Configurable connection sockets for applied part (A - K)

No.0306_000007424

Item	Designation
L	Conductor
N	Neutral conductor
PE	Protective earth

- Action
- Prepare the test setup.
 - Connect the device-specific test adapter on one end to the device under test and on the other end to the tester's configurable socket "A" for applied parts (paying attention to the configuration!).
 - Follow the instructions on the tester.

NOTE

For symmetrical mains plugs that have no preferential position in the socket-outlet, perform the test twice! Perform the second test with the plug rotated 180° in the socket. In many test devices the mains plug rotation is simulated by means of a built-in selector switch. Document the higher measured value.

NOTE

Always enter the reference value (initial value measured) in the "Test Report" or "Result Sheet" document.

NOTE

If the measured values are between 90% and 100% of the permissible limit value, apply the reference value and the previously measured values of the recurrent test to assess electrical safety!

Test Measure leakage current at SpO₂ connection. The reference value must not exceed **5000** µA.

Result **SpO₂ connection reference value**

[_____ µA]

Test Measure leakage current at SpO₂ connection. The recurrent test value must not exceed **5000** µA.

Result **SpO₂ connection recurrent test**

[_____ µA]

3.2 Electrical safety according to IEC 60601-1

NOTE

The device conforms to the conditions of protection class I, type BF.

NOTE

The following presents a description of the electrical safety tests according to IEC 60601-1. Whether the standard is applicable or not depends on national regulations and its use must be decided on site under consideration of applicable national regulations.

Evita 4 / Evita 4 edition
Electrical safety

NOTE

When testing according to IEC 60601-1, systems must be tested and not individual devices.

Systems should be treated as devices.

A system is a combination of several devices of which at least one is a medical electrical device which is connected to other devices by functional connections or by a transportable multiple socket outlet.

NOTE

With devices that are connected to other devices by means of a data cable, this connection must be disconnected prior to performing the electrical safety check, in order to avoid incorrect measurements.

NOTE

The test to UL2601-1 is satisfied by testing to IEC 60601-1. Differing limit values are marked.

3.2.1 Visual check

Prerequisites The tester and the device under test are switched off.

WARNING

Hazardous voltage.

Touching live components can lead to serious injury or death.

- Disconnect the power cord from the AC outlet before checking the power fuse-links.

Action • Disconnect the power plug from the mains socket.

Test – The power fuse-links of the device under test match the specifications on the rating plate.
– The power cable and plug are not dirty or damaged.


Result **Condition checked.**

[_____ OK]

No.0306_000007424

3.2.2 Protective earth resistance

Test set-up

**Fig. 13** Protective earth conductor resistance

Item	Designation
1	Tester
2	Device under test
3	Tester probe cable
L	Conductor
N	Neutral conductor
PE	Protective earth

- Action**
- Prepare the test setup.
 - Switch the tester on.
 - Configure the tester appropriately and follow the instructions on the tester.
 - Using the tip of the probe cable, scan each of the measuring points on the device under test listed under "Protective earth resistance measuring points" one after the other, moving the power supply cord along the entire length during the measurement. The resistance must not change when you do so.

Test The protective earth resistance of single devices (including the permanently installed power supply cord, or power supply cord detachable only by use of tools) must not exceed **0.2 Ohm** in each case.

Result **Maximum measured value.**

[_____ Ω]

3.2.3 Protective earth resistance measuring points

- Action**
- Scan the following measuring points for protective conductor resistance measurement one after the other using the tip of the probe cable:
 - Power supply unit potential equalization pin
 - Gas connection and metal connectors on the pneumatic assembly
 - Inspiratory block
 - Side rails on the device under test

Result Measurement points scanned.

[_____ OK]

Evita 4 / Evita 4 edition
Electrical safety

3.2.4 Earth leakage current

Test set-up

NOTE

To avoid false measurements, install the device under test insulated.




Fig. 14 Earth leakage current

Item	Designation
1	Tester (test device)
2	Device under test
L	Conductor
N	Neutral conductor
PE	Protective earth

- Action
- Prepare the test setup.
 - Switch on the tester and the device under test.

NOTE

For symmetrical mains plugs that have no preferential position in the socket-outlet, perform the test twice! Perform the second test with the plug rotated 180° in the socket. In many test devices the mains plug rotation is simulated by means of a built-in selector switch. Document the higher measured value.

- Follow the instructions on the tester.

Test Normal condition (N.C.): The value must not exceed **500** μA .

Result **Normal condition (N.C.)**

[_____ μA]

Test Single fault condition (S.F.C.): The value must not exceed **1000** μA .

Result **Single fault condition (S.F.C.)**

[_____ μA]

3.2.5 Patient leakage current

Prerequisites The tester is switched on.

No.0306_000007424

Test set-up




Fig. 15 Patient leakage current test setup

Item	Designation
1	Tester
2	Device under test
3	Applied part, device-specific test adapters for tester
4	Ports for applied parts
L	Conductor
N	Neutral conductor
PE	Protective earth

Prerequisites The following device-specific test adapters are required for the device under test:

- 7910195, Patient adapter
- 7910364, Measuring line, leakage current, 2-pin, VDE test connector, temperature
- 7901772, VDE test connector Babylog (optional)
- 7901068, Shorting plug Nellcor, type L

- Action**
- Prepare the test setup.
 - Plug the test adapter into the relevant connection on the device under test and connect the laboratory plug to the tester.
 - Follow the instructions on the tester.

NOTE

For symmetrical mains plugs that have no preferential position in the socket-outlet, perform the test twice! Perform the second test with the plug rotated by 180° in the power socket. In many test devices the mains plug rotation is simulated by means of a built-in selector switch. Document the higher measured value.

Test Measure leakage current at temperature sensor connection. In normal condition (N.C.), the value must not exceed **100 μA**.

Result **Temperature sensor connection, normal condition (N.C.)**

[_____ μA]

Evita 4 / Evita 4 edition
Function and condition test

Test	Measure leakage current at temperature sensor connection. In single fault condition (S.F.C.), the value must not exceed 500 μA .	
Result	Temperature sensor connection, single fault condition (S.F.C.)	[<input type="text"/> μA]
Test	Measure leakage current at Neo-Flow connection (optional). In normal condition (N.C.), the value must not exceed 100 μA .	
Result	NeoFlow connection, normal condition (N.C.)	[<input type="text"/> μA]
Test	Measure leakage current at Neo-Flow connection (optional). In single fault condition (S.F.C.), the value must not exceed 500 μA .	
Result	NeoFlow connection, single fault condition (S.F.C.)	[<input type="text"/> μA]
Test	Measure leakage current at SpO ₂ connection (optional). In normal condition (N.C.), the value must not exceed 100 μA .	
Result	SpO₂ connection, normal condition (N.C.)	[<input type="text"/> μA]
Test	Measure leakage current at SpO ₂ connection (optional). In single fault condition (S.F.C.), the value must not exceed 500 μA .	
Result	SpO₂ connection, single fault condition (S.F.C.)	[<input type="text"/> μA]

4 Function and condition test

This section contains tests to establish whether the device and the accessories used conform to the stipulations of the Instructions for Use in terms of condition and function.

4.1 Condition of basic device

Prerequisites The device is fully assembled.

4.1.1 Labels and instructions for use

Test	<ul style="list-style-type: none"> Labels are complete and legible. The instructions for use are available, according to the user/owner. Only applicable in Germany: The medical product log is available (according to user/owner). 	
Result	Condition checked.	[<input type="text"/> OK]

No.0306_000007/424

4.1.2 Base unit

- Test
- The housing is neither contaminated nor damaged.
 - Markings/labelling are complete and legible.

Result **Condition and function checked.**

[_____ OK]

4.1.3 Trolley

- Test
- Trolley is undamaged, all screw connections are tight.
 - The castors are undamaged and firmly attached to the trolley.

Result **Condition checked.**

[_____ OK]

4.1.4 Compressed gas connecting hoses

- Test
- Plugs, screw connections and tubes/hoses are undamaged.
 - Plugs, screw connections and tubes/hoses correspond to country-specific coding.

Result **Condition and country-specific coding checked.**

[_____ OK]

4.2 Function check

4.2.1 PEEP/PIP valve calibration

Prerequisites The device is on standby.

- Action
- Connect the "AIR" and "O₂" compressed gas connecting tubes of the ventilator to the central gas supply.
 - Disconnect the ventilation hoses from the inspiratory socket and expiratory socket.
 - Remove the flow sensor.
 - Detach patient system.
 - Unscrew the diaphragm cover (Fig. 16/7) from the patient system and remove the diaphragm (Fig. 16/6).
 - Screw the diaphragm cover (Fig. 16/7) onto the patient system.
 - Seal off the expiratory socket (Fig. 16/5) with a rubber plug.
 - Remove the "flow sensor connection" lip seal (Fig. 16/8) and seal off the connection with a rubber plug.

Evita 4 / Evita 4 edition
Function and condition test




Fig. 16 Patient system, components

Legend to Fig. 16

Item	Designation
1	Filter
2	Non-return valve
3	Plug
4	Copper tube (observe preferred position!)
5	Expiratory socket
6	Diaphragm with sealing washer
7	Diaphragm cover with O-ring
8	Lip seal

- Install the patient system on the device.
- When using a disposable expiratory valve, a short-circuit can be created between the two lip seals into the connection housing as an alternative.
- Press the "Calib. config." key.
- Press the "Ventilation" key.
- Enter the following number code: 4 7 9 9 4 7 4 7.

Test Calibration starts.

NOTE

The calibration may take more than 1 minute!

The calibration has been completed successfully if a flow pulse of less than/equal to 0.2 seconds can be heard at the inspiratory port.

NOTE

If the calibration failed (audible flow impulse was approx. 1 s), repeat the calibration. If the calibration fails several times, replace the PEEP/PIP valve and repeat the calibration.

Action • Restore the patient system to its correct condition. Note: Make sure that the sealing to the flow sensor is fitted correctly and that the non-return valve D 5.1 is not broken.


Result **Calibration completed.**

[_____]OK

No.0306_000007424

4.2.2 Power-on test

Test set-up

**Fig. 17** Device ready for use

Item	Designation
1	Tubing system according to instructions for use
2	Breathing bag
3	Cable from temperature sensor

- Action**
- Connect the ventilator to compressed gas supply and to the mains power supply.
 - Connect an adult tubing system without humidifier.
 - Connect the "AWT temperature sensor", if fitted.
 - Connect the "CO₂ sensor", if fitted.
 - Connect the "NeoFlow sensor", if fitted.
 - Connect the "Remote Pad", if fitted
 - Switch on the lung ventilator.

- Test**
- The ventilator carries out its self-test.
 - At the end of the self-test all LEDs are lit for approx. 2 seconds (also on the optional Remote Pad, if fitted).
 - The information on the display must not be corrupted by faulty pixels.

Result Power-on test checked.

[_____]OK]

4.2.3 Test according to internal device checklist as per instructions for use**NOTE**If the "CO₂" option is fitted, wait 3 minutes!

- Action**
- Call up device check and perform all test steps in succession. Follow the instructions on the display.

No.0306_0000007424

Evita 4 / Evita 4 edition
Function and condition test

Test At the end of the device check, a result sheet will be displayed containing the following information:

- "Check symbol" = correct result
- "F" = faulty result
- "(- -)" = check not performed

If there are leaks, the maximum permissible leak rate is 100 mL/min, as testing is being performed with a tubing system with no humidifier.

NOTE

Perform the leak test with all expiratory valves present.

Result All tests successfully completed.

[] OK

4.2.4 Non-return valve (D5.1) in expiratory valve

NOTE

Check all the expiratory valves registered in the "Test Report" or on the "Result Sheet".

Prerequisites The ventilator is switched off. The expiratory valve is adapted to the ventilator.

Test set-up




Fig. 18 Test setup for non-return valve D5.1 in expiratory valve

Item	Designation
1	Pressure gauge, class 1.6, measuring range e.g. +12 mbar to –12 mbar
2	Silicone tube, e.g. 7 x 2.5 Si, natural colour
3	Syringe, e.g. 60 mL
4	Tube clamp
5	Catheter connector 5.5 mm (outer)
6	LUER-LOCK with LUER-LOCK adapter, connected to expiratory socket

No.0306_000007/24


- Action
- Prepare the test setup.
 - Use the syringe to generate a negative pressure of -7 mbar, and then reduce it to -4 mbar.
- Test
- After about 5 seconds, the vacuum for each expiratory valve tested is still at least -1 mbar.
- Result **Function checked.**

[_____ OK]

4.2.5 Non-return valve D3.2 (10 mbar)

Prerequisites The ventilator is switched off.

Test set-up

**Fig. 19** Test setup for non-return valve D3.2

Item	Designation
1	Pressure gage, class 1.6, measuring range e.g. $+12$ mbar to -12 mbar
2	Tube, e.g. 7 x 2.5 Si, natural colour
3	Syringe, e.g. 60 mL
4	Tube clamp
5	Catheter connector 5.5 mm (outer)
6	LUER-LOCK with LUER-LOCK adapter, connected to inspiratory socket

- Action
- Prepare the test setup.
 - Use the syringe to feed a volume of 50 mL into the inspiratory socket within 2 seconds to 4 seconds.
- Test
- The reading on the comparative manometer is **5 mbar to 10 mbar**.
- Result **Enter the measured value.**

[_____ mbar]

- Action
- Withdraw the syringe.

4.2.6 Emergency air valve Y3.1/non-return valve D3.1

Prerequisites The ventilator is switched off.

No.0306_0000007424

Evita 4 / Evita 4 edition
Function and condition test

Test set-up

NOTE

Use an adult ventilation tube (0.9 m to 1.2 m long) between the inspiratory socket of the ventilator and the Y-piece (pressure measurement point).




Fig. 20 Test setup; emergency air valve Y3.1/non-return valve D3.1

Item	Designation
1	Adjustable compressed gas source, e.g. test pressure regulator
2	Ejector
3	Adult corrugated tube, e.g. 0.35 m long
4	Flowmeter; class 2.5, measuring range e.g. 10 L/min to 120 L/min
5	Manometer, class 1.6, measuring range e.g. +12 mbar to -12 mbar
6	Adult corrugated tube, 0.9 m to 1.2 m long
7	Y-piece
8	Adult corrugated hose, 0.9 m to 1.2 m long, plugged onto inspiratory socket

- Action
- Prepare the test setup.
 - With a test pressure regulator and an injector, set a flow of 55 L/min to 60 L/min.

Test – The measured value on the reference manometer is **0** to **-6** mbar.

Result **Enter the measured value.**


[_____ mbar]

- Action
- Remove the test setup.

No.0306_000007424

4.2.7 Safety valve D3.3 (120 mbar)

Test set-up

**Fig. 21** Test setup; safety valve D3.3

Item	Designation
1	Pressure gage, class 1.6, measuring range e.g. 0 mbar to 120 mbar
2	Silicone hose, e.g. 7 x 2.5 Si, natural color
3	Syringe, at least 60 mL
4	Tube clamp
5	Catheter connector 5.5 mm (outer)
6	LUER-LOCK with LUER-LOCK adapter, connected to inspiratory socket

NOTE

Applies to Evita 4/2 dura as from software version 4.2n:

If the ventilator is switched to Standby mode immediately after completion of the self-test, it will enter the safety switchoff/gas saving mode, and no pressure will be built up (Y1.3 is not actuated).

Remedy: start the "Leakage" device test and then interrupt it at once.

NOTE

In this test, the "Pressure measurement inop" alarm may be delivered; ignore this alarm. This alarm must disappear again when the test setup is removed.

Pressure may also build up as a result of leakage at the mixer if this was not zeroed at power-on. The test can be performed immediately, however.

- Action
- Prepare the test setup.
 - Set the device to "Standby" mode and confirm. From the "Check" menu, just start the leakage test using the "Check" softkeys and then stop it immediately using the "Abort" softkey.
 - Use the syringe to feed a volume of 50 mL into the inspiratory socket within 2 seconds to 4 seconds.

Evita 4 / Evita 4 edition
Function and condition test

- Test – The reading on the comparative manometer is **101** mbar to **110** mbar.
- Result **Enter the measured value.** [_____mbar]
- Action • Remove the test setup.

4.2.8 DC power supply unit with LEDs

NOTE

This test relates only to the new DC power supply unit 8419616.

- Prerequisites The ventilator is connected to the mains power supply.
The ventilator is switched off.
- Test The mains power LED on the rear of the power supply unit is lit green.
The LED for the optional external rechargeable battery on the rear of the power supply unit is lit green (rechargeable battery fully charged) or yellow (rechargeable battery charging). If no external battery is connected, the LED remains unlit.
The LED for the internal rechargeable battery on the rear of the power supply unit is lit green (rechargeable battery fully charged) or yellow (rechargeable battery charging). An internal battery must always be present.
- Action • Switch on the device.
- Test All LEDs on the rear of the power supply unit are out.
After the self-test the mains power supply connection is indicated by a yellow dot next to the power icon in the bottom right corner of the screen.
The state of charge of the internal and external rechargeable batteries is indicated in the bottom right corner of the screen as follows:
– Green battery icon = rechargeable battery fully charged.
– "Transparent" battery icon = rechargeable battery charging.
- Action • Cut the mains power supply.
- Test If the "external battery" option is connected, the "External battery active!!!" alarm is generated.
The yellow dot next to the on-screen mains power icon disappears.
The on-screen icon for the external rechargeable battery becomes transparent. A yellow dot is displayed next to the external rechargeable battery icon.
- Test If the "external battery" option is not present, or if the connector from the external rechargeable battery to the power supply unit is unplugged, the "Internal battery active!!!" alarm is generated.
The yellow dot next to the on-screen mains power icon has disappeared.
The on-screen icon for the internal rechargeable battery becomes transparent. A yellow dot is displayed next to the internal rechargeable battery icon.
- Result **Functions checked.**

[_____OK]

No.0306_000007424

4.2.9 DC module (if present)**NOTE**

This test item applies to the following power supply units:
 AC power supply unit 8306520 with retrofitted DC module 8306530
 For the above variant standard ex factory
 For the new DC power supply unit 8419616 with LEDs

NOTE

The internal battery may take up to 20 minutes to discharge. If the battery life and data backup tests are started here and run in parallel with the following test steps, test time can be considerably reduced.

- Action
- Carry out the test step [Internal rechargeable battery](#) in the test "DC module (if present), additional tests".

NOTE

It is possible that the device may switch off during one of the following test steps due to the internal battery discharging.

In this case first perform the following step and then repeat the interrupted test item.

- Action
- Carry out the test steps [Mains power failure and data backup](#) and ["External rechargeable battery" option](#) in the test "DC module (if present), additional tests".

Result **Function checked.**

[_____ OK]

4.3 Testing ventilation in adult mode

- Prerequisites
- The ventilator is in adult mode.
 - An adult tubing system is connected.

Evita 4 / Evita 4 edition
Function and condition test

Test set-up




Fig. 22 Sensor calibration

Item	Designation
1	External manometer e.g. DPI 705, 7910722
2	Silicone tube, e.g. 1198343
3	External O ₂ meter, e.g. Miniox 3000, 7911610
4	O ₂ connecting cable Miniox 3000
5	O ₂ measuring adapter (scope of supply, Miniox 3000)
6	Adult ventilation tube
7	LUER-LOCK adapter, 8290285 LUER-LOCK, 8600214
8	Test lung, 8403201

NOTE

In the following tests, the "IPPV" and "BIPAP/PCV+" ventilation modes are used in adult and paediatric mode.

If one of these modes has not been programmed by the user, this can be configured for the tests in accordance with the instructions for use.

Before handing the ventilator back to the user, restore the original user configuration.

No.0306_000007424

- Action
- Make the following settings:
 - Flow trigger = off
 - Autoflow = off
 - Sigh = off
 - Mode = IPPV
 - O₂ concentration = 21 %
 - VT = 500 mL
 - Tinsp = 2 s
 - f = 10 1/min
 - Flow = 25 L/min
 - Pmax = max.
 - PEEP = 5 mbar
 - On the "Alarm Limits" menu:
 - Alarm limits: Vti = max.; Paw = max., MVlow < 5L/min, etCO₂ = 0
 - Start ventilation.

4.3.1 Calibrating the sensors

Calibrate the sensors during ventilation.

- Action
- Press the "Calibration" key.
 - Press the "Flow" softkey.
 - Press the "Start cal." softkey.
- Test
- The device calibrates the flow sensor.
- Action
- Press the "O₂" softkey.
- Test
- The device calibrates the O₂ sensor.
- Result
- Calibration of sensors successful.**

[_____]OK]

4.3.2 O₂ comparative measurement

- Action
- Set the individual O₂ presets on the device and after a few seconds compare the measured values.
- Test
- Compare the device measured value on the "Values" menu with the external measured value.

O ₂ preset	O ₂ measurement from device/external tester
21 vol. %	19 to 23 vol. %
30 vol. %	27 to 33 vol. %
90 vol. %	85 to 95 vol. %

- Result
- O₂ comparative values within tolerance.**

[_____]OK]

Evita 4 / Evita 4 edition
Function and condition test

4.3.3 Volume metering/flow measurement

NOTE

If the ambient pressure is above 965 mbar, perform the "Measurement at high ambient pressure" test.

If the ambient pressure is below 965 mbar, perform the "Measurement at low ambient pressure" test.

Measurement at high ambient pressure

- Action
- Open the "Service diagnosis" start screen and enter the following codes:
 - 3999 and 3948 (leakage compensation off)
 - 3999 and 3958 (conversion off)
 - 7299 and 7235 and 0000 (set the tube compliance correction to "0").

NOTE

All settings remain active until the ventilator is switched off.

NOTE

Any fault 00.71.002 which occurs can be ignored.

- Action
- Start ventilation.
- Test
- Compare the specified "VT" value with the measured "VTe" value.
- Action
- Set the O₂ concentration to 21%.
 - Calibrate the flow sensor.
 - Read out "VTe".
- Test
- VTe = 450 to 550 mL
- Action
- Set the O₂ concentration to 100%.
 - Calibrate the flow sensor.
 - Read out "VTe".
- Test
- VTe = 450 to 550 mL (Note: The oxygen concentration rises slowly!).
- Result
- Measurement at high ambient pressure. Function checked.**

[_____]OK

Measurement at low ambient pressure

- Action
- Open the "Service diagnosis" start screen and enter the following codes:
 - 3999 and 3948 (leakage compensation off)
 - 7299 and 7235 and 0000 (set the tube compliance correction to "0").
 - 4655 (Start service mode).
- Action
- Under "Pneumatic", "Sensors" read the two factors "NTPS -> BTPS" and "BTPS -> NTPD".
 - The VTe readout must be corrected on all ventilators by the humidity error according to the following formula:

No.0306_000007424

$$V_{Te} = \frac{V_{Te \text{ display value}}}{(\text{Factor NTPS} \rightarrow \text{BTPS}) \times (\text{Factor BTPS} \rightarrow \text{NTPD})}$$

18/35

Fig. 23 Formula for correction of humidity error

- Set the O₂ concentration to 21%.
- Calibrate the flow sensor.
- Read out "VTe".

Test VTe = 450 to 550 mL

- Action
- Set the O₂ concentration to 100%.
 - Calibrate the flow sensor.
 - Read out "VTe".

Test VTe = 450 to 550 mL

Result **Measurement at low ambient pressure. Function checked.**

[_____ OK]

- Action
- Switch the device to standby mode.

4.3.4 BIPAP/PCV+ ventilation mode

- Action
- Select "BIPAP/PCV+" ventilation mode.
 - Make the following entries, confirming each with the rotary knob:
 - O₂ concentration = 21%
 - T_{insp} = 2 s
 - f = 10 1/min
 - Ramp = 0.2 s
 - P_{insp} = 25 mbar
 - PEEP = 5 mbar
 - PASB = 0 mbar
 - Tube compensation option off.
 - Alarm limits: set V_{ti} and P_{aw} to "max".
 - Start the ventilation.

Plausibility checking of pressure measurement and control

- Action
- Read out measured values.

Test Read the PEEP measurement from the device and compare it with the external measured value.

Device measured values	External measured values
PEEP = 4 to 6 mbar (during expiratory phase)	PEEP = 4 to 6 mbar (during expiratory phase)
P _{peak} = 23 to 27 mbar (during inspiratory phase)	P _{peak} = 23 to 27 mbar (during inspiratory phase)

- Action
- Set "PEEP" to 1 mbar.
 - Set "P_{insp}" to 50 mbar.
 - Read out measured values.

Test Compare measured value with external measuring device:

No.0306_000007424

Evita 4 / Evita 4 edition
Function and condition test

Device measured values	External measured values
PEEP = 0 to 2 mbar (during expiratory phase)	PEEP = 0 to 2 mbar (during expiratory phase)
Ppeak = 47 to 53 mbar (during inspiratory phase)	Ppeak = 47 to 53 mbar (during inspiratory phase)

Result **Functions checked.**

[☐ OK]

Action • Switch the device to standby mode.

4.4 Additional function tests

4.4.1 Power failure and data backup (ventilators without DC module)

- Action • Start the ventilation.
• Disconnect the ventilator power plug from the mains socket outlet.
- Test – The acoustic power failure warning sounds.
– The signal can be clearly heard distinguished from the general noise level.
- Action • Connect the power plug of the ventilator to the mains power supply.
- Test – Following a shortened version of the self-test, the ventilator restarts operation with the previously set ventilation mode.

Result **Function checked.**

[☐ OK]

4.5 DC module or Benning power supply unit (if fitted), additional tests

NOTE

The tests for the "DC module" are run in parallel with other function tests.

NOTE

This test item applies to the following power supply units:

AC power supply unit 8306520 with retrofitted DC module 8306530

For the above variant standard ex factory

For the new DC power supply unit 8419616 with LEDs

4.5.1 Internal rechargeable battery

NOTE

Power supply units with a DC module are always fitted with internal rechargeable batteries.

Prerequisites The ventilator is connected to the compressed gas and mains power supply.

No.0306_000007424

- Action
- Unplug the connector of the "external rechargeable battery" option from the power supply unit.
 - Switch on the device.
 - Set controlled ventilation for adults.
 - Disconnect the ventilator power plug from the mains socket-outlet.
- Test
- The DC connector is undamaged.
 - The ventilator continues operating without interruption.
- Action
- Continue operating the ventilator from the internal battery supply until the battery is depleted (with a fully charged battery this will take approx. 20 minutes).
- Result **Internal rechargeable battery, state of charge and function tested.**

[_____ OK]

4.5.2 Mains power failure and data backup

- Test
- The audible power failure alarm sounds.
 - The signal can be clearly heard from a distance of 3 m and distinguished from the general noise level.
- Action
- Connect the power plug of the ventilator to the mains power supply.
- Test
- Following a shortened version of the self-test, the ventilator must restart operation with the previously selected ventilation mode.
- Result **Power failure and data backup, function tested.**

[_____ OK]

4.5.3 "External rechargeable battery" option

- Action
- Connect the cable connector of the external rechargeable batteries to the power supply unit.
 - Disconnect the ventilator power plug from the mains socket-outlet.
- Test
- The ventilator continues operating without interruption.
 - The power supply from the external rechargeable batteries is indicated on the display.
- Result **External rechargeable battery, function tested.**

[_____ OK]

4.6 CO₂ sensor option

- Prerequisites The ventilator is connected to the compressed gas and mains power supply and switched on. The internal device check has been carried out.

4.6.1 Condition

- Prerequisites The ventilator is connected to the compressed gas and mains power supply and switched on. The internal device check has been carried out.

No.0306_0000007424

Evita 4 / Evita 4 edition
Function and condition test

- Test
- The housing and the anti-kink sleeve are undamaged.
 - The cable and the connector are undamaged.
 - The window in the CO₂ sensor is neither contaminated nor damaged.
 - The mounting of the test filter on the connecting cable is undamaged. The test filter is neither contaminated nor scratched.
 - The cuvette windows are neither contaminated nor damaged.
 - The windows in the standby holder for the CO₂ sensor are neither contaminated nor damaged.

Result **Condition checked.**

[_____]OK]

4.6.2 Measurement accuracy with test filter

NOTE

Before testing with the test filter, perform a zero calibration of the CO₂ sensor on the standby holder or on a clean reusable cuvette.

NOTE

Calibration data is stored in the CO₂ sensor, not in the ventilator.

NOTE

The value stated on the test filter applies to the given CO₂ sensor and not to other sensors! The labeling on the test filter can be changed after calibration with test gas.

- Action
- Connect the breathing tube system without a humidifier to the ventilator.
 - Select and start "IPPV" mode.
 - If the "CO₂" function is not active, press the "Alarm limits" key.
 - Press the "Monitoring" softkey and push the rotary knob.
- Test
- The lamp in the CO₂ sensor comes on.
- Action
- Press the "Calibration" key.
 - Press the "CO₂" softkey.
 - Slot the CO₂ sensor onto the test filter.
 - Press the "Filter test" softkey.
- Test
- The value for "FCO₂" is indicated on the display. The value matches the figure on the test filter with a tolerance of +/-0.3 vol. %.
- Action
- Fit the CO₂ sensor in the standby holder.
- Result **Measurement accuracy with test filter. Function checked.**

[_____]OK]

4.7 Options

4.7.1 Evita SAT (SpO₂, hardware)

- Test The SpO₂ sensor and cable are undamaged.

No.0306_000007424

- Action
- Activate the SpO₂ measurement function on the device.
 - Connect the SpO₂ sensor to the device.
 - Clip the SPO₂ cable by the finger clip to your own finger.

- Test
- The displayed SpO₂ value is above 90%
 - The pulse rate is plausible.

Result **Condition checked.**

[_____]OK]

4.7.2 NeoFlow (hardware)

NOTE

The option is identifiable by the Flow PCB with NeoFlow sensor connection.
For this option an expiratory valve with water trap must be used!

Condition of accessories

- Prerequisites
- The ventilator is connected to the compressed gas and mains power supply and switched on.
 - The ventilator is in neonatal mode.
 - The water trap of the expiratory valve is undamaged.

- Test
- The flow sensor cable is undamaged. The socket on the flow sensor cable engages securely with the flow sensor insert.
 - The ISO 15 flow sensor and the flow sensor insert are undamaged.

Function of accessories

- Action
- Press the "Calibration" key.
 - Press the "NeoFlow" softkey.
 - Pull the NeoFlow sensor out of the Y-piece, seal it off, and confirm.

- Test
- The message "Calibration" will be displayed after a brief period.

Result **Condition and function checked.**

[_____]OK]

4.7.3 Breathing Support Package (software)

- Action
- Press the "Mode settings" key.

- Test
- The "Tubus Comp" and "PPS" softkeys are installed.

Result **Breathing Support Package (software)**

[_____]OK]

4.7.4 ATC (Automatic Tubus Compensation, software)

- Action
- Press the "Mode Settings" key and the "BIPAP", "Additional Setting" softkeys.

- Test
- If the "Tubus Comp" softkey is installed, the option is available.

Result **The "ATC" option is present.**

[_____]OK]

No.0306_0000007424

Evita 4 / Evita 4 edition
Function and condition test

4.7.5 CapnoPlus (hardware)

- Test – The option is present if the port for the CO₂ connection is installed on the CO₂ Carrier PCB.

Result **The "CapnoPlus" option is present.**

[_____]OK]

4.7.6 SW 4.00 Plus (software)

- Action • Press the "Calibration" key, then the "Flow" softkey.
Test – The "External flow source" softkey is installed.

Result **SW 4.00 Plus (software)**

[_____]OK]

4.7.7 NIV (Non-Invasive Ventilation, software)

- Test – If the "Tube/Mask" softkey is displayed in standby, the option is available.

Result **The "mask ventilation NIV" option is present.**

[_____]OK]

4.7.8 IFCO Carrier/NurseCall (hardware)

Prerequisites The ventilator is connected to the compressed gas and mains power supply and switched on.

Nurse call

- Action • Switch the device to standby mode.
Test The standby alarm is activated (do not deactivate with the Reset key!).
Action • Measure the ohmic resistance at the nurse call connector on the back of the ventilator.

Test set-up




Fig. 24 Contacts on the nurse call connector (top view)

- Test – Continuity between pin 3 and pin 5.
– No continuity between pin 1 and pin 3.
Action • Press the Reset key, and then measure the ohmic resistance at the nurse call connector on the back of the ventilator.

- Test – No continuity between pin 3 and pin 5.
– Continuity between pin 1 and pin 3.

Result **Nurse call. Condition checked.**

[_____]OK]

No.0306_000007424

4.7.9 Remote Pad (hardware)**Condition**

Prerequisites The ventilator is connected to the compressed gas and mains power supply and switched on.

- Test
- The housing is undamaged.
 - The keypad is undamaged.
 - The cable and the connector are undamaged.

Function

Action • Start up ventilator.

- Press the alarm silence key.

Test – The LED for this key comes on in the Remote Pad and on the ventilator.

Action • Slide the flow sensor to the left.

Test – The error message "Flow Sensor?" is displayed on red background at the top of the screen.

Action • Slide the flow sensor to the right.

Test – The background in the display changes from red to blue.

Action • Press the Alarm Reset key on the Remote Pad.

Test – The error message disappears from the display.

Action • Press the drug nebulizer key for approx. 2 s.

Test – The ventilator displays the message that drug nebulization is on.
– The LED for this key in the Remote Pad comes on.

Action • Press the "O₂ suction" key.

Test – The O₂ value rises (max. 100%).

Action • Press the "Insp. Hold" key.

Test – The airway pressure curve is held on the upper pressure level until this key is released.
– The LED on the Remote Pad is on as long as the key is pressed.

Action • Press the "Exp. Hold" key.

Test – The airway pressure curve is held on the lower pressure level until this key is released.
– The LED on the Remote Pad is on as long as the key is pressed.

Result **Condition and function checked.**

[_____]OK]

4.7.10 Communication between ventilator and monitor (hardware)**NOTE**

This test relates to the Infinity monitors SC 7000, SC 8000, SC 9000 XL, Delta, Delta XL and Kappa.

No.0306_0000007424

Evita 4 / Evita 4 edition Test equipment

- Action
- Connect between ventilator and monitor via MIB converter, preferentially at COM2 or COM3 if the "Evita Link" option is present.
 - Switch on the ventilator and the monitor.
- Test
- Communication between the two devices is present.
- Result
- Function checked.**

[_____ OK]

4.8 Final action

- Prerequisites
- The test instructions have been performed as specified.
 - All tests performed were passed successfully.

4.8.1 Device handover

CAUTION

- ▶ External battery not connected!
- ▶ If the red plug of the external battery is not connected to the back of the power supply unit, the ventilator's runtime will be reduced considerably in the event of a mains power failure. Check that the plug of the external battery is securely connected to the back of the power supply unit.

CAUTION

- ▶ Batteries not fully charged!
- ▶ If the batteries are not fully charged, the runtime of the ventilator in the event of a mains power failure may be reduced considerably. If applicable, inform the user that the rechargeable batteries still need to be recharged.

- Action
- Attach a test label to the device.
 - Supply the user/owner with a fully functioning device.

Result

Device handover

[_____ OK]

5 Test equipment

This section sets out the test equipment required for the tests in this test procedure.

5.1 Test equipment

5.1.1 Test equipment subject to mandatory calibration






NOTE

Tests must be performed using test equipment calibrated in a valid manner!


No.0306_000007424

NOTE

Use the following test equipment or equivalent aids.

Designation	Part number	Comments
Electrical safety tester, e.g. GMC Secutest	7910594	
Multimeter	7901021	
Manometer	7910722	
Flowmeter 120 L/min	7900718	
Adjustable compressed gas source (O ₂ test pressure regulator; 300 bar manometer, 20 bar manometer)	7901482	

Result **Test equipment calibrated in a valid manner used**[☐ OK]**5.1.2 Test equipment not subject to mandatory calibration**

Designation	Part number	Remarks
Tube 7 x 2.5 Si nat.	1198343	

No.0306_0000007424

Evita 4 / Evita 4 edition
Test equipment







Designation	Part number	Remarks
Rubber plug	1294113	
Tube clamp	7900484	
Ejector	7900930	
Shorting cable Nellcor, type L	7901068	
VDE test connector Babylog	7901772	
Patient adapter	7910195	
Measuring line, leakage current, 2-pin, temp.	7910364	
Syringe, 60 mL	7910745	

No.0306_000007424

Designation	Part number	Remarks
Probe cable	7911196	
USB adapter	7911718	
MX 300-i oxygen monitor	7911955	
LUER-LOCK adapter	8290285	
LUER-LOCK	8600214	
Breathing bag with 7 mm catheter connector, set	8403201	

No.0306_0000007424

Evita 4 / Evita 4 edition
Test equipment

Designation	Part number	Remarks
Catheter connector 5.5. mm (outer), set	8403684	
O ₂ adapter	8405807	
Y-piece	2M12754	
Corrugated tube	8402041 or equivalent	
Connecting sleeve	M13506	
Y-piece	M25650	

5.1.3 Additionally required items

Designation	Part number	Comments
Test label		No illustration available.
Tubing system	As specified in instructions for use	
"AIR" pipeline supply connecting tube	As specified in instructions for use	
"O ₂ " pipeline supply connecting tube	As specified in instructions for use	

No.0306_000007424

6 Annex

This section contains further information which may be helpful in testing the device.

6.1 Evita 4 service mode start screen

6.1.1 Access

Service mode can be called up at any time in the course of operation. This mode only involves data readout; settings are not altered.

1. Press the "Configuration" key.
2. Select the "Basic setting" menu.
3. Select the "Service Diagnosis" menu.




Fig. 25 Detail view of display; Service mode start screen

4. Enter the following number code: 4 6 5 5.
5. Select the required menu.

5.2 Result Sheet Test instructions / Service Card IPM



Result Sheet Test Instructions / Service Card IPM
Evita 4 / Evita 4 edition

Order number:

Location: _____

Serial no. (basic unit): _____

Department: _____

Cust. invent. no.: _____

Maintenance interval: _____

Other / Delivery date: _____

Key

- ✓ / OK = OK
- + = Spare part used
- ! = Error / Report
- / = Accessory not available
- = Not applicable

Applies to Test Instructions / Service Card IPM Revision 13.0

Test	Result
1 Device configuration	
1.1 Device configuration	
1.1.1 Serial number (if not otherwise recorded)	
<input type="checkbox"/> 1.1.1.1 basic device	txt
<input type="checkbox"/> 1.1.1.2 disposable expiratory valve	txt
<input type="checkbox"/> 1.1.1.3 expiratory valve 1	txt
<input type="checkbox"/> 1.1.1.4 expiratory valve 2	txt
<input type="checkbox"/> 1.1.1.5 control unit	txt
<input type="checkbox"/> 1.1.1.6 CO ₂ sensor	txt
<input type="checkbox"/> 1.1.2 Software version	txt
<input type="checkbox"/> 1.1.3 Operating hours	h
<input type="checkbox"/> 1.1.4 Device ID	txt
2 Maintenance parts	
2.1 Maintenance intervals, overview	
2.1.1 Maintenance intervals and required sets and parts	
2.2 Maintenance parts by specified interval	
<input type="checkbox"/> 2.2.1 Service set, once a year, MX08220	dat
2.2.2 Option for Canada and USA, filter element MP03903	
<input type="checkbox"/> 2.2.2.1 AIR filter	dat
<input type="checkbox"/> 2.2.2.2 O ₂ filter	dat
2.2.3 Gas inlet block	
<input type="checkbox"/> 2.2.3.1 Service set, gas inlet block FAS, MX08225	dat
<input type="checkbox"/> 2.2.3.2 Service set, gas inlet block Dräger, MX08226	dat
<input type="checkbox"/> 2.2.3.3 Service set, gas inlet block FAS, MX08227	dat
<input type="checkbox"/> 2.2.3.4 Service set, gas inlet block Dräger, MX08228	dat
<input type="checkbox"/> 2.2.4 Optional rechargeable battery (DC module), 8421988	dat
<input type="checkbox"/> 2.2.5 Optional rechargeable battery (trolley), 8422125	dat
2.2.6 Replacing the LCD backlight after 25000 operating hours	
<input type="checkbox"/> 2.2.6.1 Replacement of no longer available displays and backlights	txt
<input type="checkbox"/> 2.2.6.2 Maintenance-free display	OK
2.2.7 Maintenance parts as required	
<input type="checkbox"/> 2.2.7.1 Maintenance parts as per instructions for use	OK
3 Electrical safety	
3.1 Electrical safety to DIN EN 62353 (IEC 62353)	
<input type="checkbox"/> 3.1.1 Visual check	OK
3.1.2 Protective earth resistance	
<input type="checkbox"/> 3.1.2.1 Maximum measured value of device with power supply cord	Ω
<input type="checkbox"/> 3.1.2.2 Maximum measured value of optional power supply cord	Ω
<input type="checkbox"/> 3.1.3 Protective earth resistance measuring points	OK
<input type="checkbox"/> 3.1.4 Test adapters for equipment leakage current	OK

Test	Result
3.1.5 Equipment leakage current	
<input type="checkbox"/> 3.1.5.1 Reference value	µA
<input type="checkbox"/> 3.1.5.2 Recurrent test	µA
<input type="checkbox"/> 3.1.6 Test adapters for patient leakage current	OK
3.1.7 Leakage current, mains on applied part with airway temperature sensor AWT 01, if installed	
<input type="checkbox"/> 3.1.7.1 Reference value	µA
<input type="checkbox"/> 3.1.7.2 Recurrent test	µA
3.1.8 Leakage current on applied part with test adapter (normal condition)	
<input type="checkbox"/> 3.1.8.1 Initial measured value "IAC", airway temperature sensor	µA
<input type="checkbox"/> 3.1.8.2 Initial measured value "IDC", airway temperature sensor	µA
<input type="checkbox"/> 3.1.8.3 Recurrent test "IAC", airway temperature sensor	µA
<input type="checkbox"/> 3.1.8.4 Recurrent test "IDC", airway temperature sensor	µA
<input type="checkbox"/> 3.1.8.5 Initial measured value "IAC", NeoFlow sensor	µA
<input type="checkbox"/> 3.1.8.6 Initial measured value "IDC", NeoFlow sensor	µA
<input type="checkbox"/> 3.1.8.7 Recurrent test "IAC", NeoFlow sensor	µA
<input type="checkbox"/> 3.1.8.8 Recurrent test "IDC", NeoFlow sensor	µA
3.1.9 Optional SpO₂, leakage current, mains on applied part	
<input type="checkbox"/> 3.1.9.1 SpO ₂ connection reference value	µA
<input type="checkbox"/> 3.1.9.2 SpO ₂ connection recurrent test	µA
3.2 Electrical safety according to IEC 60601-1	
<input type="checkbox"/> 3.2.1 Visual check	OK
<input type="checkbox"/> 3.2.2 Protective earth resistance	Ω
<input type="checkbox"/> 3.2.3 Protective earth resistance measuring points	OK
3.2.4 Earth leakage current	
<input type="checkbox"/> 3.2.4.1 Normal condition (N.C.)	µA
<input type="checkbox"/> 3.2.4.2 Single fault condition (S.F.C.)	µA
3.2.5 Patient leakage current	
<input type="checkbox"/> 3.2.5.1 Temperature sensor connection, normal condition (N.C.)	µA
<input type="checkbox"/> 3.2.5.2 Temperature sensor connection, single fault condition (S.F.C.)	µA
<input type="checkbox"/> 3.2.5.3 NeoFlow connection, normal condition (N.C.)	µA
<input type="checkbox"/> 3.2.5.4 NeoFlow connection, single fault condition (S.F.C.)	µA
<input type="checkbox"/> 3.2.5.5 SpO ₂ connection, normal condition (N.C.)	µA
<input type="checkbox"/> 3.2.5.6 SpO ₂ connection, single fault condition (S.F.C.)	µA
4 Function and condition test	
4.1 Condition of basic device	
<input type="checkbox"/> 4.1.1 Labels and instructions for use	OK
<input type="checkbox"/> 4.1.2 Base unit	OK

Revision 13.0

Evita 4 / Evita 4 edition

1 / 3

Test	Result
<input type="checkbox"/> 4.1.3 Trolley	OK
<input type="checkbox"/> 4.1.4 Compressed gas connecting hoses	OK
4.2 Function check	
<input type="checkbox"/> 4.2.1 PEEP/PIP valve calibration	OK
<input type="checkbox"/> 4.2.2 Power-on test	OK
<input type="checkbox"/> 4.2.3 Test according to internal device checklist as per instructions for use	OK
<input type="checkbox"/> 4.2.4 Non-return valve (D5.1) in expiratory valve	OK
<input type="checkbox"/> 4.2.5 Non-return valve D3.2 (10 mbar)	mbar
<input type="checkbox"/> 4.2.6 Emergency air valve Y3.1/non-return valve D3.1	mbar
<input type="checkbox"/> 4.2.7 Safety valve D3.3 (120 mbar)	mbar
<input type="checkbox"/> 4.2.8 DC power supply unit with LEDs	OK
<input type="checkbox"/> 4.2.9 DC module (if present)	OK
4.3 Testing ventilation in adult mode	
<input type="checkbox"/> 4.3.1 Calibrating the sensors	OK
<input type="checkbox"/> 4.3.2 O ₂ comparative measurement	OK
4.3.3 Volume metering/flow measurement	
<input type="checkbox"/> 4.3.3.1 Measurement at high ambient pressure.	OK
<input type="checkbox"/> 4.3.3.2 Measurement at low ambient pressure.	OK
<input type="checkbox"/> 4.3.4 BiPAP/PCV+ ventilation mode	OK
4.4 Additional function tests	
<input type="checkbox"/> 4.4.1 Power failure and data backup (ventilators without DC module)	OK
4.5 DC module or Benning power supply unit (if fitted), additional tests	
<input type="checkbox"/> 4.5.1 Internal rechargeable battery	OK
<input type="checkbox"/> 4.5.2 Mains power failure and data backup	OK
<input type="checkbox"/> 4.5.3 External rechargeable battery option	OK
4.6 CO₂ sensor option	
<input type="checkbox"/> 4.6.1 Condition	OK
<input type="checkbox"/> 4.6.2 Measurement accuracy with test filter	OK
4.7 Options	
<input type="checkbox"/> 4.7.1 Evita SAT (SpO ₂ , hardware)	OK
<input type="checkbox"/> 4.7.2 NeoFlow (hardware)	OK
<input type="checkbox"/> 4.7.3 Breathing Support Package (software)	OK
<input type="checkbox"/> 4.7.4 ATC (Automatic Tubus Compensation, software)	OK
<input type="checkbox"/> 4.7.5 CapnoPlus (hardware)	OK
<input type="checkbox"/> 4.7.6 SW 4.00 Plus (software)	OK
<input type="checkbox"/> 4.7.7 NIV (Non-Invasive Ventilation, software)	OK
<input type="checkbox"/> 4.7.8 IFCO Carrier/NurseCall (hardware)	OK
<input type="checkbox"/> 4.7.9 Remote Pad (hardware)	OK
<input type="checkbox"/> 4.7.10 Communication between ventilator and monitor (hardware)	OK
4.8 Final action	
<input type="checkbox"/> 4.8.1 Device handover	OK
5 Test equipment	
5.1 Test equipment	
<input type="checkbox"/> 5.1.1 Test equipment subject to mandatory calibration	OK
5.1.2 Test equipment not subject to mandatory calibration	
5.1.3 Additionally required items	
6 Annex	
6.1 Evita 4 service mode start screen	
6.1.1 Access	

Report:

Test has been performed according to the test instructions.

Name: :

Date/Signature:..... :

This page has been left blank intentionally.



Directive 93/42/EEC concerning medical
devices



Manufacturer



Drägerwerk AG & Co. KGaA

Moislinger Allee 53 – 55

D-23542 Lübeck

Germany



+49 451 8 82-0

FAX

+49 451 8 82-2080



<http://www.draeger.com>

9036758 – 5664.500 en

© Drägerwerk AG & Co. KGaA

Revision: 6.0 –

Subject to change without notice.

Will not be replaced in the event of
modifications.

